## Container Filter with Indicator

### EN - INTENDED USE:

Key Surgical® Container Filters are accessory parts for a sterile barrier system (sterilization container) intended to be used to form a reliable barrier against microbes during the sterilization of devices (instruments) of suitable dimensions in sterilization containers.

### PRECAUTIONS:

Used container filters must be removed prior to the processing (cleaning) of the sterile container systems. A new filter must be inserted prior to a new sterilization of the container.

Single use only; reuse may lead to patient cross contamination.

### CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

### INSTRUCTIONS FOR USE:

1. Container filters may only be handled by qualified personnel to prevent damage. Damaged filters should not be used.
2. Select appropriate size. To ensure a safe germ barrier, check the container filter used on the respective container system for correct fit.
3. When inserting the container filters, the instructions for use of the manufacturer of the container systems must be observed to ensure correct handling of the various holder systems and their maintenance during filter change.
4. Consult the container manufacturer's IFU for container reprocessing information.
5. Verify the indicator changes from pink to dark brown/black as diagramed below.

### STERILIZATION:

Compatible with steam sterilization processes up to 134°C. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Color change process indicator from pink to dark brown/black

### STORAGE:

Store at +42°F (+5°C) to 86° F (+30°C), 30%rH bis 60%rH, avoid direct sunlight.

### DISPOSAL:

Dispose after use according to facility protocol.

### WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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**KEY SURGICAL LLC**
MDSS UK RP LTD