## Tray Liners and Soaker Sheets

### EN - INTENDED USE:

Key Surgical® Tray Liners and Soaker Sheets are intended for placement on the bottom of a sterilisation tray to absorb wetness and humidity during processing.

### CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

### PRECAUTIONS:

Store away from a source of ignition.  
Single use only; reuse may lead to patient cross contamination.

### PREPARATIONS FOR USE:

- Select appropriate size product for tray  
Note: Tray Liners may be cut to the required size to fit the selected tray
- Product may be placed inside and/or under the tray
- Place product inside the wrapped tray or rigid container system

TLP-xxx-xxx and TL-xxx-xxx products can be used in the following configurations:

Pre-Vacuum Steam/100% Ethylene Oxide (EtO)

- Up to 2 tray liners inside and 2 tray liners outside
- Up to 8 folded layers (4 tray liners)

### INSTRUCTIONS FOR USE:

- Sterilise according to the information below
- Place product inside or under tray
- Place instruments in tray and wrap or close container according to manufacturer's instructions

### STERILISATION:

Provided non-sterile and are single use.

Not suitable for vaporized hydrogen peroxide/plasma sterilisation.

Validated for sterilisation efficacy using the following cycle parameters:

Steam Sterilisation			
Item Number	Cycle	Temperature	Minimum Cycle
TLP-012-018 TLP-020-025 TL-101-200 TL-202-500 TL-203-500 TL-204-500 TL004A TL005A TL006A TL007A TL008A TL009A TL010A TL011A TL012A TL013A SOK001A SOK002A SOK003A SOK006A SOK008A SOK009A SOK010A	Pre-Vacuum	132°C (270°F)	3 minutes
		132°C (270°F)	4 minutes
		134°C (273°F)	3 minutes
TL-101-200 TL-202-500 TL-203-500 TL-204-500 TL004A TL005A TL006A TL007A TL008A TL009A TL010A TL011A TL012A TL013A SOK001A SOK002A SOK003A SOK006A SOK008A SOK009A SOK010A	Gravity	121°C (250°F)	30 minutes

100% Ethylene Oxide (EtO)				
Item Number	Concentration	Temperature	Exposure Time	Humidity
TLP-012-018 TLP-020-025 TL-101-200 TL-202-500 TL-203-500 TL-204-500 TL004A TL005A TL006A TL007A TL008A TL009A TL010A TL011A TL012A TL013A SOK001A SOK002A SOK003A SOK006A SOK008A SOK009A SOK010A	735 mg/L	55° C	60 minutes	50-70%

Items 1136011 and 1137011 are compatible for steam sterilization up to 134°C for 3.5 minutes or 121°C for 15 minutes.

All other products are compatible for steam sterilization up to 134°C for up to 18 minutes.

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

### DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

### WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
UK REP MDSS UK RP LTD