## Tray Mats

### **EN** - INTENDED USE:

Key Surgical® Tray Mats are intended to be placed under the tray inside the sterilization wrap for organization and storage.

### **CONTRAINDICATIONS:**

No known contraindications and/or adverse effects.

### **PRECAUTIONS:**

Store away from any sources of ignition.

Single use only; reuse may lead to patient cross contamination.

### **INSTRUCTIONS FOR USE:**

1. Does not need to be cleaned or disinfected before initial use
2. Choose a mat which is slightly larger than the underside of the tray or basket being wrapped.
3. Place the mat on a flat surface.
4. Place the mat on the sterilization wrap in the place where you intend to put the basket.
5. Place the basket or tray on the mat.
6. When using DRYMATxxx or DRYMATGRNxxx, ensure that the white side is uppermost, so it contacts the tray or basket.
7. Load basket or tray and wrap as normal.

### **STERILIZATION:**

Mats are provided non-sterile and are single use. Sterilization compatibility below is referenced and determined by the operator manuals for the sterilisation devices.

<b>STERILIZATION COMPATILITY</b>	
MATxxx	Steam V-PRO® STERRAD®
DRYMATxxx DTRYMATGRNxxx	Steam

Note – Consult sterilization device operator manual for material compatibility.

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

### **DISPOSAL:**

Used: Dispose after use according to facility protocol. Do not reuse.

Unused: Can be disposed of in the plastic recycling stream.

### **WARNING:**

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
UK REP MDSS UK RP LTD