

Name: See Through Point Protectors

Product Codes: PPT007C and PPT008C

Name of Manufacturer: Clinipak Ltd

Address: Beech House

Knives Beech Business Centre

Loudwater

High Wycombe

Bucks

HP10 9SD

Tel: +44 (0)1628 810626

Fax: +44 (0)1628 810760

Website: clinipak.co.uk

Intended Use

Plastic Point Protectors are simple-to-assemble Plastic-based products which are designed for protecting chisels and other sharp instruments during sterilisation and during transport and storage.

They are particularly useful when the user wishes to identify the instrument being protected.

See Through Point Protectors are intended to be used by professional healthcare personnel working in sterile services (or similar) departments.

Sterilisation Compatibility

See Through Point Protectors are suitable for use in steam (up to 137°C), ethylene oxide and plasma (including the Sterrad Process) sterilisation processes.

Storage Conditions

It is recommended that See Through Point Protectors are stored under the conditions stated on the label of the box. If this is not possible, they should be prevented from getting wet or from being exposed to high temperatures (>40degC).

Single Use

Plastic Point Protectors are single use only.

Precautions

See Through Point Protectors should be stored away from any sources of ignition.

Disposal

Unused See Through Point Protectors are manufactured from Polypropylene and are recyclable in the plastic waste recycling stream. They are also suitable for disposal to landfill.

Used See Through Point Protectors should be disposed of by considering possible microbial contamination. Users should reference their own procedures in this matter.

How to use See Through Point Protectors:

Follow the steps below:

