Instrument Tip Caps – TPE

EN - INTENDED USE:

Key Surgical® Instrument Tip Caps are designed to protect the tips of surgical instruments when not in use or during instrument sterilization.

PRECAUTION:

Does not exhibit a sharp melting point but may soften gradually during sterilization. Avoid letting caps touch or using caps in a way that will exert pressure on the cap during sterilization as this may cause them to melt.

Should not be altered and are not intended for implantation.

Must be removed from the instrument prior to use.

Single use only; reuse may lead to patient cross contamination.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

MATERIAL:

Evoprene™ G970 (TPE) material that is free of BPA and DEHP. Not made with natural rubber latex. The colorant added varies only in color, not in the basic chemical composition.

INSTRUCTIONS FOR USE:

1. Does not need to be cleaned or disinfected before use. If desired, the caps can be put through a washer or alcohol bath.
2. Place the cap on the tip(s) of an appropriately sized clean instrument before sterilization and storage. Hinged instruments should be sterilized in an open position.
3. Package instruments in sterile packaging device (e.g., pouch, wrapped tray, closed container).
4. Place package in sterilizer. The instrument surface under the cap will be sterile, regardless of use of vented or non-vented caps.
5. Single use and can be sterilized only once.

STERILIZATION:

Validated for sterilization efficacy according to the following parameters:

Steam Sterilization			
Cycle Type	Temperature	Minimum Time	Minimum Validated Dry Time
Pre-Vacuum	134° C (273° F)	3 minutes	20 minutes
Pre-Vacuum	121° C (250° F)	15 minutes	20 minutes

STORAGE:

Store according to the information on the device labelling and/or this IFU.

DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred

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KEY SURGICAL LLC
MDSS UK RP LTD