Cannulas – Non-Sterile

EN - INTENDED USE:

Key Surgical® Cannulas are a non-sterile device intended to be used for topical wound irrigation. Cannulas are not intended for direct contact with wounds.

PRECAUTIONS:

Single use only; reuse may lead to patient cross-contamination.

Visually inspect the cannulas before use to ensure that they are free from visible contaminants.

CONTRAINDICATION:

No known contraindications and/or adverse effects.

INSTRUCTIONS FOR USE:

1. Inspect for usability and surface integrity.
2. Clean and sterilize before use.
3. This product is single use. Upon completion of procedure, dispose of used cannula according to facility protocol.

CLEANING AND DISINFECTION:

Manual or mechanical processing with or without ultrasonic treatment.

Chemicals and temperatures to be used for cleaning and disinfecting: acidic / neutral / alkaline with or without surfactant additives, chemically at max. 60° C / 140° F or with fully demineralized, thermally, at max. 93° C / 199° F.

It is assumed that commercially available products approved for the intended application are used and that the user complies with recommended concentrations, exposure times, and temperatures. It must be guaranteed also that no residues remain on the products. Fully demineralized water must be used for the final rinse.

Drying at max. 120° C / 248° F.

STERILIZATION:

Validated for sterilization according to the following parameters:

Steam Sterilization			
Cycle Type	Temperature	Minimum Time –Full Cycle	Minimum Dry Time
Pre-Vacuum	121° C (250° F)	15 minutes	20 minutes
Pre-Vacuum	134° C (273° F)	3.5 minutes	20 minutes

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

STORAGE REQUIREMENTS:

Store according to the information on the device labelling and/or this IFU.

DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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MDSS UK RP LTD