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**Hersteller: Key Surgical GmbH**

**Produktname:**

Silikongitterabdeckung

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**Product name**

**Silicone Mesh Cover**

*Please observe the reprocessing instructions in the following table. These instructions refer to cleaning, disinfection and sterilization. The manifold reprocessing possibilities available derive from the material compatibility of this product. The operator bears sole responsibility for assuring a successful reprocessing outcome.*

<b>Preparation at location of use</b>	
dry / wet Not applicable (single use product)	Wet means the products are immersed in a non-fixating active cleaning disinfection solution immediately after the operation. Dry means the products are transferred to CSSD without any disinfectant or other additional fluids immediately after the operation.
<b>Cleaning and Disinfection</b>	
Manual or automated reprocessing	
<b>Chemicals and Temperatures for Cleaning and Disinfection</b>	
Acidic / neutral / alkaline with/without added surfactants, chemical at max 60 °C or with deionised water thermally at max. 93 °C	It is assumed that commercially available products that have been approved for cleaning and disinfection will be used. Likewise, it is assumed that the recommended concentrations, exposure times and temperatures will be observed. It must be ensured that no residues remain on the products. Demineralized water must be used for the final rinse.
<b>Drying</b>	
Max. 100 °C	
<b>Packaging</b>	
Packaging materials in accordance with the EN 868 and ISO 11607 series of standards, which are approved by the manufacturer for the specified sterilization process	
<b>Sterilization</b>	
One or more procedures to choose from: Validated steam sterilization process in the 134 °C program or validated steam sterilization process in the 121 °C program.	Sterilisation and hold times are subject to national regulations and guidelines and, as such, cannot be stipulated in general terms. The operator bears responsibility for ensuring that reprocessing and sterilisation, as actually conducted with the equipment, materials and personnel used for reprocessing and sterilisation, achieve the required results. Validation and routine monitoring of the process are needed to that effect.
<b>Alternative sterilisation processes*</b>	*Note on alternative sterilisation processes: Steam sterilisation has become established worldwide as a very safe and reliable sterilisation method and is thus the method of choice for temperature- and humidity-resistant sterile supplies. Here reference is generally made to steam sterilisation using a validated steam sterilisation process (see also DIN EN ISO 17665). There is therefore no need to use alternative sterilisation processes, e.g. formaldehyde or ethylene oxide, for steam-tolerant products. However, each operator of a sterilisation unit is free to have sterilisation validation conducted using an alternative sterilisation process.
<b>Storage</b>	
No special requirements	Please note the general principles and requirements in dealing with sterile goods and sterile packaging.
<b>Other important notes</b>	
<i>The preparation instructions given here cannot replace detailed process descriptions, as the large number of preparation processes used worldwide cannot be described in detail. All information is without guarantee. The aluminium trays are not marketed as medical products in terms of their intended use. The area of application and the consequences arising from this are the sole responsibility of the operator.</i>	
<b>Contact:</b>	Key Surgical GmbH, Zum Windpark 1, 23738 Lensahn Tel. 04363/905900, Fax 04363/90590590

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**Manufacturer: Key Surgical GmbH**