Hot Air Film

EN - INTENDED USE:

Key Surgical® Hot Air Film is used as sterilization packaging for dry heat sterilization.

CAUTION:

Cannot be sealed with impulse or continuous heat-sealing machines. Alternative sealing methods are not known. Dry heat tape is recommended for sealing.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

INSTRUCTIONS FOR USE:

Pouch:

1. Select the appropriate pouch size.
2. Place item(s) into the pouch.
3. Fold the open side of the pouch 4-5 times.
4. The folded side can be fixed with a piece of adhesive tape so that the side does not roll up again.
5. Stick an adhesive tape suitable for hot-air sterilization 2 times around the folded side.

Tube:

1. Cut the tubular film to the desired size.
2. Fold one side of the tube 4-5 times.
3. The folded side can be fixed with a piece of adhesive tape so that the side does not roll up again.
4. Stick an adhesive tape suitable for hot-air sterilization 2 times around the folded side.
5. Place item(s) into the pouch.
6. Seal the second side following the same procedure as before.

STERILIZATION:

Validated for sterilization according to the following parameters:

	Temperature	Minimum Time – Full Cycle
Cycle Type (Dry Heat)	140° C (284° F)	60 minutes
	160° C (320° F)	45 minutes
	180° C (356° F)	30 minutes

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

STORAGE:

Store in a cool, dry place.

DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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MDSS UK RP LTD