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with Process challenge device (PCD)

Type 2

Chemical Indicator for Checking Steam Sterilization by User

The *PVP Steam-Control Indicator* changes colour, as a function of the time, from different shades of light-brown through grey to black on being exposed to hot steam between 121° C (15 min) and 134°C (3.5 min).

Pursuant to DIN EN ISO 11140-1, the indicator thus complies with the requirements addressed to indicators for Type 5.

It is recommended that indicators be used when conducting validation or microbiological investigation of the sterilizer. The black colour tone obtained can be used as a typical colour comparison for the respective procedure (121 °C and 134 °C).

The indicators are available as a pad comprising 25 perforated indicator strips. Individual strips can be easily torn off at the perforated narrow end before use and then enclosed with the items to be sterilized.

If the indicator is to be used in a process challenge device the strips have to be folded at the seam so that the laminated indicator side lies on top of itself.

In this combination (indicator strip and PCD) this product is applicable for class B sterilizers with pulsed flow vacuum process according to DIN EN 13060 and large sterilizers according to DIN EN 285 for monitoring steam sterilization for hollow instruments (A).

At the end of sterilization the indicator strips will have assumed a uniform black colour similar to the reference colour already obtained.

If the indicator has changed colour in the manner required, it can be assumed that the prerequisites for orderly and effective sterilization have been assured.

If this is not the case, the sterilization process should be critically reviewed bearing in mind all the instructions and guidelines given by the manufacturer.

Microbiological investigation of the sterilizer may be needed.

The indicator strip is slotted at the front side so that by removing the protective layer, the indicator strip can be affixed to documentation.


Please note absolutely:

To prevent as far as possible any effect that light or other environmental influences might have on the indicator, the latter must be removed from the packaging only immediately prior to use.

The batch number is specified on the rear of each indicator pad, outside the perforated indicator strip. This must always be quoted for product identification in the event of any queries.

Process challenge device (PCD)

1. Unscrew the cap of the PCD.
2. Fold the chemical indicator strip in such a way that the indicator is on the inside.


3. Insert the folded indicator strip with the closed end first in the cap of the PCD.

4. Insert the indicator strip completely.

5. Close the PCD by screwing the cap firmly.

6. Place the PCD close to the bottom of the sterilizer and sterilize it together with the load.

Indicator evaluation

7. Unscrew the PCD and take out the chemical indicator strip.

Colour Change of Steam Control Indicators

before sterilization:

after orderly sterilization:

all indicating fields are dark brown or black

failed colour change:

the fields of the indicator strip have different colours (in case of use in challenge device)

possible reason:

inadequate air removal

the fields of the indicator strip have an evenly green-brown colour.

possible reason:

insufficient saturated steam

Documentation sheet sterilization: **STEAM** 121°C / 134°C

template

autoclave:

date	Batch number	program	Indicator	Process indicator correct color change (yes/no)	Release batch (name)