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Vertrieb: Key Surgical GmbH

Produktname:

Silikonringe

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Product name:

Silicone ring

Please observe the processing instructions given in the following table, which relate to cleaning, disinfection and sterilization. While there are many processing options available, the correct choice depends on the material of the item at hand (i.e. its compatibility with the selected procedure). Note that the operator is solely responsible for the success of the selected procedure.	
Preparation at location of use	
Recommendation: The silicone rings must be replaced each time before use	
Cleaning and disinfection	
manual or mechanical reprocessing without ultrasound treatment	
Cleaning and disinfecting chemicals and temperatures	
acidic / neutral / alkaline >pH 10 with or without tenside additives, chemically at max. 60° C / 140° F or with VE water, thermally, at max. 93° C / 199° F	It is assumed that commercially available products approved for the intended application are used and that the user complies with recommended concentrations, exposure times and temperatures. It must be guaranteed also that no residues remain on the products. Fully demineralized water must be used for the final rinse.
Drying	
drying at max. 100° C / 212° F	
Checks, maintenance/care and inspection	
Inspection for usability and surface integrity	
Sterilization	
One or several procedures available: Validated steam sterilization process 134°C (273 °F) / 3.5 min, 5 min respectively 18 min, or validated steam sterilization process 121°C (250°F) / 15 min.	Sterilization and holding times are subject to national regulations and guidelines and therefore cannot be set universally. It is the operator's responsibility to ensure that the intended results are actually achieved with the processing and sterilization methods, equipment, materials and personnel used. This, in turn, requires validation and routine monitoring of all relevant processes.
Other sterilization processes	
Alternative sterilization methods*	*Notes on alternative sterilization methods: Steam sterilization has established itself worldwide as a very safe and highly reliable sterilization method and, therefore, represents the method of choice with regard to temperature- and humidity-insensitive goods. As a rule, steam sterilization in the form of a validated process is required (see also DIN EN ISO 17665-1). Therefore, it is not necessary to use alternative sterilization methods – such as low-temperature, formaldehyde or ethylene oxide for steam-sterilizable medical devices. Nonetheless, sterilizer operators are free to validate sterilization processes using alternative methods for sterilizing their products.
Storage	
Requirements:	Stored at room temperature. Do not expose to direct sunlight.
Other important notes	
<i>The processing information given herein cannot replace detailed process descriptions because it is not possible to give a detailed account of the numerous processing methods used all over the world. All information is given without guarantee.</i>	
Contact:	Key Surgical GmbH, Zum Windpark 1, 23738 Lensahn Tel. 04363/905900, Fax 04363/90590590

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Distributor: Key Surgical GmbH