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**Hersteller: Key Surgical GmbH**

**Produktname:**

Silikon-Noppenmatten

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**Product name:**

Silicone Nub Mat

Please observe the processing instructions given in the following table, which relate to cleaning, disinfection and sterilization. While there are many processing options available, the correct choice depends on the material of the item at hand (i.e. its compatibility with the selected procedure). Note that the operator is solely responsible for the success of the selected procedure.	
<b>Preparation at location of use</b>	
wet or dry	Wet means the products are immersed in a non-fixating active cleaning disinfection solution immediately after the operation. Dry means the products are transferred to CSSD without any disinfectant or other additional fluids immediately after the operation.
<b>Cleaning and disinfection</b>	
manual or mechanical reprocessing without ultrasound treatment	
<b>Cleaning and disinfecting chemicals and temperatures</b>	
acidic / neutral / alkaline >pH 10 with or without tenside additives, chemically at max. 60° C / 140° F or with VE water, thermally, at max. 93° C / 199° F	It is assumed that commercially available products approved for the intended application are used and that the user complies with recommended concentrations, exposure times and temperatures. It must be guaranteed also that no residues remain on the products. Fully demineralized water must be used for the final rinse.
<b>Drying</b>	
drying at max. 100° C / 212° F	
<b>Checks, maintenance/care and inspection</b>	
Inspection for usability and surface integrity	
<b>Sterilization</b>	
One or several procedures available: Validated steam sterilization process 134°C (273 °F) / 3.5 min, 5 min respectively 18 min, or validated steam sterilization process 121°C (250°F) / 15 min.	Sterilization and holding times are subject to national regulations and guidelines and therefore cannot be set universally. It is the operator's responsibility to ensure that the intended results are actually achieved with the processing and sterilization methods, equipment, materials and personnel used. This, in turn, requires validation and routine monitoring of all relevant processes.
<b>Other sterilization processes</b>	
Alternative sterilization methods*	*Notes on alternative sterilization methods: Steam sterilization has established itself worldwide as a very safe and highly reliable sterilization method and, therefore, represents the method of choice with regard to temperature- and humidity-insensitive goods. As a rule, steam sterilization in the form of a validated process is required (see also DIN EN ISO 17665-1). Therefore, it is not necessary to use alternative sterilization methods – such as low-temperature, formaldehyde or ethylene oxide for steam-sterilizable medical devices. Nonetheless, sterilizer operators are free to validate sterilization processes using alternative methods for sterilizing their products.
<b>Storage</b>	
Requirements:	Stored at room temperature. Do not expose to direct sunlight.
<b>Other important notes</b>	
<i>The processing information given herein cannot replace detailed process descriptions because it is not possible to give a detailed account of the numerous processing methods used all over the world. All information is given without guarantee.</i>	
<b>Contact:</b>	Key Surgical GmbH, Zum Windpark 1, 23738 Lensahn Tel. 04363/905900, Fax 04363/90590590

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**Manufacturer: Key Surgical GmbH**