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### Interlock BD*pro* **STEAM**

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#### Content

- 1 PCD with stainless steel holder
- 250 indicators in a plastic sliding box
- 1 O-ring
- Screw with O-ring for leakage test

#### Advantages

- Easy handling
- Robust construction
- Logical-constructive structure following hollow instruments
- Processing of leakage tests
- Simple O-ring replacement
- Long-lasting operating times
- Reliable results

#### Pre-Caution

- To eliminate unnecessary environmental impact on the indicator, take the indicator out of its packaging just before application. Do not store the indicators close to heat sources, e.g. on the steriliser

#### General notes:

- The indicator strip and the PCD form the indicator system.
- A visual inspection on the integrity and the fit of the O-ring has to be performed, before using the PCD.
- Involve the indicator system in the assessment of performance, for validation and routine controls of sterilisation processes for medical devices with saturated steam.
- The leakage test needs to be undertaken every 250 cycles. Therefore, place the screw (M3) with the O-ring vertically into the steam inlet opening at the bottom of the PCD as shown on figure 3. Then replace the O-ring of the indicator mount and place an indicator strip into it. Run a sterilisation process. The result should be no colour change of the indicator. If the indicator colour changes, repeat the procedure. If the colour of the indicator strip changes again, contact Interlock.
- It is recommended to use heat protection gloves when taking the PCD out of the sterilizer


#### Application:

The Interlock BD*pro* is an air removal system for steam sterilisers according to the standard DIN EN 285 in the healthcare. That means that air gets removed by fractional pre-vacuum steps. According to DIN EN ISO 11140-4 the function of the test system is the complete air removal and therewith the steam penetrations as an alternative to the daily execution of the required BD test. The indicator's colour change depends on the steam impact. If enough saturated steam condensates, the colour change to black takes place within this temperature-time-combination of 134°C and 3.5 min or 121°C and 15 min.

**Attention: Before first use, the screw must be removed!**

#### Equipping of the PCD to the indicator system

Screw off the indicator mounting capsule. Bend the indicator strip at the fold so the printed sides are facing each other. Slide it into the gap like shown in the figure below. Screw the PCD together and place it close to the door in the lower part of the chamber.

When the program has finished, remove the indicator from the PCD and evaluate the colour change. The indicator is self-adhesive and can be paste into the documentation. For this purpose remove the backing layer at the backside perforation.

**Information on indicator evaluation and possible sources of errors**

The indicator changes colour gradually from cream to grey to black when being affected by saturated steam depending on time.

| Appearance of the indicator   | Notes on evaluation   |
|---|---|
| <br>cream            | Initial state of the indicator  |
| <br>black            | Proper test result, indicated value of colour density = 1.2   |
| <br>colour gradation | The gradient can occur in different colour intensities. It occurs when the indicator was influenced by no or less steam, i.e. the PCD is not completely free from inert gas, e.g. air. Different causes are possible: e.g. leakage, lack of air removal, insufficient steam quality, overheating. |
| <br>grey             | The air was removed completely from the PCD but lack of steam condensation on the indicator. Leads to an insufficient colour depth and indicates overheated steam or short exposure time.   |
| <br>cream            | No steam has penetrated the indicator. Possible cause: No or less air removal, too big amount of inert gases.   |

**Figure 3**

**Storage conditions:** see product label

**PCD Holder**

The holder can be mounted on the PCD by slightly extending both wire ends and pushing the holder over the PCD

Catalog number

Serial number

Batch code

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Humidity limit, limitation

Temperature limitation

Keep away from direct sunlight

Keep dry