Stainless Steel Baskets

EN - INTENDED USE:

Key Surgical® Stainless Steel Baskets (tray) are designed to contain and protect surgical instruments during the sterilization process.

PRECAUTIONS:

Do not overload.
Overloading may cause damage.
Visually inspect prior to use. Check for loose wires or burrs that may puncture, or tear sterilization wrap or pouch.
Inspect handles to ensure they are secure prior to use and transport.
Avoid prolonged exposure to saline to minimize the chance of corrosion.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

INSTRUCTIONS FOR USE:

1. Reusable and should be decontaminated, cleaned, and inspected prior to each use according to the instructions below.
2. Place unpouched instruments in tray. An approved liner, mat, or other positioning device can be used to accommodate the instruments to be sterilized.
3. Wrap with an approved non-woven textile wrap or place in any other FDA cleared packaging system (i.e. sterilization pouch or container). Do not wrap more than one tray per wrap.
4. The tray should always be placed flat on the sterilizer rack prior to sterilization. Stacking of trays inside the sterilization chamber is not recommended.
5. Sterilize the loaded wrapped tray in accordance with the instrument device manufacturer's IFU.
6. After sterilization, move wrapped tray to a sterile storage area for storage or transport to the point of use location.
7. Clean as soon as possible after use. If cleaning must be delayed, immerse in an enzymatic solution or water to prevent drying and encrustation of surgical soil. Remove excessive soil with a disposable wipe.

MANUAL CLEANING:

1. Pre-rinse under cold tap water for one (1) minute to remove gross debris
2. Soak for a minimum of two (2) minutes in a pH neutral detergent, prepared in accordance with the manufacturer's IFU.
3. Rinse under cold tap water for one (1) minute.
4. Ultrasonically clean for a minimum of five (5) minutes in a neutral pH detergent, prepared in accordance with the manufacturer's IFU.
5. Rinse under cold tap water for one (1) minute.

AUTOMATED CLEANING:

It may be necessary to manually clean prior to automated processing to improve the removal of adherent soil. Follow the previous instructions for manual cleaning.

1. Run the automatic wash cycle – minimum cycle parameters:
 - 1 minute cold pre-rinse
 - 5 minute enzyme wash at 43° C (109° F) minimum temperature
 - 1 minute cold rinse
 - 7 minute dry at 90° C (194° F) minimum temperature

CLEANING INSPECTION:

Visually inspect before sterilization or storage to ensure the complete removal of soil from surfaces. If soil is still present, re-clean.

STERILIZATION:

Validated for sterilization according to the following parameters:

Steam Sterilization		
Cycle Type	Temperature	Minimum Time – Full Cycle
Gravity	121° C (250° F)	30 minutes
Pre-Vacuum	132° C (270° F)	4 minutes
Pre-Vacuum	134° C (275° F)	3 minutes
Immediate Use Pre-Vacuum	132° C (270° F)	3 minutes

100% Ethylene Oxide (EtO)	
Concentration	725 mg/L
Temperature	55° C

Exposure Time	60 minutes
Humidity	50-80%

STERRAD® System and Cycle	
STERRAD® System and Cycle within the United States	STERRAD® System and Cycle outside United States
STERRAD® 50	STERRAD® 50
STERRAD® 200	STERRAD® 200 Short cycle
STERRAD® 100S	STERRAD® 100S Short cycle
STERRAD® NX Advanced cycle	STERRAD® NX Advanced cycle
STERRAD® NX Standard cycle	STERRAD® NX Standard cycle
STERRAD® 100NX Standard cycle	STERRAD® 100NX Standard cycle

STERIZONE® System	STERIZONE® VP4
-------------------	----------------

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL:

Requires replacing if there are any signs of deterioration or loss of functionality. Dispose according to facility protocol.

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