Transportation Wrap

EN - INTENDED USE:

Key Surgical® Clinipak Transportation Wrap and Transportation Wrap Dry is intended for use as an outer wrap that must be used in combination with one or more inner barrier sterilisation wraps.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

PRECAUTIONS:

Store away from any sources of ignition.

Single use only; reuse may lead to patient cross contamination.

Transportation Wrap Dry must not be used in H2O2/Plasma sterilisation processes.

INSTRUCTIONS FOR USE:

Sequential Wrapping: Envelope Fold/Square Fold

This wrapping technique allows for sequential opening when items are unwrapped after sterilisation.

1. Select the appropriate size and material weight for inner and outer layer.
2. Place item(s) into a basket/tray on the inner wrap.
3. Wrap using aseptic wrapping technique.
4. Place the wrapped item on the outer layer
5. Wrap using aseptic wrapping technique.
6. Seal the wrap in a closed position using a suitable closure method.
7. Label the package.

Simultaneous Double Wrapping: Envelope Fold/Square Fold

This wrapping technique allows items to be opened when items are unwrapped after sterilisation.

1. Select the appropriate size and material weight for inner and outer layer.
2. Align the inner wrap and outer layer, with the inner wrap on top of the outer wrap.
3. Place item(s) into a basket/tray then place the item in the middle of the wrap and wrap both sheets together.
4. Seal the wrap in a closed position using a suitable closure method.
5. Label the package.
6. Wrap using aseptic wrapping technique.
7. Seal the wrap in a closed position using a suitable closure method.
8. Label the package.

STERILISATION:

Compatible with steam and ETO sterilisation processes.

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

Unused: Can be disposed of in the plastic recycling stream.

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MANUFACTURER: KEY SURGICAL LLC