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REF 08500

Table of Contents:

Content of test-equipment according to DIN 96298-3	13
Preparation	13
Testing the cutting efficiency of bone splitters, hollow chisel forceps and bone punches	13
through-cutting conchotomes	14
scalpels and knives	15
scissors	15
chisels, raspatories, curettes and sharp spoons	16
Testing the elasticity of ring-nosed forceps	16
Testing the functionality of the profile of ring-nosed forceps with atraumatic perforation	17
Testing the leak tightness of ring-nosed forceps with atraumatic perforation	17
Testing the elasticity of needle holders	18
Testing the functionality of the profile of needle holders	18
Testing the gripping strength of non-cutting conchotomes	19
Testing the lenses of rigid endoscopes	19
List of orders	20

Test-equipment according to DIN 96298-3

Instruction for use

Test-equipment

- for bone splitters, hollow chisel forceps and bone punches: cardboard 160 g/m²
- for cutting conchotome: synthetic film with a strength of 100 µm
- for scalpels and knives: synthetic film with a strength of 100 µm / test drum made of synthetic material
- for surgical scissors: red latex strap, styrene-butadiene rubber
- for chisels, raspatories, bone curettes and sharp spoons: plastic round profile
- for ring-nosed forceps: plates made of synthetic material 1,0 / 1,5 / 2,0 / 2,5 / 3,0 / 4,0 mm, paper 30 g/m², polybag 50 µm
- for needle holder: wire 0,4 / 0,6 / 0,8 / 1,0 mm, paper 30 g/m², magnifying glass
- for non-cutting conchotome: synthetic film with a strength of 300 µm
- for rigid endoscopes: monocular magnifying glass

Additional test-equipment (not according to DIN 96298-3):

- for rigid endoscopes: magnifying glass for testing the inner lenses

1.0 Preparation

Clean all instruments before tests.

2.1 Testing the cutting efficiency of bone splitters, hollow chisel forceps and bone punches

Test material: cardboard 160 g/m² (refill pack REF 95172)


- a) Bone splitters and hollow chisel forceps:
To test the cutting efficiency, make three cuts through the cutting material. The test cutting length is at least 1/3 of the cutting length in the anterior part of the working end. It must be possible to separate the test material smoothly, and without it slipping.
- b) Bone punches:
To test the cutting efficiency, make three cuts through the cutting material. The test cutting length is at least 2/3 of the cutting length in the anterior part of the working end. It must be possible to separate the test material smoothly, and without it slipping.

Bone splitters

Bone punches

Hollow chisel forceps

2.2 Cutting efficiency of through-cutting conchotomes

Test material: synthetic foil with a strength of 100 μm (refill pack REF 95141)

Test: To test the cutting efficiency, make three cuts through the cutting material. It must be possible to separate the test material smoothly, and without it slipping.


Through-cutting conchotomes

2.3 Cutting efficiency of scalpels and knives

Test material: plastic foil with a strength of 100 µm (refill pack 95141) / test drum (refill pack 95140)

Test: Scalpels and knives, with the exception of sharp eye instruments, are tested by making a cut through the foil, while guiding the cut across the anterior half of the blade. The foil must be separated with a smooth cut.

Sharp eye instruments such as e.g. cataract knives and lancets are tested by making a cut through the foil stretched across the test drum.

Knives

Sharp eye instruments

2.4 Cutting efficiency of scissors

Test material: red latex tape (refill pack REF 08502), styrene-butadiene rubber NR-SBR (refill pack REF 08534)

Test: Use the test material listed in table 1. To test the cutting efficiency, make three continuous cuts across 2/3 of the cutting length of the scissors through the test material, without exerting any lateral pressure. It must be possible to separate the test material smoothly, and without it slipping.

FM	Strength or layers and type of test material
	Styrene-butadiene rubber and latex tape
≥ 0,70	3,00 mm NR-SBR
0,71 - 1,00	2,00 mm NR-SBR
1,01 - 2,00	1,00 mm NR-SBR
micro or spring scissors	0,2 mm latex tape red

Table 1

2.5 Cutting efficiency of chisels, raspatories, curettes and sharp spoons

Test material: plastic round profiles with a diameter of 6 mm (refill pack 95169)

Test: The instrument is pulled or pushed over the test material at an angle appropriate to the application with as little force as possible. The test is considered to have been passed if material is removed in the process.

Curettes


3.0 Testing the elasticity - Ring-nosed forceps

Test accessories: plastic plates 1,0 / 1,5 / 2,0 / 2,5 / 3,0 / 4,0 mm (refill pack REF 95154)

Test: Place a piece of test material listed in table 2 between the functional surfaces in the anterior third of the working end. The ring-nosed forceps must be fully engaged (all ratchet teeth engaged) and stored at room temperature 3h +/- 5 minutes after being subjected to this stress. On completion of the test, there must not be any deformations or lasting changes noticeable.

Total length	Length of the mouth piece	Strength of the test material
< 100 mm	each	1,0 mm
100 mm to < 160 mm	< 25 mm	1,5 mm
100 mm to < 160 mm	25 mm to 30 mm	1,5 mm
100 mm to < 160 mm	> 30 mm	2,0 mm
160 mm to 220 mm	< 25 mm	2,0 mm
160 mm to 220 mm	25 mm to 30 mm	2,5 mm
160 mm to 220 mm	> 30 mm	3,0 mm
> 220 mm	< 25 mm	2,5 mm
> 220 mm	25 mm to 40 mm	3,0 mm
> 220 mm	> 40 mm	4,0 mm

Table 2


Ring-nosed forceps

3.1 Testing the fillet of the tooth profil - Ring-nosed forceps with atraumatic perforation

Test material: paper 30 g/m² (refill pack REF 95143)

Test: The test material is placed between the functional surfaces of the clamp and then the clamp is fully closed. After holding for at least 2 seconds, the clamp is opened again and the test material withdrawn. The impression of the longitudinal profile must be uniformly visible and must not be perforated in any place.

Atraumatic clamp

3.2 Testing leak tightness - Ring-nosed forceps with atraumatic perforation

Test material: polybag 50 µm (refill pack REF 09912)

Test: The polybag is filled half full with water and sealed. Then the polybag is gripped with the clamp jaw on the top diagonally to the corners. The clamp lock must be fully closed at this point. The bag is then turned over and no water is allowed to flow into the clamped corner. Even after removing the clamp, no water may flow out of the bag.

4.0 Testing the elasticity – Needle holder

Test material: wire 0,4 / 0,6 / 0,8 / 1,0 mm (refill pack REF 95155), hand magnifying glass (5x zoom) (REF 1904180)

Test: Place a piece of wire as listed in Table 3 between the functional surfaces at the tip of the needle holder. Then the needle holder must be fully engaged (all ratchet teeth engaged) and stored at room temperature (23°C +/- 2°C) 3h +/- 5 min. after being subjected to this stress. On completion of the test, there must not be any deformations or lasting changes noticeable. No tears should be noticeable during visual inspection with a magnifying glass.

Total length [mm]	Barrel-Surfaces profile	Diameter test material [mm]
bis 160	flat or cross groove, division 0,2	0,4
	flat or cross groove division 0,4 bis 0,6	0,6
über 160	flat or cross groove division 0,2	0,8
	cross groove division 0,4 bis 0,6	1

Table 3

Needle holder

4.1 Testing the functionality of the profile – Needle holder

Test material: Paper 30 g/m² (refill pack REF 95143)

Test: The paper is placed between the functional surfaces of the needle holder and the latter is closed to the last notch. The needle holder is then opened again and the test material must be evaluated visually. Two results of such a test are shown in Figure 1 and 2. Figure 1 shows that the position of the functional surfaces is not in order, as the front half of the impression is not sufficiently imaged. Figure 2, on the other hand, shows the impression over the entire surface and thus the position of the functional surfaces is in order.


Figure 1

Figure 2

5.0 Testing the gripping strength - non-through-cutting conchotomes

Test material: synthetic foil with a strength of 300 μm (refill pack REF 95139)

Test: To test the gripping strength, the test material must continue to be held by the jaws even using strong counterforce.


non-through-cutting conchotomes

6.0 Testing the lenses of rigid endoscopes

Test material: (monocular) magnifying glass (REF 08509)

Test: Place the (monocular) magnifying glass before one eye and direct the endoscope towards a bright, reflecting surface and look with the magnifying glass into the ocular end of the endoscope. The optical axes (of the eye using the magnifying glass as well as of the longitudinal axis of the endoscope) must be perfectly aligned. While varying the distance between the endoscope and the eye (with the axis orientation unchanged), look through the individual lenses within the endoscope into the optics. By focusing on a lens surface contaminants can be detected or ruled out. By focusing on a lens, lens breakages can be seen.

7.0 Refill package for test-equipment according to DIN 96298-3

SKU	Description	UOM
95172	cardboard, 160 g/m ²	50 pieces
95141	synthetic foil with a strength of 100 µm	100 pieces
95140	test drum	1 piece
08502	red latex tape, 100% latex, powdered with talc	1 tape
08534	styrene-butadiene rubber NR-SBR	1 piece
95169	round profile with a diameter of 6 mm	1 piece
95154	plates made of synthetic material 1,0 / 1,5 / 2,0 / 2,5 / 3,0 / 4,0 mm	each 3 pieces
95143	paper pack, 30 g/m ²	1 pack
95155	wire 0,4 / 0,6 / 0,8 / 1,0 mm	each 10 pieces
95139	synthetic foil with a strength of 300 µm	50 pieces
1904180	hand magnifying glass, 4x zoom with led	1 piece
08509	monocular magnifying glass, 10x zoom	1 piece
09912	polybag 50 µm	20 pieces