Tray Liners and Soaker Sheets

EN - INTENDED USE:

Key Surgical® Tray Liners and Soaker Sheets are intended for placement on the bottom of a sterilisation tray to absorb wetness and humidity during processing.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

PRECAUTIONS:

Store away from a source of ignition.
Single use only; reuse may lead to patient cross contamination.

PREPARATIONS FOR USE:

- Select appropriate size product for tray
Note: Tray Liners may be cut to the required size to fit the selected tray
- Product may be placed inside and/or under the tray
- Place product inside the wrapped tray or rigid container system

TLP-xxx-xxx and TL-xxx-xxx products can be used in the following configurations:

- Pre-Vacuum Steam/100% Ethylene Oxide (EtO)
- Up to 2 tray liners inside and 2 tray liners outside
 - Up to 8 folded layers (4 tray liners)

INSTRUCTIONS FOR USE:

- Sterilise according to the information below
- Place product inside or under tray
- Place instruments in tray and wrap or close container according to manufacturer's instructions

STERILISATION:

Provided non-sterile and are single use.
Not suitable for vaporized hydrogen peroxide/plasma sterilisation.
Validated for sterilisation efficacy using the following cycle parameters:

Steam Sterilisation				
Item Number	Cycle	Temperature	Minimum Cycle	
TLP-012-018 TL-202-500	Pre-Vacuum	132°C (270°F)	3 minutes	
TLP-020-025 TL-203-500		132°C (270°F)	4 minutes	
TL-101-200 TL-204-500		134°C (273°F)	3 minutes	
TL-101-200 TL-202-500	Gravity	121°C (250°F)	30 minutes	
TL-203-500 TL-204-500				
100% Ethylene Oxide (EtO)				
Item Number	Concentration	Temperature	Exposure Time	Humidity
TLP-012-018 TLP-020-025 TL-101-200 TL-105-250 TL-202-500 TL-203-500 TL-204-500	735 mg/L	55° C	60 minutes	50-70%

Items 1136011 and 1137011 are compatible for steam sterilization up to 134°C for 3.5 minutes or 121°C for 15 minutes.
All other products are compatible for steam sterilization up to 134°C for up to 18 minutes.

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

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MANUFACTURER: KEY SURGICAL LLC