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Key Surgical GmbH

Produktname:

MagnoMat (Magnetische OP-Matte)

REF 26300

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Processing Instruction

Produktname: **Key Surgical GmbH**
MagnoMat (Magnetic OR-Mat)
 REF 26300

Please observe the reprocessing instructions in the following table. These instructions refer to cleaning, disinfection and sterilisation. The manifold reprocessing possibilities available derive from the material compatibility of this product. The operator bears sole responsibility for assuring a successful reprocessing outcome.

Preparation at the site of use	
wet or dry	Wet means the instruments are immersed in a non-fixating active cleaning disinfection solution immediately after the operation. Please observe the manufacturers instructions. Dry means the instruments are transferred to CSSD without any disinfectant or other additional fluids immediately after the operation.
Cleaning and Disinfection	
Manual or automated reprocessing without ultrasonic treatment	
Chemicals and Temperatures for Cleaning and Disinfection	
Acidic / neutral / alkaline with/without addition of surfactants, chemically at max. 60 °C or thermally with demineralised water at max. 93 °C	It is assumed that commercially available products that have been approved for cleaning and disinfection will be used. Likewise, it is assumed that the recommended concentrations, exposure times and temperatures will be observed. It must be ensured that no residues remain on the products. Demineralised water must be used for the final rinse.
Drying	
Max. 100 °C	
Inspection, Maintenance and Testing	
Visual inspection for suitability for use.	
Sterilisation	
A choice of one or several processes: Validated steam sterilisation process with 134 °C program or validated steam sterilisation process with 121°C program.	Sterilisation and hold times are subject to national regulations and guidelines and, as such, cannot be stipulated in general terms. The operator bears responsibility for ensuring that reprocessing and sterilisation, as actually conducted with the equipment, materials and personnel used for reprocessing and sterilisation, achieve the required results. Validation and routine monitoring of the process are needed to that effect.
Alternative sterilisation processes*	*Note on alternative sterilisation processes: Steam sterilisation has become established worldwide as a very safe and reliable sterilisation method and is thus the method of choice for temperature- and humidity-resistant sterile supplies. Here reference is generally made to steam sterilisation using a validated steam sterilisation process (see also DIN EN ISO 17665). There is therefore no need to use alternative sterilisation processes, e.g. formaldehyde or ethylene oxide, for steam-tolerant products. However, each operator of a sterilisation unit is free to have sterilisation validation conducted using an alternative sterilisation process. The sterilisation processes H2O2 and dry heat must not be used.
Storage	
Requirements:	The MagnoMat must be stored at room temperature. Do not expose to direct sunlight. Store in dry conditions.
Further Instructions	
<i>The reprocessing instructions given here are not a substitute for a detailed process description since it is not possible to elaborate in detail on the myriad reprocessing methods used worldwide. All information without guarantee.</i>	
Contact details:	Key Surgical GmbH, Zum Windpark 1, D-23738 Lensahn Tel. +49 4363/905900, Fax +49 4363/90590590

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Key Surgical GmbH, Zum Windpark 1, 23738 Lensahn