Spreaders made of Stainless Steel

EN - INTENDED USE:

Key Surgical Spreaders made of stainless steel are designed to hold non-disassemble clamps open during the cleaning and sterilization process.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

INSTRUCTIONS FOR USE:

1. Open instrument to appropriate width for cleaning.
2. Place the spreader over the clamp handle.
3. Clean/sterilize instrument according to instrument manufacturer guidelines.
4. Clean spreader as soon as possible after use.

CLEANING AND DISINFECTION:

Manual or mechanical reprocessing with/without ultrasonic treatment

1. Chemicals and temperatures for cleaning and disinfection*.
 - Acidic / neutral / alkaline with/without added surfactants, chemical at max. 60 °C or with deionized water thermal at max. 93 °C
2. Drying
 - Max. 100°C

* It is assumed that commercially available products approved for the application are used for cleaning and disinfection. Likewise, that the recommended concentrations, exposure times and temperatures are observed. It must be ensured that no residues remain on the products. Deionized water must be used for final rinsing.

CLEANING INSPECTION:

Visually inspect before sterilization or storage to ensure the complete removal of soil from surfaces. If soil is still present, re-clean the band clamps.

STERILIZATION:

Material Compatible with sterilization according to the following parameters:

Steam Sterilization		
Cycle Type	Temperature	Cycle Time
Pre-Vacuum	Up to 121° C (250° F)	Up to 20 minutes
Pre-Vacuum	Up to 134° C (273° F)	Up to 18 minutes

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL:

Requires replacing if there are any signs of deterioration or loss of functionality. Dispose according to facility protocol.

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KEY SURGICAL LLC