Heat Protective Gloves (silicone) 09484, 09368

EN - INTENDED USE:

Key Surgical® heat protective gloves provide high protection against moderate temperature up to 240° C/464° F (intended for use during removal of products from autoclaves/wash-disinfector).

PRECAUTION:


Perform visual inspection for usability and cleanliness of the heat protective gloves before use

Note: Damaged gloves no longer provide protection and must be replaced.

APPLIED STANDARDS:

Tested according to:
EN 407:2004-11, DIN EN 420:2010-03 (EN 420:2003+A1:2009)

Tested against thermal hazards & meet minimum requirements x2xxxx
Protection against contact heat Performance level 2: 250°C/ 482°F / 15 Sek

The marking "X" instead of a number means that the glove has not been tested for the corresponding risk.

CLASSIFICATION:

Category II per Regulation (EU) 2016/425 on personal protective equipment

EU Type – Examination Certificate issued by notified body HOHENSTEIN LABORATORIES GMBH & CO. KG (NB 0555). Schloß Hohenstein 74357 Bönningheim, Germany

CLEANING AND MAINTENANCE:

Mechanical reprocessing with or without ultrasound treatment:

Cleaning and disinfecting chemicals and temperatures:

acidic / neutral / alkaline with or without surfactant additives, chemically at max. 60°C / 140°F or with demineralized water, thermally, at max. 93°C / 199°F. It is assumed that commercially available products approved for the intended application are used and that the user complies with recommended concentrations, exposure times and temperatures. Demineralized water must be used for the final rinse to ensure no residues are present post-reprocessing.

Drying: max. 100°C / 212°F

STORAGE CONDITIONS:

Store in a dry area away from direct sunlight.

DISPOSAL:

Dispose according to facility protocol.

Declaration of Conformity can be found at keysurgical.com/resources.

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KEY SURGICAL, LLC