
InkCheck® (REF 08101) Product information

Seal integrity test according to DIN EN ISO 11607-1, Annex B / ASTM F1929 for the routine monitoring of your sealing device.

Execution of the test:

1. Turn on the sealing device, select the desired program if necessary and wait until the corresponding temperature has been reached.
2. Seal a bag or tube section approx. 20 - 25 cm wide.
3. Cut the bag off approx. 5 cm above the sealed seam.
4. Remove an InkCheck from the ampoule block. Hold the ampoule with the cap facing up and unscrew the cap.
5. Empty the contents of the ampoule (2 ml) into the opened bag just above the sealed seam.
6. Allow the InkCheck to work for 5 seconds (according to ASTM F1929:2012 Annex X1*) and visually check the sealed seam for channel formation. The exposure time of 5 seconds should not be significantly exceeded as this may lead to misinterpretations (see below).
7. Channels are displayed within a few seconds by the InkCheck. Optionally, you can use the magnifying glass provided to look more closely at possible sources of error.
8. Document the result on the attached documentation sheet (copy template).

Example

Completed test

failed test

Possible misinterpretations:

- 1) The test liquid penetrates porous materials such as paper and Tyvek® very quickly. I.e. if the exposure time is too long, the liquid can creep through the porous layer behind the sealing seam, which can possibly be misinterpreted as a channel. In contrast, an actual channel would be created almost immediately after contact with the test liquid across the entire seal width.
- 2) Wrinkles and kinks in the sterile packaging may occur if the flexible porous material of the sterile packaging has been bent, folded or kinked, which may damage the fine paper fibres. This can occur, for example, during storage or transport and usually does not pose a problem for the use of the sterile packaging. However, even the smallest gaps can form inside the paper layer, which usually have no negative influence on the use of the packaging, but are detected by the test liquid.


The test liquid makes its way through the microcracks, which can possibly create a channel-like appearance. This is therefore not a defect of the sealing seam, but an occasional phenomenon which could be misinterpreted as a channel. It should be avoided if possible to bend or strongly bend sterile packaging in order not to damage the fibres of the porous side, but this often cannot be avoided. If the phenomenon is detected by the InkCheck, the packaging can still be used, as this is not critical material damage.

*** Standard requirement according to ASTM F1929:2012 Annex X1*): "... 5 seconds is sufficient as time, since channel formation is typically detected within the first few seconds. The guide value for a complete package with four sealing seams is 20 seconds (4 x 5 seconds = 20 seconds). If the reaction time of the test liquid is 20 seconds for a sealed seam, this can lead to misinterpretations because the penetration of the porous layer (Tyvek® or paper) can be too massive.**

Regulatory symbols used:

Catalog number

Storage life

Batch code

Keep dry

Keep away from direct sunlight

KEY SURGICAL, LLC
8101 WALLACE ROAD / EDEN PRAIRIE, MN 55344 USA
(US) TEL: 800.541.7995 (UK) TEL: +44(0)1628 810626
(EU) TEL: +49 4363 905900

MANUFACTURER FOR US: MANUFACTURED FOR

InkCheck (REF 08101) Documentation sheet

Sealing seam tightness test according to DIN EN ISO 11607-1 Appendix B (ASTM F1929:2012)

General Particulars:

Hospital/ department: _____

Sealing device no. / manufacturer: _____

Information on sterile packaging

Bag Tube Side fold: yes / no Material: Paper / Tyvek®

Manufacturer: _____

Test performed on/ at	set sealing temp. / program if necessary	REF-Nr.	Lot-Nr.	Format (Length x width)	Exp. Date	Test completed	if failed, description of mistakes	signature examiner
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		
						yes <input type="checkbox"/> / no <input type="checkbox"/>		

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
MANUFACTURER
FOR US: MANUFACTURED FOR

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