Instrument Tip Caps

EN - INTENDED USE:

Key Surgical® Instrument Tip Caps are designed to protect the tips of surgical instruments when not in use or during instrument sterilization.

PRECAUTION:

Does not exhibit a sharp melting point but may soften gradually during sterilization. Avoid letting caps touch or using caps in a way that will exert pressure on the cap during sterilization as this may cause them to melt. Should not be altered and are not intended for implantation. Must be removed from the instrument prior to use. Single use only; reuse may lead to patient cross contamination.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

MATERIAL:

Made of medical-grade polyvinyl chloride (PVC) material that is radiopaque and free of BPA and DEHP.

INSTRUCTIONS FOR USE:

- Does not need to be cleaned or disinfected before use. If desired, can be put through a washer or alcohol bath.
- Place the cap on the tip(s) of an appropriately sized clean instrument before sterilization and storage. Hinged instruments should be sterilized in an open position.
- Package instruments in sterile packaging device.
- Place package in sterilizer. The instrument surface under the cap will be sterile, regardless of use of vented or non-vented caps.
- Single use and can be sterilized only once.

STERILIZATION:

Validated for sterilization according to the following parameters:

Steam Sterilization			
Cycle Type	Temperature	Minimum Time – Full Cycle	Minimum Validated Dry Time
Gravity	121° C (250° F)	30 minutes	30 minutes
Pre-Vacuum	132° C (270° F)	4 minutes	20 minutes
Pre-Vacuum	134° C (273° F)	18 minutes	Not Applicable
Pre-Vacuum	134° C (273° F)	3 minutes	Not Applicable
Immediate Use Pre-Vacuum	132° C (270° F)	3 minutes	Not Applicable

100% Ethylene Oxide (EtO)	
Concentration	6% Ethylene Oxide, 94% CO ₂
Temperature	47° C
Exposure Time	90 min
Pressure	1.7 bar

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

STORAGE:

Store according to the information on the device labelling and/or this IFU.

DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
**UK REP** MDSS UK RP LTD