## Instrument Sleeves

### EN – INTENDED USE:

Key Surgical® Instrument Sleeves are designed to protect the tips and edges of surgical instruments when not in use, during transportation and during sterilization.

### CONTRAINDICATIONS:


No known contraindications and/or adverse effects.

### PRECAUTIONS:

Single use only; reuse may lead to patient cross contamination.


### INSTRUCTIONS FOR USE:

1. Does not need to be cleaned or disinfected before initial use.
2. Fold at perforation:


a.

- b. For sleeves with corner perforations, Fold the bottom corners at perforation inward to create a pocket.

3. Tuck sleeve under slit


a.


b.

4. Slide instrument into top perforated slit and weave through to bottom slit and into bottom of sleeve.

a.

b.

5. Package instruments in sterile packaging device
6. Place package in sterilizer. The instrument surface under the sleeve will be sterile.

### STERILIZATION:

Provided non-sterile and are single use.

CS-1010W, CS-2025W, CS-3050W: Validated at the following cycle parameters:

Steam Sterilization		
Cycle Type	Temperature	Minimum Time – Full Cycle
Gravity	121° C (250° F)	30 minutes
Pre-Vacuum	132° C (270° F)	4 minutes
Pre-Vacuum	134° C (273° F)	3 minutes

CS-1010W, CS-2025W, CS-3050W: Evaluated for sterilization compatibility and found to be compatible with the following parameters:

STERRAD® System and Cycle
STERRAD® 100S
STERRAD® 100S Short
STERRAD® NX Standard
STERRAD® 100NX Standard
STERRAD® 100NX Express
STERRAD® 100NX DUO
STERRAD® 100NX Flex

All other products: Validated at the following cycle parameters:

Steam Sterilization		
Cycle Type	Temperature	Minimum Time – Full Cycle
Pre-Vacuum	132° C (270° F)	3 minutes
Pre-Vacuum	134° C (273° F)	3 minutes

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

### STORAGE:

Do not store in temperatures greater than 40°C (104°F) degrees.

### DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

### WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
UK REP MDSS UK RP LTD