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Interlock LCpro **STEAM**

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Instructions for use

Interlock LCpro **STEAM**

Testing-system for the routine process control of autoclaves

DIN EN ISO 11140-1 TYP 2

REF 08120 – Set

REF 08121 – Refill 500 pcs. 4-stage indicator

REF 08125 – Process challenging device (PCD)

REF 08130 – Refill 200 pcs. 4-stage indicator for steam sterilization, 134°/18 min.

REF 08131 – Plastic spare part screw M3

Content

- 1 PCD with stainless steel holder
- 500 (2 x 250) indicators in a plastic sliding box
- 2 O-rings
- Screw with O-ring for leakage test

Advantages

- Easy handling
- Robust construction
- Logical-constructive structure following hollow instruments
- Processing of leakage tests
- Simple O-ring replacement
- Long-lasting operating times
- Reliable results

Pre-Caution

- To eliminate unnecessary environmental impact on the indicator, take the indicator out of its packaging just before application. Do not store the indicators close to heat sources, e.g. on the steriliser

General notes

- The 4 phase indicator strip and the PCD form the indicator system.
- The PCD must have the same temperature as the sterile load.
- A visual inspection on the integrity and the fit of the O-ring has to be performed, before using the PCD.
- The leakage test needs to be undertaken every 250 cycles. Therefore place the screw (M3) with the O-ring into the steam inlet opening at the bottom of the PCD. Take attention that you put the screw vertical into the PCD because the sealing ring has to be into the notch (Figure 3). Then replace the O-ring of the indicator mount and place an indicator strip into it. Run a sterilisation process. The result should be no colour change of the indicator. If the indicator colour changes, repeat the procedure. If the colour of the indicator strip changes again, contact Interlock.

Application:


The testing-system is aimed at the process monitoring in steam sterilisers in the healthcare according to the standard DIN EN 285. That means that air gets removed by fractional pre-vacuum steps. Within the German healthcare the general rule in terms of temperature/time combination is 134°C and 5 min.

The PCD's function is to examine the complete air removal and so the steam penetration (analogue DIN EN 867-5). The indicator's colour change depends on the steam impact. A colour change to black takes place when the amount of saturated steam was impacting on the indicator that complies with a steam temperature of 134°C for 5 minutes pure exposure time. The device can also be used at 121°C with a holding time of min. 15 minutes.

Attention: Before first use, the screw must be removed!

Equipping of the PCD to the indicator system

Screw off the indicator mounting capsule. Bend the indicator strip at the fold so the printed sides are facing each other. Slide it into the gap like shown in the figure below. Screw the PCD together and place it close to the door in the lower part of the chamber.

When the program has finished, remove the indicator from the PCD and evaluate the colour change. The indicator is self-adhesive and can be paste into the documentation. For this purpose remove the backing layer at the backside perforation.

Information on indicator evaluation and possible sources of errors

The indicator changes colour gradually from cream to grey to black when being affected by saturated steam depending on time.

Appearance of the indicator	Notes on evaluation
 cream	Initial state of the indicator
 black	Proper sterilisation process operation
 colour change light to dark	The gradient can occur in different colour intensities. It occurs when the indicator was influenced by no or less steam, i.e. the PCD is not completely free from inert gas, e.g. air. Different causes are possible: e.g. leakage, lack of ventilation, insufficient steam quality, overheating
 equal greyish brown	The air was removed completely from the PCD but lack of steam condensation on the indicator. Leads to an insufficient colour depth and indicates overheated steam or short exposure time.
 irregular black with grey zones	Irregular steam condensation Possible causes could be: <ul style="list-style-type: none"> • The steriliser was not preheated or overloaded • The PCD and/or load not conditioned • The speed of pressure change is too high.

Figure 3

Storage conditions: see product label
PCD Holder

The holder can be mounted on the PCD by slightly extending both wire ends and pushing the holder over the PCD

Catalog number

Serial number

Batch code

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Humidity limit,
limitation

Temperature limitation

Keep away from
direct sunlight

Keep dry