## Plastic Sterilization Trays

### EN - INTENDED USE:

Key Surgical® Plastic Sterilization Trays are designed to store and protect delicate instruments during sterilization. These trays are not intended for use with immediate-use steam sterilization (IUSS). The trays feature holes and additional slots in the lid and base to increase sterilant and air circulation. Holes in the silicone mat inserts align with the base to ensure adequate drying and prevention of wet packs. The use of the plastic trays may prolong the life of surgical instruments by protecting the items during storage, transportation, and sterilization.

### CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

### INSTRUCTIONS FOR USE:

1. Reusable and should be decontaminated, cleaned, and inspected prior to each use according to the instructions below.
2. Configure the tray and mat so the mat pins point upward and the holes in the tray base and mat align to accommodate the instruments to be sterilized. The mat is used to help stabilize, separate, and protect instruments. Instruments should be arranged in a single layer with adequate spacing to prevent them from contacting each other.
3. An approved sterility assurance device (i.e. chemical indicator) should be placed inside and on each level of the tray in a location easily visible to the user.
4. Cover the tray with the lid and press down until the side tabs click shut, securing the lid.
5. Wrap with an approved non-woven textile wrap or place in any other FDA cleared packaging system (i.e. sterilization pouch or container). Do not wrap more than one tray per wrap.
6. The tray should always be placed flat on the sterilizer rack prior to sterilization. Stacking of trays inside the sterilization chamber is not recommended.
7. Sterilize the loaded wrapped tray in accordance with the instrument device manufacturer's IFU.
8. After sterilization, move wrapped tray to a sterile storage area for storage or transport to the point of use location.
9. Clean as soon as possible after use. If cleaning must be delayed, immerse in an enzymatic solution or water to prevent drying and encrustation of surgical soil. Remove excessive soil with a disposable wipe.

### MANUAL CLEANING:

1. Pre-rinse under cold tap water for one (1) minute to remove gross debris.
2. Soak for a minimum of two (2) minutes in a pH neutral detergent, prepared in accordance with the manufacturer's IFU.
3. Rinse under cold tap water for one (1) minute.
4. Ultrasonically clean for a minimum of five (5) minutes in a neutral pH detergent, prepared in accordance with the manufacturer's IFU.
5. Rinse under cold tap water for one (1) minute.
6. Visually inspect to ensure complete removal of soil from surfaces. No visible soil should be observed.
7. If soil is still visible, repeat above steps until free from visible soil.

### AUTOMATED CLEANING:

It may be necessary to manually clean prior to automated processing to improve the removal of adherent soil. Follow the previous instructions for manual cleaning.

1. Run the automatic wash cycle – minimum cycle parameters:
  - 1 minute cold pre-rinse
  - 5 minute enzyme wash at 43° C (109° F) minimum temperature
  - 1 minute cold rinse
  - 7 minute dry at 90° C (194° F) minimum temperature
2. Visually inspect to ensure complete removal of soil from surfaces. No visible soil should be observed.
3. If soil is still visible, repeat above steps until free from visible soil.

### INSPECTION:

Inspect for damage: Requires replacing if there are any signs of deterioration or loss of functionality.

### STERILIZATION:

Validated for sterilization efficacy according to the following parameters:

Steam Sterilization		
Cycle Type	Temperature	Minimum Time – Full Cycle
Gravity	121° C (250° F)	30 minutes
Pre-Vacuum	132° C (270° F)	4 minutes
Pre-Vacuum	134° C (273° F)	3 minutes

STERRAD® System and Cycle	
STERRAD® System and Cycle within the United States	STERRAD® System and Cycle outside United States
STERRAD® 100S	STERRAD® 100 Short cycle
STERRAD® NX Standard cycle	
STERRAD®100NX Standard cycle	

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

### DISPOSAL:

Requires replacing if there are any signs of deterioration or loss of functionality. Dispose according to facility protocol.

### WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
MDSS UK RP LTD