Silicone Finger Mat

EN - INTENDED USE:

Key Surgical® Silicone Finger Mat is designed to be placed in the bottom of a tray to hold, organize, and protect surgical instruments during sterilization.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

INSTRUCTIONS FOR USE:

1. Reusable and should be decontaminated, cleaned, and inspected prior to each use according to the instructions below.
2. May cut mat to the desired length to fit tray.
3. Configure the tray and mat so the mat pins point upward. The mat is used to help stabilize, separate, and protect instruments. Instruments should be arranged in a single layer with adequate spacing to prevent them from contacting each other.
4. Clean as soon as possible after use. If cleaning must be delayed, immerse in an enzymatic solution or water to prevent drying and encrustation of surgical soil. Remove excessive soil with a disposable wipe.

MANUAL CLEANING:

1. Pre-rinse under cold tap water for one (1) minute to remove gross debris.
2. Soak for a minimum of two (2) minutes in a pH neutral detergent, prepared in accordance with the manufacturer's IFU.
3. Rinse under cold tap water for one (1) minute.
4. Ultrasonically clean for a minimum of five (5) minutes in a neutral pH detergent, prepared in accordance with the manufacturer's IFU.
5. Rinse under cold tap water for one (1) minute.
6. Visually inspect to ensure complete removal of soil from surfaces. No visible soil should be observed.
7. If soil is still visible, repeat above steps until free from visible soil.

AUTOMATED CLEANING:

It may be necessary to manually clean prior to automated processing to improve the removal of adherent soil. Follow the previous instructions for manual cleaning.

1. Run the automatic wash cycle – minimum cycle parameters:
 - 1 minute cold pre-rinse
 - 5 minute enzyme wash at 43° C (109° F) minimum temperature
 - 1 minute cold rinse
 - 7 minute dry at 90° C (194° F) minimum temperature
2. Visually inspect to ensure complete removal of soil from surfaces. No visible soil should be observed.
3. If soil is still visible, repeat above steps until free from visible soil.

INSPECTION:

Inspect for damage: Requires replacing if there are any signs of deterioration or loss of functionality.

STERILIZATION:

Validated for sterilization efficacy according to the following parameters:

Steam Sterilization		
Cycle Type	Temperature	Minimum Time – Full Cycle
Gravity	121° C (250° F)	30 minutes
Pre-Vacuum	132° C (270° F)	4 minutes
Pre-Vacuum	134° C (273° F)	3 minutes

STERRAD® System and Cycle	
STERRAD® System and Cycle within the United States	STERRAD® System and Cycle outside United States
STERRAD® 100S	STERRAD® 100 Short cycle
STERRAD® NX Standard cycle	
STERRAD® 100NX Standard cycle	

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL:

Requires replacing if there are any signs of deterioration or loss of functionality. Dispose according to facility protocol.

WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
MDSS UK RP LTD