## Choice Dry Green bags, TCOVER and TCOVERDRY Bags

### **EN** - INTENDED USE:

Key Surgical® Choice Dry Green, TCOVER bags, TCOVERDRY bags are intended for use in conjunction with sterile barrier wraps during sterilisation, transportation, and storage to prevent damage to the sterilisation wrap.

### **CONTRAINDICATIONS:**

No known contraindications and/or adverse effects.

### **PRECAUTIONS:**

Visually inspect before use to ensure that they are free from visible contaminants or any damage such as rips or tears.

Single use only; reuse may lead to patient cross contamination.

### **INSTRUCTIONS FOR USE:**

1. Product does not need to be cleaned or disinfected prior to initial use.
2. Wrap the basket or tray containing the medical devices to be sterilised in inner barrier sterilisation wrap.
3. Place the wrapped basket or tray inside the bag.
4. Fold the lip over three times.
5. Secure with autoclave tape.

### **STERILISATION:**

Evaluated for sterilisation compatibility and found to be compatible with steam sterilisation.

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials, and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

### **DISPOSAL:**

Used: Dispose after use according to facility protocol. Do not reuse.

Unused: Can be disposed of in the plastic recycling stream.

### **WARNING:**

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
MDSS UK RP LTD