

Indicator Dot Tags

EN - INTENDED USE:

Indicator Dot Tags are used to identify instruments, trays, carts, and other equipment as they move from the sterilization to staging areas of surgery. Indicator dots on tags change color to indicate exposure to steam sterilization.

PRECAUTIONS:

Single use only; reuse may lead to patient cross contamination.

CAUTION:

Store at 20-25°C (68-77°F), excursions 15-30°C (59-86°F) allowed. Relative humidity range is 30-70%. Significant changes in storage conditions for prolonged periods can have an adverse effect on the product. Extreme storage conditions such as exposure to direct sunlight and/or storage on top of or near a heat source should be avoided. Store away from hydrogen peroxide and other sterilant.

INSTRUCTIONS FOR USE:

1. Does not need to be cleaned or disinfected before use.
2. Tags can be sterilized according to the parameters detailed below. Tags are single-use and can be sterilized only once.
3. Verify that the indicator dot changes from blue to dark. Typical color change will result in dark brown to black after exposure to steam sterilization.
4. Tags can be broken and removed by hand or with the use of cutters.


STERILIZATION:

Provided non-sterile. Validated for steam sterilization at the following cycle parameters:

Steam Sterilization		
Cycle Type	Temperature	Minimum Time – Full Cycle
Pre-Vacuum	132° C (270° F)	4 minutes

DISPOSAL:

Dispose after use according to facility protocol. Do not reuse.

 MANUFACTURER: KEY SURGICAL
8101 WALLACE ROAD / EDEN PRAIRIE, MN 55344 USA
(US) PH: 800.541.7995 KEYSURGICAL.COM
(UK) TEL: +44(0)1628 810626 KEYSURGICAL.CO.UK
(EU) TEL: +49 4363 905900 KEYSURGICAL.DE



CONSULT
INSTRUCTIONS
FOR USE



CATALOG
NUMBER



LOT
NUMBER



USE-BY DATE