Acetabular Reamer Brushes

EN - INTENDED USE:

Key Surgical® Acetabular Reamer Brush is intended to clean dirty acetabular bone reamers.

CONTRAINDICATIONS:

No known contraindications and/or adverse effects.

PREPARATION:

Must be cleaned and disinfected prior to reuse.

CLEANING & DISINFECTION:

Clean after each use and disinfect between cases or daily, at a minimum. Can be cleaned manually or using an automatic washer/disinfector.

INSPECTION:

Inspect between uses and replaced when worn, frayed, bent or otherwise damaged. Damaged brushes cannot be cleaned effectively. Worn or damaged bristles and handles are ineffective in cleaning and may damage the device. Requires replacing if there are any signs of deterioration or loss of functionality

STERILIZATION:

This sterilization is intended for one time use of brushes in a sterile setting and is not for brushes that have already been used.

New, unused brushes are provided non-sterile and have been validated for sterilization efficacy using gravity and pre-vacuum steam sterilization according to the following parameters:

Steam Sterilization*		
Cycle Type	Temperature	Minimum Time –Full Cycle
Gravity	121° C (250° F)	30 minutes
Pre-Vacuum	132° C (270° F)	4 minutes
Pre-Vacuum	134° C (273° F)	3 minutes

*The following cannot be sterilized: 26023

Due to the steam sterilization process, there may be a slight discoloration of the nylon brush bristles; however, this does not impair the fit, form, or function of the brush.

STORAGE:

Store brushes in a manner that will reduce cross-contamination.

DISPOSAL:

Dispose according to facility protocol.

WARNING:

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
MDSS UK RP LTD