Surgical Tray Corner Protectors (Silicone)

EN - INTENDED USE: Key Surgical® Tray Corner Protectors are designed for use on the bottom and/or top of surgical instrument trays to provide protection during sterilization, transportation, and/or storage.

CONTRAINDICATIONS: No known contraindications and/or adverse effects.

PRECAUTIONS: Visually inspect the devices before use to ensure that they are free from visible contaminants or any damage such as rips or tears.

INSTRUCTIONS FOR USE:

- Does not need to be cleaned or disinfected before initial use. If desired, can be cleaned and disinfected prior to initial use.
- Place the corner protectors on the top and bottom corners of the tray before wrapping with sterilization wrap.
- Wrap the tray according to standard wrapping procedures. Appropriate sterilant contact will occur on the surface of the corner protectors as well as the tray surface under the corner protectors.
- Reusable and should be cleaned and decontaminated prior to reuse according to the instructions below.

MANUAL CLEANING:

- Pre-rinse under cold tap water for one (1) minute to remove gross debris.
- Soak for a minimum of two (2) minutes in a pH neutral detergent, prepared in accordance with the manufacturer's IFU.
- Rinse under cold tap water for one (1) minute.
- Ultrasonically clean for a minimum of five (5) minutes in a neutral pH detergent, prepared in accordance with the manufacturer's IFU.
- Rinse under cold tap water for one (1) minute.
- Visually inspect to ensure complete removal of soil from surfaces. No visible soil should be observed.
- If soil is still visible, repeat above steps until free from visible soil.

AUTOMATED CLEANING: It may be necessary to manually clean prior to automated processing to improve the removal of adherent soil. Follow the previous instructions for manual cleaning.

- Load in a mesh wire tray with concave surface facing downward to facilitate draining.
- Run the automatic wash cycle – minimum cycle parameters:
 - 1 minute cold pre-rinse
 - 5 minute enzyme wash at 43° C (109° F) minimum temperature
 - 1 minute cold rinse
 - 7 minute dry at 90° C (194° F) minimum temperature
- Visually inspect to ensure complete removal of soil from surfaces. No visible soil should be observed.
- If soil is still visible, repeat above steps until free from visible soil.

INSPECTION: Inspect for damage: Requires replacing if there are any signs of deterioration or loss of functionality.

STERILIZATION: Evaluated for sterilization compatibility and found to be compatible with the following parameters:

Steam Sterilization			
Cycle Type	Temperature	Minimum Time – Full Cycle	Minimum Validated Dry Time
Gravity	121° C (250° F)	30 minutes	20 minutes
Pre-Vacuum	132° C (270° F)	4 minutes	20 minutes
Pre-Vacuum	134° C (273° F)	3 minutes	20 minutes
STERRAD® System and Cycle			
STERRAD® 100S/100S Short			
STERRAD® NX Standard/Advanced			
STERRAD® 100NX Standard/Express/DUO/Flex			
V-PRO System (CS9113)			
V-PRO® maX Lumen/Non Lumen/Flexible			
V-PRO® maX 2 Lumen/Non Lumen/Flexible/Fast Non Lumen			
V-PRO® 60 Lumen/Non Lumen/Flexible			
V-PRO® s2 Lumen/Non Lumen/Flexible			
V-PRO® 1/1 Plus			

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL: Tested to withstand a minimum of 15 sterilization cycles with no impact to the fit, form or function. Requires replacing if there are any signs of deterioration or loss of functionality. Dispose according to facility protocol.

WARNING: Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
UK REP MDSS UK RP LTD