## Stainless Steel Baskets

**EN - INTENDED USE:** Key Surgical® Stainless Steel Baskets (tray) are designed to contain and protect surgical instruments during the sterilization process.

**PRECAUTIONS:** Do not overload.

Overloading may cause damage.

Visually inspect prior to use. Check for loose wires or burrs that may puncture, or tear sterilization wrap or pouch.

Inspect handles to ensure they are secure prior to use and transport.

Avoid prolonged exposure to saline to minimize the chance of corrosion.

**CONTRAINDICATIONS:** No known contraindications and/or adverse effects.

### INSTRUCTIONS FOR USE:

- Reusable and should be decontaminated, cleaned, and inspected prior to each use according to the instructions below.
- Place unpouched instruments in tray. An approved liner, mat, or other positioning device can be used to accommodate the instruments to be sterilized.
- Wrap with an approved non-woven textile wrap or place in any other FDA cleared packaging system (i.e. sterilization pouch or container). Do not wrap more than one tray per wrap.
- The tray should always be placed flat on the sterilizer rack prior to sterilization. Stacking of trays inside the sterilization chamber is not recommended.
- Sterilize the loaded wrapped tray in accordance with the instrument device manufacturer's IFU.
- After sterilization, move wrapped tray to a sterile storage area for storage or transport to the point of use location.
- Clean as soon as possible after use. If cleaning must be delayed, immerse in an enzymatic solution or water to prevent drying and encrustation of surgical soil. Remove excessive soil with a disposable wipe.

### MANUAL CLEANING:

- Pre-rinse under cold tap water for one (1) minute to remove gross debris
- Soak for a minimum of two (2) minutes in a pH neutral detergent, prepared in accordance with the manufacturer's IFU.
- Rinse under cold tap water for one (1) minute.
- Ultrasonically clean for a minimum of five (5) minutes in a neutral pH detergent, prepared in accordance with the manufacturer's IFU.
- Rinse under cold tap water for one (1) minute.
- Visually inspect to ensure complete removal of soil from surfaces. No visible soil should be observed.
- If soil is still visible, repeat above steps until free from visible soil.

**AUTOMATED CLEANING:** It may be necessary to manually clean prior to automated processing to improve the removal of adherent soil. Follow the previous instructions for manual cleaning.

- Run the automatic wash cycle – minimum cycle parameters:
  - 1 minute cold pre-rinse
  - 5 minute enzyme wash at 43° C (109° F) minimum temperature
  - 1 minute cold rinse
  - 7 minute dry at 90° C (194° F) minimum temperature
- Visually inspect to ensure complete removal of soil from surfaces. No visible soil should be observed.
- If soil is still visible, repeat above steps until free from visible soil.

**INSPECTION:** Inspect for damage: Requires replacing if there are any signs of deterioration or loss of functionality.

STERILIZATION: Validated for sterilization according to the following parameters: Steam Sterilization		
Cycle Type	Temperature	Minimum Time – Full Cycle
Gravity	121° C (250° F)	30 minutes
Pre-Vacuum	132° C (270° F)	4 minutes
Pre-Vacuum	134° C (275° F)	3 minutes
Immediate Use Pre-Vacuum	132° C (270° F)	3 minutes
100% Ethylene Oxide (EtO)		
Concentration		725 mg/L
Temperature		55° C
Exposure Time		60 minutes
Humidity		50-80%

Evaluated for sterilization compatibility and found to be compatible with the following parameters:

STERRAD® System and Cycle	
STERRAD® 100S/100S Short	
STERRAD® NX Standard	
STERRAD® 100NX Standard/Express/DUO/Flex	
STERIZONE® System	STERIZONE® VP4
V-PRO System (MTXXXX)	
V-PRO® maX Lumen/Non Lumen/Flexible	
V-PRO® maX 2 Lumen/Non Lumen/Flexible/Fast Non Lumen	
V-PRO® 60 Lumen/Non Lumen/Flexible	
V-PRO® s2 Lumen/Non Lumen/Flexible	

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

**DISPOSAL:** Requires replacing if there are any signs of deterioration or loss of functionality. Dispose according to facility protocol.

**WARNING:** Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

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KEY SURGICAL LLC
MDSS UK RP LTD