Band Clamps

EN - INTENDED USE: Key Surgical band clamps are designed to keep the joints and jaws of bone forceps and other instruments open during the cleaning and sterilization process.

CONTRAINDICATIONS: No known contraindications and/or adverse effects.

INSTRUCTIONS FOR USE:

1. Open instrument to appropriate width for cleaning.
2. Place band clamp over instrument handles.
3. Clean/sterilize instrument according to instrument manufacturer guidelines.
4. Clean band clamps as soon as possible after use. If cleaning must be delayed, immerse in an enzymatic solution or water to prevent drying and encrustation of surgical soil. Remove excessive soil with a disposable wipe.

MANUAL CLEANING:

1. Pre-rinse under cold tap water for one (1) minute to remove gross debris
2. Soak for a minimum of two (2) minutes in a pH neutral detergent, prepared in accordance with the manufacturer's instructions for use.
3. Rinse under cold tap water for one (1) minute.
4. Ultrasonically clean for a minimum of five (5) minutes in a neutral pH detergent, prepared in accordance with the manufacturer's instructions for use.
5. Rinse under cold tap water for one (1) minute.

AUTOMATED CLEANING: It may be necessary to manually clean prior to automated processing to improve the removal of adherent soil. Follow the previous instructions for manual cleaning.

1. Run the automatic wash cycle – minimum cycle parameters:
 - 1 minute cold pre-rinse
 - 5 minute enzyme wash at 43° C minimum temperature
 - 1 minute cold rinse
 - 7 minute dry at 90° C minimum temperature

CLEANING INSPECTION: Visually inspect before sterilization or storage to ensure the complete removal of soil from surfaces. If soil is still present, re-clean the band clamps.

STERILIZATION: Material Compatible with sterilization according to the following parameters:

Steam Sterilization		
Cycle Type	Temperature	Cycle Time
Gravity	Up to 121° C (250° F)	Up to 30 minutes
Pre-Vacuum	Up to 121° C (250° F)	Up to 20 minutes
Pre-Vacuum	Up to 134° C (273° F)	Up to 18 minutes
100% Ethylene Oxide (EtO)		
Concentration	735 mg/L	
Temperature	55° C	
Exposure Time	60 minutes	
Humidity	50-70%	
V-PRO System		
V-PRO® maX Lumen/Non Lumen/Flexible V-PRO® maX 2 Lumen/Non Lumen/Flexible/Fast Non Lumen V-PRO® 60 Lumen/Non Lumen/Flexible V-PRO® s2 Lumen/Non Lumen/Flexible		

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL: Requires replacing if there are any signs of deterioration or loss of functionality. Dispose according to facility protocol.

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KEY SURGICAL LLC