Reusable Tray Tags

[Metal Tray Tag MTXXXX, MTXXXX10;

Plastic Tray Tag PTXXXX], Plastic Insertion Tag, Plastic Mesh Basket Labels

(Note: See BROC3138 for Key Surgical® Indicator Dot Tags, see BROC3151 for Key Surgical® Single-Use Tags)

EN - INTENDED USE: The Key Surgical® Tray Tags are intended to be attached to the side of a tray or container for purposes of color coding or affixing a label.

CONTRAINDICATIONS: No known contraindications and/or adverse effects.

NOTE: Stacking heavy objects on plastic tags may cause cracking/breaking of tag.

INSTRUCTIONS FOR USE:

- Attach tag to tray by:
 - Sliding the tray clip over the edge of the tray, OR:
 - Alternatively, some tag designs may require a fastener such as a ring, clip, or tie to connect the tag to the tray.
- Self-adhesive labels may be affixed to the tag.
- Clean as soon as possible after use. If cleaning must be delayed, immerse in an enzymatic solution, apply pre-treatment product, or water to prevent drying and encrustation of surgical soil. Remove excessive soil with a disposable wipe.

MANUAL CLEANING:

- Pre-rinse under cold tap water for one (1) minute to remove gross debris
- Soak for a minimum of two (2) minutes in a pH neutral detergent, prepared in accordance with the manufacturer's instructions for use.
- Rinse under cold tap water for one (1) minute.
- Ultrasonically clean for a minimum of five (5) minutes in a neutral pH detergent, prepared in accordance with the manufacturer's instructions for use.
- Rinse under cold tap water for one (1) minute.

AUTOMATED CLEANING: It may be necessary to manually clean prior to automated processing to improve the removal of adherent soil. Follow the previous instructions for manual cleaning.

- Run the automatic wash cycle – minimum cycle parameters:
 - 1 minute cold pre-rinse
 - 5 minute enzyme wash at 43° C minimum temperature
 - 1 minute cold rinse
 - 7 minute dry at 90° C minimum temperature

CLEANING INSPECTION: Visually inspect before sterilization or storage to ensure the complete removal of soil from surfaces. If soil is still present, re-clean.

STERILIZATION: Provided non-sterile and are compatible with sterilization according to the following parameters:

Tag Type	Steam Sterilization		
	Cycle Type	Temperature	Maximum Time – Full Cycle
Metal Tray Tag (MTXXXX, MTXXXX10) Plastic Tray Tag (PTXXXX)	Gravity	121° C (250° F)	30 minutes
	Pre-Vacuum	132° C (270° F)	4 minutes
Metal Tray Tag (01444, 01445, 01449, 01452-01454, 01493, 01495, 01520, 01559, 01588, 01590)	Gravity	121° C (250° F)	30 minutes
	Pre-Vacuum	132° C (270° F)	4 minutes
Plastic Tray/Insertion Tag/Mesh Basket Labels (01423-01440, 01547, 01550-01556, 01561-01564, 01566-01571, 01573, 01575-01576,	Gravity	121° C (250° F)	30 minutes

01578-01587, 015970-01599, 01900-01947, 01956-01968, 01970, 01976-01979, 01980-01989, 01992-01994, 01996-01999, 02200-02226, 02233-02237, 02239-02247, 02249-02252)	Pre-Vacuum	132° C (270° F)	4 minutes
100% Ethylene Oxide (EtO)			
Metal Tray Tag (MTXXXX, MTXXXX10)	Concentration	725 mg/L	
	Temperature	55° C	
	Exposure Time	60 minutes	
	Humidity	50-70%	
STERIZONE® System			
Metal Tray Tag (MTXXXX, MTXXXX10)	STERIZONE® VP4		
STERRAD® System and Cycle			
Metal Tray Tag (MTXXXX, MTXXXX10) Plastic Tray Tag (PTXXXX)	STERRAD® 100S/100S Short		
	STERRAD® NX Standard/Advanced		
Plastic Tray Tag (PTXXXX)	STERRAD® 100NX Standard/Express/DUO/Flex		
	V-PRO® System		
Plastic Tray Tag (PTXXXX)	V-PRO® maX Lumen/Non Lumen/Flexible V-PRO® maX 2 Lumen/Non Lumen/Flexible/Fast Non Lumen V-PRO® 60 Lumen/Non Lumen/Flexible V-PRO® s2 Lumen/Non Lumen/Flexible		

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL: Requires replacing if there are any signs of deterioration or loss of functionality. Dispose according to facility protocol.

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KEY SURGICAL, LLC