

## Magnetic Mats

### INTENDED USE:

Key Surgical® Magnetic Mats are intended to be used in a hands-free transfer of stainless steel surgical sharps and instruments during procedures.

### PRECAUTIONS:

- These instructions are for general use and application of the magnetic mats. These instructions do not represent instructions for any procedures during which the products may be associated.
- The magnetic mats must be sterilized prior to use and before entering the sterile field.

### CAUTION:

- Do not use product in MRI controlled environments. Product is magnetic and not MRI safe.

### CONTRAINDICATION:

- Caution should be used during procedures with materials which may be adversely affected by magnetic fields.

### CLEANING AND DISINFECTION:

Clean the magnetic mats prior to re-use.

1. Prepare a mild, neutral pH enzymatic cleaning solution in accordance with the manufacturer's instructions.
2. Soak soiled mat according to the instructions for use for the enzymatic solution to loosen or remove soil. Abrasive cleaning chemicals, brushes or pads should not be used on the mat.
3. After cleaning with detergent, the mat should be rinsed with demineralized water.
4. Visually inspect the mat before sterilization or storage to ensure the complete removal of soil. Inspection should verify cleanliness, lack of damage and proper function.

### STERILIZATION:

Magnetic mats are provided non-sterile and have been validated for sterilization according to applicable international process standards and guidance for the following methods and parameters.

<i>Cycle Type (Steam)</i>	<i>Temperature</i>	<i>Minimum Time – Full Cycle</i>
Gravity	121 ° C (250° F)	30 minutes
Pre-Vacuum	132° C (270° F)	4 minutes
Pre-Vacuum	134° C (273° F)	3 minutes

### STORAGE:

The magnetic mats must be stored at room temperature in dry conditions. Do not expose the mats to direct sunlight.

### DISPOSAL:

The magnetic mats have been tested to withstand a minimum of 15 sterilization cycles with no impact to the fit, form or function; however, the magnetic mats have a life span and will require replacing if there are any signs of deterioration or loss of functionality. Dispose of product following facility policy for disposal.

 **MANUFACTURER: KEY SURGICAL LLC**  
8101 WALLACE ROAD / EDEN PRAIRIE, MN 55344 / USA  
PH: 800.541.7995 KEYSURGICAL.COM

**REF**  
CATALOG  
NUMBER

**LOT**  
LOT  
NUMBER

