

Reusable Tray Tags

[Metal Tray Tag MT-XXXX, MT-XXXX-10;
Plastic Tray Tag PT-XX-XXX]

(Note: See BROC-3138 for Key Surgical® Indicator Dot Tags, see BROC-3151 for Key Surgical® Single-Use Tags)

INTENDED USE:

The Key Surgical® Tray Tags are intended to be attached to the side of a tray or container for purposes of color coding or affixing a label.

CONTRAINDICATIONS:

There are no known contraindications and/or adverse effects.

NOTE: Stacking heavy objects on plastic tags may cause cracking/breaking of tag.

INSTRUCTIONS FOR USE:

1. Attach tag to tray by:
 - a. Sliding the tray clip over the edge of the tray, OR:
 - b. Alternatively, some tag designs may require a fastener such as a ring, clip, or tie to connect the tag to the tray.
2. Self-adhesive labels may be affixed to the tag.

POINT OF USE CARE:

Clean tray tag as soon as possible after use. If cleaning must be delayed, immerse in an enzymatic solution, apply pre-treatment product, or water to prevent drying and encrustation of surgical soil. Remove excessive soil with a disposable wipe.

MANUAL CLEANING:

1. Pre-rinse under cold tap water for one (1) minute to remove gross debris
2. Soak for a minimum of two (2) minutes in a pH neutral detergent, prepared in accordance with the manufacturer's instructions for use.
3. Rinse under cold tap water for one (1) minute.
4. Ultrasonically clean for a minimum of five (5) minutes in a neutral pH detergent, prepared in accordance with the manufacturer's instructions for use.
5. Rinse under cold tap water for one (1) minute.

AUTOMATED CLEANING:

It may be necessary to manually clean prior to automated processing to improve the removal of adherent soil. Follow the previous instructions for manual cleaning.

1. Run the automatic wash cycle – minimum cycle parameters:
 - 1 minute cold pre-rinse
 - 5 minute enzyme wash at 43° C minimum temperature
 - 1 minute cold rinse
 - 7 minute dry at 90° C minimum temperature

CLEANING INSPECTION:

Visually inspect before sterilization or storage to ensure the complete removal of soil from surfaces. If soil is still present, re-clean the tray tag.

STERILIZATION:


Tray tags are provided non-sterile and may be sterilized using the following cycle parameters:

| Tag Type | Steam Sterilization | | |
|---|---|-----------------|---------------------------|
| | Cycle Type | Temperature | Minimum Time – Full Cycle |
| Metal Tray Tag (MT-XXXX, MT-XXXX-10) Plastic Tray Tag (PT-XX-XX) | Gravity | 121° C (250° F) | 30 minutes |
| | Pre-Vacuum | 132° C (270° F) | 4 minutes |
| Metal Tray Tag (MT-XXXX, MT-XXXX-10) | Pre-Vacuum | 132° C (270° F) | 3 minutes |
| | Immediate Use Pre-Vacuum | 132° C (270° F) | 3 minutes |
| 100% Ethylene Oxide (EtO) | | | |
| Metal Tray Tag (MT-XXXX, MT-XXXX-10) Plastic Tray Tag (PT-XX-XX) | Concentration | | 725 mg/L |
| | Temperature | | 55° C |
| | Exposure Time | | 60 minutes |
| | Humidity | | 50-70% |
| STERIZONE® System | | | |
| Metal Tray Tag (MT-XXXX, MT-XXXX-10) | STERIZONE® VP4 | | |
| STERRAD® System and Cycle | | | |
| Metal Tray Tag (MT-XXXX, MT-XXXX-10) Plastic Tray Tag (PT-XX-XX) | STERRAD® 50/100S/100S Short/200/200 Short STERRAD® NX Standard/Advanced STERRAD® 100NX Standard | | |

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL:

Tray Tags have a life span and will require replacing if there are any signs of deterioration or loss of functionality. Dispose of product following facility policy for disposal.

 **MANUFACTURER: KEY SURGICAL LLC**
8101 WALLACE ROAD / EDEN PRAIRIE, MN 55344 / USA
PH: 800.541.7995 KEYSURGICAL.COM

LOT LOT NUMBER **REF** CATALOG NUMBER