

Silicone Finger Mat

INTENDED USE:

Key Surgical® Silicone Finger Mat is designed to be placed in the bottom of a tray to hold, organize, and protect surgical instruments during sterilization.

CONTRAINDICATIONS

There are no known contraindications and/or adverse effects.

INSTRUCTIONS FOR USE:

1. The Silicone Finger Mat should be cleaned prior to first use in accordance with the cleaning instructions below.
2. Cut mat to the desired length.
3. Configure the mat properly to where the fingers point upward to accommodate the instruments that will be sterilized. The mat should be used to help hold, organize, and protect instruments. Instruments should be arranged in a single layer with adequate spacing to prevent them from contacting each other.
4. The Silicone Finger Mat is reusable and should be cleaned and decontaminated prior to reuse.

CLEANING AND DISINFECTION:

1. Prepare a mild, neutral pH enzymatic cleaning solution in accordance with the manufacturer's instructions.
2. Soak soiled mat according to the instructions for use for the enzymatic solution to loosen or remove soil. Hard to reach places such as between the pins on the mat can be cleaned using a nylon bristled brush (i.e., Key Surgical part N-2100, N-3000, N-4000). Abrasive cleaning chemicals, brushes or pads should not be used on the mat.
3. After cleaning with detergent, the mat should be rinsed thoroughly with warm tap water for at least one minute. A soft absorbent towel (i.e., Key Surgical part CR-21214) should be used to dry the mat or allow adequate time to air dry.
4. Visually inspect the mat before sterilization or storage to ensure the complete removal of soil. Inspection should verify cleanliness, lack of damage and proper function.

STERILIZATION:

The Silicone Finger Mat has been validated for sterilization according to applicable international process standards and guidance for the following methods and parameters:


Cycle Type (Steam)	Temperature	Minimum Time – Full Cycle
Pre-Vacuum	132° C (270° F)	4 minutes

STERRAD® System and Cycle within the United States	STERRAD® System and Cycle outside United States
STERRAD® 100S	STERRAD® 100 Short cycle
STERRAD® NX Standard cycle	
STERRAD®100NX Standard cycle	

It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process.

DISPOSAL:

The silicone finger mat has a life span and will require replacing if there are any signs of deterioration or loss of functionality. Dispose of product following facility policy for disposal.

 MANUFACTURER: KEY SURGICAL LLC
8101 WALLACE ROAD / EDEN PRAIRIE, MN 55344 / USA
PH: 800.541.7995 KEYSURGICAL.COM

 REF CATALOG
NUMBER

 LOT
NUMBER

 CONSULT
INSTRUCTIONS
FOR USE



 EC  REP
MDSS GMBH
SCHIFFGRABEN 41, D-30175 / HANNOVER, GERMANY