


2200 RENAISSANCE BASING VIEW, BASINGSTOKE,
HAMPSHIRE, ENGLAND, RG21 4EQ

Tel: +44 (0) 1628 810626
Fax: +44 (0) 1628 810760

sales@keysurgical.co.uk keysurgical.co.uk



STERILE PROCESSING
EDITION 31



Ordering Information

PLACING ORDERS AND PRODUCT ADVICE

If you require advice or wish to place an order, you can contact Key Surgical from Monday to Friday between 8:30 am to 5:00 pm as follows

- 1

TELEPHONE

+44 (0) 1628 810626
- 2

FAX

+44 (0) 1628 810760
- 3

E-MAIL

sales@keysurgical.co.uk
- 4

WEBSITE

keysurgical.co.uk
- 5

MAIL

Beech House, Knaves Beech Business Centre,
Loudwater, High Wycombe, Bucks,
HP10 9SD

FURTHER INFORMATION

If you need any more information, our knowledgeable Customer Services Team are on hand to answer any of your questions about product characteristics or use of products.

Alternatively, you can contact your local Account Manager directly to take advantage of their knowledge and training.

ORDERING INFORMATION AND GUIDELINES

Please contact our Customer Services Team for information about minimum order values. We aim to dispatch items within three working days if in stock.



THE KEY SURGICAL STORY

Key Surgical was founded in 1988 by three entrepreneurial women with years of experience in both the nursing profession and the medical supply industry. Motivated by the growing need for a single source supplier that would provide excellent Customer service, the partners began building an extensive product line based on market need. We are dedicated to meeting the healthcare industry's highest quality standards for patient care, safety, and efficiency in hospitals, clinics, and surgical centres.

Over 30 years later, the same Customer-driven principals and patient-focused values upon which Key Surgical was founded continue to guide everything we do. With over 3,000 of the highest quality products available to Customers worldwide, Key Surgical is poised to continue to be a compassionate industry leader for years to come.

OUR MISSION

Our mission at Key Surgical is three-fold. First, to hire and retain good people who understand that nothing happens until a Customer decides to buy something from us. Second, to grow by providing the quality products, value and service that our Customers want, when they want them. Finally, to be good citizens by providing a rewarding and fun place for employees to work and by sharing the benefits of our growth with the community in which we operate.

A NOTE ABOUT QUALITY

There are no shades of grey in the medical marketplace when it comes to quality. Key Surgical is committed to manufacturing and distributing the highest quality medical products, and to maintaining the highest level of compliance to FDA and ISO regulations. Nobody can claim to be perfect, but we strive to continuously improve our processes and systems so that we will always be closer to perfect tomorrow than we are today.

LETTER TO OUR CUSTOMERS

Thank you for being our Customer and for choosing us to be your trusted supplier for products and solutions for your sterile processing department. Our Customers are at the forefront of everything we do, we are dedicated to meeting the healthcare industry's highest quality standards for patient care, safety, and efficiency in hospitals, clinics, and surgical centres around the world. As your partner, we promise to continue providing a comprehensive selection of products and personalised service you won't find anywhere else.

Our Customer base of hospitals, theatre centres, and distributors continues to expand throughout the UK and in 70 other countries around the world. We're committed to continuously making investments in our products, quality system, regulatory and compliance efforts, and operations with the goal of being the easiest company with which to do business.

Although seemingly simple, focusing on delivering an exceptional experience with us every step of the way is, and will continue to be, the foundation of why we do what we do. Our entire team works hard every day to be recognised by you as a company of dedication, honesty, and integrity. In striving toward that recognition, we promise to keep our feet planted firmly on the ground and to make sure we brighten both your day and ours in the process.

If we can ever be of assistance, please call or email us. Businesses grow because they do the right things for their Customers and their employees. We intend to keep doing those things and to keep growing, and we truly appreciate your business and thank you for your continued support.

Sincerely,



Arnaud Corsat
Managing Director of EMEA and Exports



Ross Tustain
Managing Director, UK













TABLE OF CONTENTS

-  00 – Instrument Cleaning and Maintenance
-  48 – Sterilisation Packaging
-  80 – Trays and Accessories
-  128 – Containers and Accessories
-  134 – Instrument Identification
-  148 – Instrument Protection
-  178 – Tags
-  212 – Transport and Disposal
-  224– Index



PLEASE NOTE: Unless specified, all products are not made with natural rubber latex.

INSTRUMENT CLEANING AND MAINTENANCE

The Key Surgical Instrument Cleaning line is designed to provide you with the tools needed to clean your surgical instruments. In addition to our wide variety of cleaning brushes, we also offer enzymatic solutions to help improve the decontamination process.

With bristles available in nylon, brass, polypropylene or stainless steel, we also provide an incredible variety of lengths and diameters to fit your needs.

Key Surgical Endoscope Brushes are available in a wide range of diameters and lengths designed for manually cleaning the channels of a variety of endoscopes such as colonoscopes, gastroscopes, ureteroscopes, and duodenoscopes. Our selection of brushes to clean valves, suction wells, or elevator mechanisms as well as brush kits will have your department set up for success.

With our vast assortment of brushes - from channel cleaning to endoscope cleaning to steriliser cleaning - we can help you match the right brush to the right cleaning task at hand.

Channel Cleaning Brushes

Pipe Cleaners

Toothbrush-Style Cleaning Brushes

General Instrument Cleaning Brushes


Steriliser Cleaning Brushes

Brush Assortments

Endoscope Cleaning Brushes

Please note: Unless specified, all products are not made with natural rubber latex.





BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.


	Product code	Diameter	Bristle area length	Length	Qty
A	BR-06-039	1mm	6.35mm	152mm	2 pcs
B	BR-08-039	1mm	12.7mm	203.2mm	2 pcs
	BR-18-039	1mm	6.35mm	457mm	2 pcs
C	BR-06-062	1.58mm	12.7mm	152mm	2 pcs
	BR-18-062	1.58mm	12.7mm	457mm	2 pcs
D	BR-08-070	1.77mm	12.7mm	203.2mm	2 pcs
	BR-18-070	1.77mm	12.7mm	457mm	2 pcs
	BR-24-070	1.77mm	12.7mm	609.6mm	2 pcs
E	BR-06-079	2mm	12.7mm	152mm	2 pcs
F	BR-08-079	2mm	19.05mm	203.2mm	2 pcs
	BR-12-079	2mm	19.05mm	305mm	2 pcs
G	BR-24-080	2mm	25.4mm	609.6mm	2 pcs
H	BR-08-087	2.2mm	12.7mm	203.2mm	2 pcs
	BR-18-087	2.2mm	12.7mm	457mm	2 pcs
I	BR-08-092	2.3mm	12.7mm	203.2mm	2 pcs
J	BR-24-098	2.49mm	12.7mm	609.6mm	2 pcs
K	BR-08-105	2.67mm	25.4mm	203.2mm	2 pcs
	BR-18-105	2.67mm	25.4mm	457mm	2 pcs
L	BR-08-118	3mm	38.1mm	203.2mm	2 pcs
M	BR-12-118	3mm	69.85mm	305mm	2 pcs

Instrument Cleaning – Channel Cleaning Brushes


CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-18-125	3.18mm	50.8mm	457mm	2 pcs
B	BR-08-131	3.3mm	19.05mm	203.2mm	2 pcs
C	BR-08-144	3.66mm	19.05mm	203.2mm	2 pcs
D	BR-08-158	4mm	19.05mm	203.2mm	2 pcs
	BR-12-158	4mm	19.05mm	305mm	2 pcs
	BR-18-158	4mm	19.05mm	457mm	2 pcs
E	BR-12-197	5mm	76.2mm	305mm	2 pcs
	BR-16-197	5mm	76.2mm	406mm	2 pcs
	BR-18-197	5mm	76.2mm	457mm	2 pcs
	BR-24-197	5mm	76.2mm	609.6mm	2 pcs
F	BR-24-200	5.08mm	50.8mm	609.6mm	2 pcs
G	BR-12-236	6mm	50.8mm	305mm	2 pcs
	BR-16-236	6mm	50.8mm	406mm	2 pcs
	BR-24-236	6mm	50.8mm	609.6mm	2 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE



BRUSH TIP DETAILS SHOWN ACTUAL SIZE


Instrument Cleaning – Channel Cleaning Brushes

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-06-250	6.35mm	50.8mm	152mm	2 pcs
	BR-08-250	6.35mm	50.8mm	203.2mm	2 pcs
	BR-18-250	6.35mm	50.8mm	457mm	2 pcs
B	BR-16-276	7mm	76.2mm	406mm	2 pcs
	BR-24-276	7mm	76.2mm	609.6mm	2 pcs
C	BR-24-312	7.9mm	76.2mm	609.6mm	2 pcs
D	BR-06-375	9.5mm	38.1mm	152mm	2 pcs
E	BR-24-375	9.5mm	76.2mm	609.6mm	2 pcs
F	BR-16-393	10mm	76.2mm	406mm	2 pcs
	BR-18-393	10mm	76.2mm	457mm	2 pcs
	BR-24-393	10mm	76.2mm	609.6mm	2 pcs
G	BR-24-425	10.8mm	76.2mm	609.6mm	2 pcs

Instrument Cleaning – Channel Cleaning Brushes



CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-08-440	11.18mm	76.2mm	203.2mm	2 pcs
	BR-24-440	11.18mm	76.2mm	609.6mm	2 pcs
B	BR-16-443	11.25mm	76.2mm	406mm	2 pcs
	BR-24-443	11.25mm	76.2mm	609.6mm	2 pcs
C	BR-16-472	12mm	76.2mm	406mm	2 pcs
D	BR-06-500	12.7mm	50.8mm	152mm	2 pcs
E	BR-16-591	15mm	50.8mm	406mm	2 pcs
	BR-18-591	15mm	50.8mm	457mm	2 pcs

BRUSH TIP DETAILS SHOWN ACTUAL SIZE



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-08-630	16mm	50.8mm	203mm	2 pcs
B	BR-16-787	20mm	50.8mm	406mm	2 pcs
C	BR-18-880	22.35mm	88.9mm	457mm	2 pcs
D	BR-24-880	22.35mm	107.95mm	609.6mm	2 pcs

K

Instrument Cleaning

Instrument Cleaning – Channel Cleaning Brushes

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end (pull through style excluded).

	Product code	Description	Diameter	Bristle area length	Length	Qty
A	BR-3044	nylon	0.5mm	6.35mm	76mm	10 pcs
B	BR-3045	nylon	0.5mm	6.35mm	203.2mm	10 pcs
C	BR-3047	nylon	0.9mm	12.7mm	101.6mm	10 pcs
D	BR-4120	nylon, centered (pull through)	6.35mm to 15.24mm	50.8mm	457.2mm	2 pcs
E	BR-4140	nylon, centered (pull through)	12.7mm	76.2mm	609.6mm	6 pcs



A



B



C



D

(PULL THROUGH)



E

(PULL THROUGH)


BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a twisted stainless steel handle with loop end. Quantities also available in 50 pcs.

	Product code	Description	Diameter	Bristle area length	Length	Qty
A	BR-8110	curved bristle	10mm	44.45mm	116.84mm	2 pcs
B	BR-8115	curved bristle	13.99mm	92.075mm	203.2mm	2 pcs
C	BR-8210	tapered fan bristle	4.75mm to 12.7mm	76.2mm	177.8mm	2 pcs
D	BR-8310	straight bristle	12.7mm	63.5mm	304.8mm	3 pcs
E	BR-8315	straight bristle	12.7mm	76.2mm	304.8mm	3 pcs
F	BR-8410	fan bristle	12.7mm	76.2mm	152.4mm	3 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

BONE REAMER BRUSHES

Channel cleaning brushes featuring antimicrobial nylon bristles on a rigid stainless steel handle with loop end.

	Product code	Diameter	Bristle area length	Length	Qty
A	BR-31-118	3mm	76mm	787mm	3 pcs
B	BR-31-197	5mm	76mm	787mm	3 pcs
C	BR-31-276	7mm	76mm	787mm	3 pcs




BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes

POLYPROPYLENE BRUSHES

Channel cleaning brushes featuring rigid polypropylene bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	PP-16-039	1mm	25.4mm	406mm	2 pcs
	PP-22-039	1mm	25.4mm	558.8mm	2 pcs
B	PP-16-062	1.58mm	25.4mm	406mm	2 pcs
	PP-22-062	1.58mm	25.4mm	558.8mm	2 pcs
C	PP-16-087	2.2mm	25.4mm	406mm	2 pcs
	PP-22-087	2.2mm	25.4mm	558.8mm	2 pcs
D	PP-16-105	2.67mm	50.8mm	406mm	2 pcs
	PP-22-105	2.67mm	50.8mm	558.8mm	2 pcs
E	PP-16-125	3.18mm	50.8mm	406mm	2 pcs
	PP-22-125	3.18mm	50.8mm	558.8mm	2 pcs
F	PP-16-158	4mm	50.8mm	406mm	2 pcs
	PP-22-158	4mm	50.8mm	558.8mm	2 pcs
G	PP-18-197	5mm	50.8mm	457mm	2 pcs
	PP-22-197	5mm	50.8mm	558.8mm	2 pcs




BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Channel Cleaning Brushes


POLYPROPYLENE BRUSHES

Channel cleaning brushes featuring rigid polypropylene bristles on a twisted stainless steel handle with loop end. Quantities available 2 pcs and 50 pcs.

	Product code	Diameter	Bristle area length	Length	Qty
A	PP-18-250	6.35mm	50.8mm	457mm	2 pcs
	PP-24-250	6.35mm	50.8mm	610mm	2 pcs
B	PP-18-393	10mm	50.8mm	457.2mm	2 pcs
C	PP-18-591	15mm	50.8mm	457.2mm	2 pcs
D	PP-18-750	19.05mm	152.4mm	457.2mm	2 pcs
E	PT-28-118 (pull through)	3mm	50.8mm	711mm	2 pcs
F	PT-28-197 (pull through)	5mm	50.8mm	711mm	2 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE



Instrument Cleaning – Pipe Cleaners

PIPE CLEANERS - REEL

Available in a range of different materials. Rolls allow you to cut to the size you need, minimising waste. Single use.

	Product code	Hardness	Material	Diameter	Length	Qty
A	26322	soft	cotton	3mm	10000mm	1 reel
B	26323	medium	cotton/nylon	3mm	10000mm	1 reel
C	26328	medium	polyester/nylon	15mm	5000mm	1 reel
D	26329	medium	polyester/nylon	15mm	5000mm	1 reel
E	26334	hard	polyester/nylon	3mm	7500mm	1 reel
F	26335	hard	polyester/nylon	6mm	7500mm	1 reel
G	26336	hard	polyester/nylon	3mm	7500mm	1 reel
H	26337	medium	polyester/nylon	6mm	7500mm	1 reel
I	PC-300P	N/A	polyester	3.18mm	15240mm	1 reel
J	PC-400P	N/A	polyester	6.35mm	15240mm	1 reel


PIPE CLEANERS - LOOSE

Available in either polyester (non-woven) or polyester and nylon. Single use.

	Product code	Hardness	Material	Diameter	Length	Qty
K	26324	soft/hard	polyester/nylon	3mm	150mm	100 pcs
L	26325	soft	polyester (non-woven)	9mm	300mm	50 pcs
M	26326	soft	polyester (non-woven)	12mm	300mm	50 pcs
N	26327	medium	polyester/nylon	15mm	300mm	20 pcs
O	26342	soft	polyester (non-woven)	3mm	300mm	50 pcs

K

Instrument Cleaning



Instrument Cleaning – Pipe Cleaners

FLEXISTEM AND MICROCLEAN

Flexistem and Microclene are unique cleaners for all small-bore cannula. Various diameters from 3mm to 15mm are available constructed from very soft to hard bristles inserted between two soft mild steel twisted wires which assist in following contours while resisting buckling. Single use. Flexistem FLEX007A and FLEX007B has a wavy bristle configuration to reduce buckling which can occur in long instruments.

FLEXISTEM

	Product code	Type	Material	Diameter	Length	Qty
A	FLEX001	soft cleaner	cotton	3mm	10000mm	1 reel
B	FLEX002	super bristle	polyester/nylon	3mm	10000mm	1 reel
	FLEX003	double ended	polyester/nylon	3mm	300mm	100 pcs
	FLEX004	fleece	polyester	9mm	300mm	50 pcs
C	FLEX005	fleece	polyester	12mm	300mm	50 pcs
D	FLEX006	bristle stem	polyester	15mm	300mm	20 pcs
E	FLEX007	bristle stem	polyester/nylon	15mm	5000mm	1 reel
	FLEX007A	bristle stem	polyester/nylon	15mm	5000mm	1 reel
	FLEX007B	bristle stem	polyester/nylon	N/A	66.5mm	10 pcs
F	FLEX033	fleece	polyester	3mm	300mm	50 pcs

MICROCLEAN

	Product code	Type	Material	Diameter	Length	Qty
G	FLEX021	super fine	polyester/nylon	3mm	7500mm	1 reel
H	FLEX022	medium fine	polyester/nylon	6mm	7500mm	1 reel
I	FLEX023	super brush	polyester/nylon	3mm	7500mm	1 reel
J	FLEX024	medium brush	polyester/nylon	6mm	7500mm	1 reel

Instrument Cleaning – Toothbrush-Style Cleaning Brushes

TOOTHBRUSH-STYLE CLEANING BRUSHES

Antimicrobial nylon, brass, or stainless bristles on a sturdy, plastic handle. Designed for external cleaning.

	Product code	Description	Total length	Qty
A	09320	extra strong synthetic bristles	165mm	10 pcs
B	09325	stainless steel bristles	180mm	1 pc
C	09485	extra strong synthetic bristles, angled handle	235mm	2 pcs
D	09478	soft synthetic bristles	215mm	1 pc
E	09488	synthetic bristles	180mm	1 pc
F	09489	brass bristles	180mm	1 pc
G	26491	soft synthetic bristles	215mm	1 pc
H	26492	hard synthetic bristles	215mm	1 pc
I	26029	nylon bristles, with sliding cover	55mm	3 pcs



A



B



C



D



E




F



G



H



I



A

B

C

D

E

F

G

H

I

Instrument Cleaning – Toothbrush-Style Cleaning Brushes

TOOTHBRUSH-STYLE CLEANING BRUSHES

Antimicrobial nylon, brass, or stainless bristles on a sturdy, plastic handle. Designed for external cleaning.

	Product code	Description	Total length	Qty
A	09084	soft bristles, double-ended	175mm	2 pcs
B	09459	extra strong synthetic bristles, double-ended	175mm	2 pcs
C	09098	nylon bristles, double-ended	175mm	2 pcs
D	09460	extra strong, synthetic bristles	215mm	3 pcs
E	B-2100	brass bristles, three rows	178mm	3 pcs
F	N-2121	nylon bristles, double-ended, three rows/one row	178mm	2 / 50 pcs
G	N-2221	nylon bristles, double-ended, three rows/one row	175mm	2 / 50pcs
H	S-2121	stainless steel bristles, double-ended, three rows/one row	178mm	2 pcs
I	S-2100	stainless steel bristles, three rows	178mm	3 pcs
J	D-2100BN	brass, nylon bristles, three rows	178mm	3 pcs
K	D-2100SN	stainless steel bristles, nylon, three rows	178mm	3 pcs



K

Instrument Cleaning



A



B



C



D



E



F

Instrument Cleaning – Toothbrush-Style and General Instrument Cleaning Brushes

TOOTHBRUSH-STYLE CLEANING BRUSHES

Antimicrobial nylon, brass, or stainless bristles on a sturdy, plastic handle. Designed for external cleaning.

	Product code	Description	Total length	Qty
A	N-2100	nylon bristles, three rows	178mm	3 pcs
B	N-2200	extra rigid nylon bristles, angled handle	235mm	2 pcs
C	S-2300	soft stainless steel bristles, three rows	178mm	3 / 50 pcs
D	S-2131	soft stainless steel bristles, one row	178mm	3 / 50 pcs
	79325	stainless steel bristles	180mm	1 pc
	96281.3	stainless steel bristles	180mm	3 pcs
	96281.4	extra strong synthetic bristles, angled handle	235mm	3 pcs

GENERAL INSTRUMENT CLEANING BRUSHES

A range of instrument cleaning brushes with strong durable nylon bristles and handles.

	Product code	Description	Total length	Qty
E	BBR001	nylon bristles	270mm	1 pc
F	BBR002	nylon bristles	395mm	1 pc
G	BBR003	nylon bristles	184mm	1 pc



G



Instrument Cleaning – General Instrument and Steriliser Cleaning Brushes

GENERAL INSTRUMENT CLEANING BRUSHES

Antimicrobial nylon or stainless steel bristles on a sturdy, plastic handle. Designed for external cleaning.

	Product code	Description	Length	Qty
A	BR-2300	stainless steel	152.4mm	2 pcs
B	BR-04-003	stainless steel	91.44mm	2 pcs
C	09099	stainless steel with nylon bristles	450mm	2 pcs

STERILISER CLEANING BRUSH

Brush allows for greater reach and cleaning ability. Bristles are at a forward angle to loosen soil without marring surface.

	Product code	Description	Bristle area	Length	Qty
D	BR-5000	double-surface, polyester bristles	152.4mm x 127mm	787.4mm	1 pc
E	BR-5200	threaded handle, plastic	N/A	914.4mm	1 pc

PISTOL-GRIP BRUSH

Soak-proof design with pistol-grip handle. Crimped nylon bristles on a block handle that resists chipping and cracking.

	Product code	Description	Bristle length	Bristle area	Length	Qty
F	BR-4000	crimped nylon bristle, plastic pistol grip	50.8mm	101.6mm x 101.6mm	203.2mm	1 pc


STERILISER DRAIN CLEANING BRUSH

White nylon bristles designed to clean out drains, lines and openings.


	Product code	Diameter	Diameter	Length	Qty
G	BR-5500	long handed brush, nylon bristles, with plastic handle	34.8mm	711.2mm	1 pc
H	BR-5600	drainline brush, nylon bristles, with stainless steel handle	50.8mm	609.6mm	1 pc




A




C




B



D



E



Instrument Cleaning – Brush Assortments

ORIGINAL BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- A** BR-4199 BR-18-062, BR-18-087, BR-18-105, BR-18-125, BR-18-250, FT-18-500, N-2100, N-2121, S-2100

SMALL DIAMETER CHANNEL BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- B** BR-4299 BR-08-158, BR-12-079, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-24-197, BR-24-200

LARGE DIAMETER CHANNEL BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- C** BR-4399 BR-06-375, BR-08-630, BR-16-276, BR-16-393, BR-18-250, BR-18-591, BR-18-880, BR-24-312, BR-24-440

POPULAR CHANNEL BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- D** BR-4499 BR-06-375, BR-08-158, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-18-250, BR-24-200

MICRO-OPHTHALMIC BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- E** BR-4599 BR-18-039, BR-18-062, BR-18-070, BR-18-087, BR-18-105, BR-18-125, BR-18-158, N-2100


STERILISER CLEANING BRUSH ASSORTMENT

Product code Includes one of each of the following brushes

- BR-6099 BR-5000, BR-5200, BR-5400, BR-5500, BR-5600

K

Instrument Cleaning



Instrument Cleaning – Endoscope Cleaning Brushes

ENDOSCOPE BRUSHES


Nylon bristles attached to flexible, coloured plastic tubing. Single-ended. Single use. Available in 1 pc or 50 pcs.

	Product code	Colour	Brush diameter	Bristle area length	Length	Qty
A	BR-075-18	blue	1.8mm	20mm	750mm	1 pc
B	BR-075-24	yellow	2.4mm	20mm	750mm	1 pc
C	BR-075-32	green	3.2mm	20mm	750mm	1 pc
D	BR-075-42	pink	4.2mm	20mm	750mm	1 pc
E	BR-100-18	blue	1.8mm	20mm	1000mm	1 pc
F	BR-100-24	yellow	2.4mm	20mm	1000mm	1 pc
G	BR-100-32	green	3.2mm	20mm	1000mm	1 pc
H	BR-100-42	pink	4.2mm	20mm	1000mm	1 pc
I	BR-150-18	blue	1.8mm	20mm	1500mm	1 pc
J	BR-150-24	yellow	2.4mm	20mm	1500mm	1 pc
K	BR-150-32	green	3.2mm	20mm	1500mm	1 pc
L	BR-150-42	pink	4.2mm	20mm	1500mm	1 pc
M	BR-200-18	blue	1.8mm	20mm	2000mm	1 pc
N	BR-200-24	yellow	2.4mm	20mm	2000mm	1 pc
O	BR-200-32	green	3.2mm	20mm	2000mm	1 pc
P	BR-200-42	pink	4.2mm	20mm	2000mm	1 pc


ENDOSCOPE BRUSHES

Tapered cleaning brushes. Nylon bristles attached to flexible plastic tubing. Single-ended. Single use. Available in 1 pc and 50 pcs.


	Product code	Brush diameter	Bristle area length	Length	Qty
Q	BR-3308	tapered 3mm to 5mm	10mm	2400mm	1 pc
R	BR-3309	tapered 3mm to 5mm	10mm	270mm	1 pc
S	BR-3310	tapered 3mm to 5mm	20mm	2400mm	1 pc



A




D




B



E



C



F


BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning – Endoscope Cleaning Brushes

ENDOSCOPE BRUSHES




A variety of single, double, and triple bristled endoscope brushes as well as valve brushes and kits.

	Product code	Description	Brush diameter	Bristle area length	Length	Qty
A	EB10100	double-ended	10mm & 5mm	30mm & 20mm	2500mm	50 pcs
B	EB10102	double-ended	3mm	20mm	2300mm	50 pcs
C	EB10103	double-ended	6mm	20mm	2500mm	50 pcs
D	EB10103.3	triple-cleaning	6mm	20mm	2500mm	50 pcs
E	EB10104	double-ended with handle	10mm & 5mm	40mm & 20mm	2300mm	25 pcs
F	EB10105	double-ended with pipe cleaner	5mm	125mm & 20mm	2300mm	50 pcs
G	EB10106	channel brush with handle	2mm	10mm	1200mm	50 pcs
H	EB10136	double-ended	3mm	20mm	1200mm	50 pcs



K

Instrument Cleaning



BRUSH TIP DETAILS
SHOWN ACTUAL SIZE

Instrument Cleaning – Endoscope Cleaning Brushes

ENDOSCOPE BRUSHES


A variety of endoscope brushes conveniently packaged with a valve brush. Single use.


Product code	Description	Qty
A EB10103.3C	includes triple-cleaning channel EB10103.3 and valve brush EB10101	50 pcs
B EB10107C	includes double-ended channel EB10103 and valve brush EB10101	50 pcs
C EB10107G	includes double-ended channel EB10102 and valve brush EB10101	50 pcs

ENDOSCOPE BRUSHES (TAPERED, SINGLE-ENDED)

Nylon bristles attached to flexible plastic tubing. Single-ended. Single use.

Product code	Description	Bristle area length	Length	Qty
BR-3711-50	tapered 4mm to 6mm	10.5mm	2400mm	50 pcs
D BR-3799-50	BR-3711-50 and BR-3201	N/A	2400mm	50 pcs





Instrument Cleaning – Endoscope Cleaning Brushes

ENDOSCOPE BRUSHES (TAPERED, DOUBLE-ENDED)

Nylon bristles attached to flexible plastic tubing. Double-ended. Single use. Available in 1 pc or 50 pcs.

	Product code	Brush diameter	Bristle area length	Length	Qty
A	BR-3320	tapered 3mm to 5mm	10mm	2400mm	1 pc
B	BR-3330	tapered 3mm to 5mm	20mm	2400mm	1 pc

ENDOSCOPE BRUSHES (NON-TAPERED, DOUBLE-ENDED)

Nylon bristles attached to flexible plastic tubing. 3mm brush on one end, 6mm brush on the other. Single use. Available in 1 pc or 50 pcs.

	Product code	Brush diameter	Bristle area length	Length	Qty
C	BR-3340	3mm & 6mm	20mm	2400mm	1 pc

PULL-THROUGH ENDOSCOPE BRUSHES

Delicate pull-through cleaning brushes. Nylon bristles attached to twisted stainless steel wire. Single-ended. Single use. Available in 1 pc or 50 pcs.

	Product code	Brush diameter	Bristle area length	Length	Qty
D	BR-3030	1.1mm	12.7mm	650mm	1 pc
	BR-3035	1.1mm	10.0mm	2300mm	1 pc
E	BR-3040	2.0mm	12.7mm	1200mm	1 pc
F	BR-3410	3.0mm	13.0mm	1200mm	1 pc
G	BR-3420	3.0mm	13.0mm	1600mm	1 pc

K

Instrument Cleaning

Instrument Cleaning – Endoscope Cleaning Brushes

A



B



C



D



E



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

ALBARRAN BRUSH

Delicate brush designed for cleaning elevator channel on a duodenoscope. Single use.

Product code	Description	Brush diameter	Bristle area length	Qty
A EB10108	albarran entrance brush	2.3mm	13mm	100 pcs


VALVE BRUSH

Short brush designed for cleaning valve ports and suction wells on an endoscope. Single use.

Product code	Description	Brush diameter	Bristle area length	Qty
B BR-3201	suction well brush	11mm	45mm	3 / 50 pcs
C BR-3202	double-ended	4.3mm & 11mm	10.2mm & 45mm	3 / 50 pcs
D EB10101	double-ended	10mm & 5mm	30mm & 20mm	100 pcs
E BR-3307	finger loop style	11mm	45mm	1 pc

DID YOU KNOW?

Key Surgical is a global leading partner to hospitals and theatres or procedural facilities, supporting the needs in endoscopy, operating room, and instrument reprocessing. Our comprehensive portfolio is designed with a focus on patient safety; allowing our Customers to efficiently source the high quality products they need with the best Customer service and overall experience with every interaction. Our UK office is based in High Wycombe but our main headquarters are in the US and our products are distributed in over 70 countries around the world. Throughout our continued growth we have remained true to our mission; to be the easiest company in our industry with which to do business.



Instrument Cleaning – Endoscope Cleaning Brushes

PREMIUM ENDOSCOPE BRUSH

Nylon bristles on twisted stainless steel wire or spring coil for hard-to-turn corners. Loop handle at one end and stainless steel beaded tip at bristle end. Reusable.

	Product code	Description	Brush diameter	Bristle area length	Qty
A	BR-3300	spring coil	tapered 4mm to 6mm	10mm	1 pc
B	BR-3305	spring coil	tapered 3mm to 5mm	10mm	1 pc
C	BR-3702	spring coil	6mm	5mm	1 pc
D	BR-3701	twisted wire	2.4mm	6mm	1 pc
E	BR-3703	twisted wire	2.4mm	6mm	1 pc

ENDOSCOPIC BRUSH ASSORTMENTS

Product code	Includes one of each of the following brushes
BR-4799	BR-12-079, BR-12-118, BR-12-197, BR-16-197, BR-16-393, BR-24-197, BR-24-393, BR-3202, N-2121
F	BR-100-99 BR-100-18, BR-100-24, BR-100-32, BR-100-42
G	BR-150-99 BR-150-18, BR-150-24, BR-150-32, BR-150-42
H	BR-200-99 BR-200-18, BR-200-24, BR-200-32, BR-200-42

K

Instrument Cleaning

STERILISATION PACKAGING

Our flexible production methods ensure that the best sizes of wrap can be produced for our Customers, even if these are not industry standard sizes, without compromising performance or protection levels.

Our Barrier Wraps provide a bacterial barrier and protect the contents of a set or pack from contamination. When classed as alcohol repellent, barrier wraps can be used within the sterile field.

Choice Wrap is designed to reduce the risk of rips, tears, and wet packs. Choice SMMS are made of polypropylene and come in three weights: 43gsm, 55gsm, 70gsm. Any two of the Choice Wraps can be interleaved in blue or green to meet the sterilisation service user requirements.

Want to trust your Transport Wrap? All Transport Wraps must be used in conjunction with one compliant Barrier Wrap to completely prevent the ingress of micro-organisms. Use of a single barrier layer with Transport Wrap can be a cost-effective solution that reduces damage in handling and transport. If adopted, changes to both wrapping and unwrapping procedures may be required to ensure aseptic opening.

Sterilisation Pouches and Reels are also available to ensure that the wrapping process is efficient and effective. Made to withstand vigorous processing and handling, pouches and reels are made from medical grade film with breathable webs and suit all common sterilisation methods.

Barrier Wrap

Crepe Wrap

Liners

Sterilisation Tape

Transport Wrap

Ancor Pouches and Reels

Wipak Pouches and Reels

Please note: Unless specified, all products are not made with natural rubber latex.





Sterilisation Packaging – Barrier Wrap

CLINIPAK CHOICE WRAP

A fully synthetic SMMS product designed to offer flexibility and enable users to pack conventionally or with 'two layers' at once. Clinipak Choice Wrap is strong and available with alcohol repellence in both green and blue for use within the sterile field. Wraps are also available with blue and green interleaved. Compatible with steam and EtO sterilisation. Single use.

A CLINIPAK CHOICE WRAP BLUE

Product code	Size	GSM	Qty
CB060	60 x 60cm	55gsm	350 sheets
CB075	75 x 75cm	55gsm	300 sheets
CB090	90 x 90cm	55gsm	300 sheets
CB100	100 x 100cm	55gsm	200 sheets
CB120	120 x 120cm	55gsm	120 sheets
CB1315	130 x 150cm	55gsm	100 sheets
CB1715	170 x 150cm	55gsm	50 sheets
CB1815	180 x 150cm	55gsm	50 sheets

B CLINIPAK CHOICE GREEN WRAP

Product code	Size	GSM	Qty
CG060	60 x 60cm	55gsm	350 sheets
CG075	75 x 75cm	55gsm	300 sheets
CG090	90 x 90cm	55gsm	300 sheets
CG100	100 x 100cm	55gsm	200 sheets

C CLINIPAK CHOICE BLUE LIGHT WRAP

Product code	Size	GSM	Qty
CBA050	50 x 50cm	43gsm	300 sheets
CBA060	60 x 60cm	43gsm	350 sheets
CBA090	90 x 90cm	43gsm	300 sheets
CBA100	100 x 100cm	43gsm	200 sheets
CBA120	120 x 120cm	43gsm	120 sheets
CBA1415	140 x 150cm	43gsm	100 sheets





Sterilisation Packaging – Barrier Wrap

A CLINIPAK CHOICE BLUE ULTRA WRAP

Clinipak Choice Blue Ultra is extra strong, drapeable, and available with alcohol repellence in both green and blue for use within the sterile field. Compatible with steam and EtO sterilisation. Single use.

Product code	Size	GSM	Qty
CBU1014	100 x 140cm	70gsm	120 sheets
CBU120	120 x 120cm	70gsm	120 sheets
CBU1315	130 x 150cm	70gsm	100 sheets
CBU140	140 x 140cm	70gsm	100 sheets

B CLINIPAK CHOICE GREEN X WRAP

A fully synthetic SMMS product designed to offer flexibility and enable users to pack conventionally or with 'two layers' at once. Available with alcohol non repellent in green for use within the sterile field. Compatible with steam and EtO sterilisation. Single use.

Product code	Size	GSM	Qty
CGX090	90 x 90cm	55gsm	300 sheets
CGX100	100 x 100cm	55gsm	200 sheets





Sterilisation Packaging – Barrier Wrap

CLINIPAK CHOICE WRAP DRY

A patented tray wrap solution is designed to heavily reduce the occurrence of wet packs. The inner surface exhibits a high degree of absorbency and the outer surface repels both water and alcohol. The absorbent surface disperses the condensation within a typical steam sterilisation process and assists with evaporation during the vacuum cycle. This wrap is manufactured with a blend of materials that are compliant with the recognised standards for toxicology, bacterial resistance, and drape. Compatible with steam and EtO sterilisation. Single use.

CLINIPAK CHOICE WRAP DRY GREEN

Product code	Size	GSM	Qty
CDG075	75 x 75cm	79gsm	300 sheets
CDG090	90 x 90cm	79gsm	300 sheets
CDG100	100 x 100cm	79gsm	200 sheets
CDG120	120 x 120cm	79gsm	120 sheets
CDG140	140 x 140cm	79gsm	100 sheets
CDG150	150 x 150cm	79gsm	100 sheets
CDG1315	130 x 150cm	79gsm	100 sheets





Sterilisation Packaging – Barrier Wrap

A RELIANCE® SOLO WRAP

A two-sheet wrap system, with two colours of latex-free wrap joined together along the sides. Each sheet is a five-layer nonwoven tray wrap, manufactured from spun-bonded and melt blown polypropylene, providing a wet and dry microbiological barrier. Bonding the wrap simplifies the wrapping process whilst increasing durability to avoid rips, tears or punctures. Wraps are also available interleaved with other compliant barrier layers, or transportation wrap. Blue and green bonded wrap. Compatible with steam and EtO sterilisation. Single use.

RELIANCE® SOLO WRAP 200B/G

Product code	Size	GSM	Qty
CT43LGT-60	60 x 60cm	43gsm	125 sets
CT43LGT-75	75 x 75cm	43gsm	150 sets
CT43LGT-90	90 x 90cm	43gsm	50 sets
CT43LGT-100	100 x 100cm	43gsm	50 sets
CT43LGT-120	120 x 120cm	43gsm	50 sets
CT43LGT-140	140 x 140cm	43gsm	25 sets
CT43LGT-1418	140 x 180cm	43gsm	25 sets

RELIANCE® SOLO WRAP 300B/G

Product code	Size	GSM	Qty
CT50MED-60	60 x 60cm	50gsm	125 sets
CT50MED75U	75 x 75cm	50gsm	50 sets
CT50MED90U	90 x 90cm	50gsm	50 sets
CT50MED-100	100 x 100cm	50gsm	50 sets
CT50MED120U	120 x 120cm	50gsm	50 sets
CT50MED140U	140 x 140cm	50gsm	25 sets
CT50MED-1418	140 x 180cm	50gsm	25 sets

RELIANCE® SOLO WRAP 400B/G

Product code	Size	GSM	Qty
CT60HVV-60	60 x 60cm	60gsm	125 sets
CT60HVV-75	75 x 75cm	60gsm	50 sets
CT60HVV-90	90 x 90cm	60gsm	50 sets
CT60HVV-100	100 x 100cm	60gsm	50 sets
CT60HVV120U	120 x 120cm	60gsm	25 sets
CT60HVV140U	140 x 140cm	60gsm	25 sets
CT60HVV-1418	140 x 180cm	60gsm	25 sets





Sterilisation Packaging – Crepe Wrap

RELIANCE® NON-WOVEN WRAP

Non-Woven Wrap is a high-performance sterilisation material that is suitable for inner and outer wrapping of medium to large trays and packs. This is a wet laid wrap with 50% cellulose and 50% synthetic fibres that improve strength and increase softness. These products are also available in interleaved green and blue wrap. Single use.

	Product code	Colour	Size	GSM	Qty
A	RELBLU090	blue	90 x 90cm	60gsm	200 sheets
	RELGRN100U	blue	100 x 100cm	60gsm	200 sheets
	RELBLU120	blue	120 x 120cm	60gsm	100 sheets
B	RELBLU100U	green	100 x 100cm	60gsm	200 sheets
	RELGRN1315U	green	130 x 150cm	60gsm	100 sheets

RELIANCE® CREPE WRAP

Crepe paper is suitable for lighter instruments and sets. Manufactured from 100% medical quality cellulose fibres and combines most economical wrapping with good drapability. A soft crepe and easier to handle due to the creping process. Compatible with steam and EtO sterilisation. Single use.

	Product code	Colour	Size	GSM	Qty
C	GRNCREPE025	green	25 x 25cm	60gsm	2000 sheets
	GRNCREPE060N	green	60 x 60cm	60gsm	500 sheets
	GRNCREPE075NU	green	75 x 75cm	60gsm	250 sheets
	GRNCREPE090U	green	90 x 90cm	60gsm	250 sheets
	GRNCREPE100U	green	100 x 100cm	60gsm	250 sheets
	GRNCREPE1014	green	100 x 140cm	60gsm	125 sheets
	GRNCREPE120U	green	120 x 120cm	60gsm	125 sheets
	GRNCREPE1315U	green	130 x 150cm	60gsm	125 sheets
D	WHCREPE030	white	30 x 30cm	60gsm	2000 sheets
	WHCREPE050N	white	50 x 50cm	60gsm	500 sheets
	WHCREPE060N	white	60 x 60cm	60gsm	500 sheets
	WHCREPE075NU	white	75 x 75cm	60gsm	250 sheets
	WHCREPE090U	white	90 x 90cm	60gsm	250 sheets
	WHCREPE100U	white	100 x 100cm	60gsm	250 sheets
	WHCREPE1014	white	100 x 140cm	60gsm	125 sheets
	WHCREPE120U	white	120 x 120cm	60gsm	125 sheets
	WHCREPE1315U	white	130 x 150cm	60gsm	125 sheets

A




B



C



D





MILLENNIUM BLUE LINER

Lightweight and easy to use liners produced with a matte blue finish to minimise glare. 100% viscose-based product. 40gsm with similar absorbency to 110gsm tray liners. Compatible with steam sterilisation. Single use.

Product code	Size	Qty
MILB002	310 x 410mm	1000 sheets
MILB003	450 x 550mm	1000 sheets
MILB004	240 x 250mm	2400 sheets
MILB005	250 x 500mm	1200 sheets
MILB006	280 x 500mm	1200 sheets
MILB007	250 x 320mm	2000 sheets
MILB008	320 x 450mm	1200 sheets
MILB009	350 x 500mm	1200 sheets
MILB011	320 x 460mm	1000 sheets
MILB012	180 x 280mm	2000 sheets
MILB016	450 x 610mm	1000 sheets
MILB017	300 x 600mm	1200 sheets
MILB019	305 x 305mm	2000 sheets
MILB020	500 x 500mm	1000 sheets
MILB021	130 x 270mm	2000 sheets
MILB022	305 x 610mm	1000 sheets
MILB025	250 x 250mm	2400 sheets
MILB027	250 x 270mm	1000 sheets
MILB029	400 x 500mm	1000 sheets
MILB056	500 x 600mm	1000 sheets
MILB120X250	120 x 250mm	4800 sheets
MILB130X230	130 x 230mm	2000 sheets
MILB170X270	170 x 270mm	1000 sheets
MILB250X300	250 x 300mm	2000 sheets
MILB300X400	300 x 400mm	1000 sheets
MILB305X450	305 x 450mm	1200 sheets
MILB350X450	350 x 450mm	1200 sheets
MILB375X600	375 x 600mm	1000 sheets
MILB500X900	500 x 900mm	400 sheets





Sterilisation Packaging – Liners

A TRAY LINERS

Our Tray Liners have high absorbency and rapid dispersion to ensure condensation is removed during autoclaving, encouraging dry packs and improving the efficacy of the sterilisation process. To be used above or beneath trays to protect against basket. Thickness, 80gsm. Compatible with steam sterilisation. Single use.

Product code	Size	Qty
TLC001	305 x 15m	30 rolls/c
TLC002	457 x 15m	20 rolls/c
TLC004	200 x 305mm	2000 sheets
TLC005	250 x 250mm	2400 sheets
TLC006	250 x 500mm	1200 sheets
TLC007	305 x 410mm	1000 sheets
TLC008	457 x 550mm	800 sheets
TLC009	305 x 305mm	2000 sheets
TLC010	305 x 457mm	1000 sheets
TLC011	305 x 610mm	1000 sheets
TLC012	457 x 457mm	1000 sheets
TLC013	457 x 610mm	600 sheets
TLC014	250 x 490mm	1200 sheets
TLC017	400 x 457mm	1000 sheets
TLC018	305 x 540mm	1000 sheets
TLC019	430 x 457mm	1000 sheets
TLC020	350 x 457mm	1000 sheets
TLC022	360 x 540mm	1000 sheets
TLC023	356 x 457mm	1000 sheets
TLC024	457 x 692mm	600 sheets
TLC027	350 x 500mm	1000 sheets
TLC029	250 x 400mm	1200 sheets
TLC030	350 x 550mm	1000 sheets
TLC031	270 x 350mm	1000 sheets
TLC034	330 x 450mm	1000 sheets
TLC036	305 x 500mm	1000 sheets
TLC045	400 x 500mm	1000 sheets
TLC050	250 x 320mm	1000 sheets
TLC051	250 x 305mm	2000 sheets
TLC052	305 x 520mm	1000 sheets
TLC055	250 x 550mm	1200 sheets
TLC058	280 x 520mm	1000 sheets
TLC060	280 x 405mm	1000 sheets

K

Sterilisation Packaging





Sterilisation Packaging – Liners

TRAY LINERS CONTINUED

Our Tray Liners have high absorbency and rapid dispersion to ensure condensation is removed during autoclaving, encouraging dry packs and improving the efficacy of the sterilisation process. To be used above or beneath trays to protect against basket. Thickness, 80gsm. Compatible with steam sterilisation. Single use.

Product code	Size	Qty
TLC066	250 x 440mm	1000 sheets
TLC069	140 x 230mm	2000 sheets
TLC071	240 x 457mm	1000 sheets
TLC072	178 x 279mm	2400 sheets
TLC078	305 x 660mm	1000 sheets
TLC079	457 x 480mm	1000 sheets
TLC081	305 x 355mm	1000 sheets
TLC085	250 x 480mm	1000 sheets
TLC086	180 x 280mm	1000 sheets
TLC087	380 x 380mm	1000 sheets
TLC095	290 x 560mm	1000 sheets
TLC096	250 x 530mm	1200 sheets
TLC098	275 x 290mm	2000 sheets
TLC100	305 x 350mm	1000 sheets
TLC102	160 x 260mm	2400 sheets
TLC103	290 x 440mm	1000 sheets
TLC104	240 x 280mm	2400 sheets
TLC107	130 x 175mm	2000 sheets
TLC108	350 x 610mm	1000 sheets
TLC109	130 x 230mm	2000 sheets
TLC110	150 x 250mm	2000 sheets
TLC114	250 x 415mm	1000 sheets
TLC115	315 x 525mm	1000 sheets
TLC116	457 x 500mm	1000 sheets
TLC120	260 x 555mm	1200 sheets
TLC122	457 x 630mm	400 sheets
TLC220	220 x 305mm	2000 sheets
TLC3044	305 x 440mm	1000 sheets
TLC340	250 x 340mm	2000 sheets
TLC3535	350 x 350mm	1000 sheets
TLC3548	350 x 480mm	1000 sheets
TLC5060	500 x 600mm	600 sheets



Sterilisation Packaging – Sterilisation Tape

STERILISATION TAPE


A wide range of unique, manufactured tapes that demonstrate excellent binding characteristics and comply with various types of sterilisation wraps.


INDICATOR AUTOCLAVE TAPES

	Product code	Colour	Adhesive	Indicator	Size	Qty
A	01130	beige	normal	steam	19mm x 50m	1 reel
B	01131	beige	normal	EtO	19mm x 50m	1 reel
C	01500	blue	strong	steam	19mm x 50m	1 reel
D	01506	blue	strong	steam	25mm x 50m	1 reel

ADHESIVE TAPES

	Product code	Description	Size	Qty
	ST130PL25	steam plain autoclave tape, (yellow label), not latex free	25mm x 50m	36 pcs
E	01502	crepe tape, brown	50mm x 50m	36 pcs





Sterilisation Packaging – Sterilisation Tape

STERILISATION TAPE

A wide range of unique, manufactured tapes that demonstrate excellent binding characteristics and comply with various types of sterilisation wraps.

LATEX-FREE TAPES

	Product code	Description	Size	Qty
A	LF19TAPE3311	regular	19mm x 50m	48 pcs
B	LF25TAPE3311	regular	25mm x 50m	36 pcs
C	LFWPTAPE19	heavy duty, waterproof paper	19mm x 50m	48 pcs
D	LFWPTAPE25	heavy duty, waterproof paper	25mm x 50m	36 pcs
E	LF25TAPEP321	extra strong adhesive	25mm x 50m	36 pcs

SMS AUTOCLAVE TAPE

	Product code	Description	Size	Qty
	ST101SMS19	steam autoclave tape, (blue label), not latex free	19mm x 50m	48 pcs





Sterilisation Packaging – Transport Wrap

A CLINIPAK TRANSPORT WRAP

The ultimate handling material used to combat damaged packs due to its strength and tear resistance. Our range is convenient when handling heavy orthopaedic sets and uses distinctive colours to easily distinguish them from barrier layers. All Transport Wraps can be supplied on their own, interleaved with any compliant barrier layer or introduced as a triple layer to guarantee a strong, high quality and cost-effective wrap combination. Available in colours purple, orange, yellow and red. Compatible with steam and EtO sterilisation. Bespoke sizes available. Single use.

Product code	Size	GSM	Qty
TWRAP050	50 x 50cm	50gsm	300 sheets
TWRAP060	60 x 60cm	50gsm	350 sheets
TWRAP075	75 x 75cm	50gsm	300 sheets
TWRAP090	90 x 90cm	50gsm	300 sheets
TWRAP100	100 x 100cm	50gsm	200 sheets
TWRAP1012	100 x 120cm	50gsm	120 sheets
TWRAP1014	100 x 140cm	50gsm	120 sheets
TWRAP120	120 x 120cm	50gsm	120 sheets
TWRAP1215	120 x 150cm	50gsm	100 sheets
TWRAP1315	130 x 150cm	50gsm	100 sheets
TWRAP140	140 x 140cm	50gsm	100 sheets
TWRAP150	150 x 150cm	50gsm	100 sheets
TWRAP1715	170 x 150cm	50gsm	50 sheets
TWRAP1815	180 x 150cm	50gsm	50 sheets

B CLINIPAK TRANSPORT WRAP DRY

Protection against damaged packs, with strength and tear resistance with a patented absorbent layer. Particularly suitable for use when handling heavy orthopaedic sets. Available in purple, to distinguish from barrier wraps and also helps to support identification during transport and storage. Compatible with steam and EtO sterilisation. Single use.

Product code	Size	GSM	Qty
TDRY075	75 x 75cm	60gsm	300 sheets
TDRY090	90 x 90cm	60gsm	300 sheets
TDRY100	100 x 100cm	60gsm	200 sheets
TDRY120	120 x 120cm	60gsm	120 sheets



Sterilisation Packaging – Amcor® Pouches and Reels

AMCOR® ULTRA POUCHES AND REELS

Designed to package a large variety of single use or reusable medical instruments creating a micro barrier. A peel open system to eliminate risk of tearing. Compatible with steam sterilisation with steam process indicators. Sterilisation indicators gives the assurance of time, temperature and protection. Available in pouches or reels. Single use.

A AMCOR® ULTRA FLAT POUCHES

Product code	Size	Qty
98ULT00830	80mm x 300mm	2000 pcs
98ULT01220	120mm x 200mm	1200 pcs
98ULT01240	120mm x 400mm	1000 pcs
98ULT01630	160mm x 300mm	1200 pcs
98ULT01650	160mm x 500mm	1000 pcs
98ULT02240	220mm x 400mm	500 pcs
98ULT02260	220mm x 600mm	400 pcs
98ULT02550	250mm x 500mm	500 pcs
98ULT02845	285mm x 450mm	500 pcs
98ULT02860	285mm x 600mm	300 pcs
98ULT03255	320mm x 550mm	300 pcs
98ULT03845	380mm x 450mm	400 pcs
98ULT03865	380mm x 650mm	300 pcs
98ULT03869	380mm x 690mm	300 pcs

B AMCOR® ULTRA FLAT REELS

Product code	Size	Qty
99ULT05070	50mm x 70m	5 pcs
99ULT08070	80mm x 70m	3 pcs
99ULT12070	120mm x 70m	3 pcs
99ULT16070	160mm x 70m	2 pcs
99ULT22070	220mm x 70m	1 pc
99ULT25070	250mm x 70m	1 pc
99ULT28570	285mm x 70m	1 pc
99ULT32070	320mm x 70m	1 pc
99ULT38070	380mm x 70m	1 pc
99ULT42070	420mm x 70m	1 pc
99ULT47070	470mm x 70m	1 pc
99ULT52070	520mm x 70m	1 pc
99ULT58070	580mm x 70m	1 pc

A



B



Sterilisation Packaging – Wipak Pouches and Reels

WIPAK POUCHES AND REELS

Wipak Steriking® transparent sterilisation Pouches and Reels for manual packaging of instruments and similar items. The range consists of heat-sealable pouches and reels, in both flat and gusseted form. In addition, flat pouches are also available in self-seal format. All the Wipak heat seal and self-seal pouches have process indicator strips for both steam and EtO gas. Single use.

WIPAK FLAT POUCHES

Product code	Size	Qty
WI-S1	50 x 250mm	2400/ctn
WI-S2	75 x 200mm	2400/ctn
WI-S4	75 x 270mm	3600/ctn
WI-S5	100 x 270mm	3000/ctn
WI-S6	150 x 270mm	1800/ctn
WI-S7	205 x 270mm	600/ctn
WI-S8	100 x 300mm	1200/ctn
WI-S9	100 x 400mm	1200/ctn
WI-S10	205 x 400mm	600/ctn
WI-S11	250 x 500mm	600/ctn
WI-S12	100 x 570mm	1800/ctn
WI-S13	150 x 400mm	1200/ctn
WI-S14	300 x 570mm	600/ctn
WI-S15	150 x 200mm	1200/ctn
WI-S16	250 x 380mm	1200/ctn
WI-S17	50 x 200mm	3600/ctn
WI-S18	75 x 300mm	4200/ctn
WI-S19	100 x 350mm	1200/ctn
WI-S20	150 x 300mm	1800/ctn
WI-S21	300 x 500mm	600/ctn
WI-S22	75 x 250mm	2400/ctn
WI-S23	75 x 230mm	1800/ctn
WI-S25	100 x 150mm	2400/ctn
WI-S27	120 x 400mm	1200/ctn
WI-S28	150 x 520mm	1200/ctn
WI-S30	270 x 440mm	600/ctn
WI-S31	420 x 500mm	500/ctn
WI-S32	420 x 600mm	500/ctn
WI-S34	160 x 600mm	1200/ctn



Sterilisation Packaging – Wipak Pouches and Reels



A WIPAK GUSSETED POUCHES

Product code	Size	Qty
WI-B35	100 x 50 x 300mm	500/ctn
WI-B31	150 x 50 x 400mm	500/ctn
WI-B36	200 x 55 x 400mm	250/ctn
WI-B37	200 x 65 x 480mm	250/ctn

B WIPAK GUSSETED POUCHES

Product code	Size	Qty
WI-RB50	75mm x 25mm x 100m	4/ctn
WI-RB51	100mm x 50mm x 100m	3/ctn
WI-RB52	150mm x 50mm x 100m	2/ctn
WI-RB53	200mm x 55mm x 100m	1/ctn
WI-RB54	250mm x 65mm x 100m	1/ctn
WI-RB55	300mm x 80mm x 100m	1/ctn
WI-RB56	360mm x 80mm x 100m	1/ctn
WI-RB57	400mm x 80mm x 100m	1/ctn

C WIPAK SELF-SEALED POUCHES

Product code	Size	Qty
WI-SS1	90 x 200mm	1000/ctn
WI-SS2	90 x 250mm	1000/ctn
WI-SS3	90 x 570mm	1000/ctn
WI-SS4	130 x 270mm	1000/ctn
WI-SS5	200 x 350mm	1200/ctn
WI-SS6	250 x 400mm	600/ctn
WI-SS7	300 x 450mm	600/ctn

A



A WIPAK FLAT REELS

Product code	Size	Qty
WI-R39	50mm x 200m	5/ctn
WI-R40	75mm x 200m	4/ctn
WI-R41	100mm x 200m	3/ctn
WI-R42	150mm x 200m	2/ctn
WI-R43	200mm x 200m	1/ctn
WI-R44	250mm x 200m	1/ctn
WI-R45	300mm x 200m	1/ctn
WI-R46	350mm x 200m	1/ctn
WI-R47	400mm x 200m	1/ctn

K

STERILISATION TRAYS AND ACCESSORIES

Our range of high-quality stainless steel washing and instrument trays and baskets are produced in line with industry standard sizes and can be made to fit customer requirements subject to a minimum order quantity.

The sterilisation tray range holds smaller items that are difficult to keep track of during sterilisation as well as larger items with small parts. Our mesh trays and instrument trays are heavy duty stainless steel for tough, clean, and corrosion-resistant surfaces.

Nic-Nac Baskets are designed to transport small items through the wash and autoclave process. All Nic-Nac Baskets are produced with a perforated stainless-steel sheet with a scalloped hinged lid.

CBM basket ranges are manufactured with AIS1304 grade stainless steel and can be used within sterilisation containers or with sterilisation wrap. Due to electro polishing, they are resistant to colour changes and have a surface treatment that increases their water drainage capacity.

Plastic Sterilisation Trays are designed to cover most requirements for protecting instruments during sterilisation and transportation. Grid holes in the tray provide drainage, fixing positions for various inserts, and buttons to raise silicone pin mats from the base. Transparent lids are strongly fitted to provide additional security and ensure all contents are highly visible. Suitable for low temperature sterilisation processes.

We also have our range of accessories such as Corner Protectors, Protective Mats, Soaker Sheets and Labels that are effective solution to reduce damage of trays but also instruments.

Trinket Baskets

Nic-Nac Baskets

Sieve Trays

Flat Base Baskets

Mesh Baskets

Micro Mesh Trays

Stainless Steel Baskets

CBM Baskets

Plastic Sterilisation Trays

Stainless Steel Scope Baskets

Silicone Inserts

Corner Protectors

Soaker Sheets

Protective Mats

Labels

Please note: Unless specified, all products are not made with natural rubber latex.





TRINKET BASKETS

Small, stainless steel trinket baskets to hold fine items such as screws and bolts. Available with lids and designed to transport small items through the washing and autoclave process. Compatible with hydrogen peroxide, EtO and steam sterilisation up to 137°C.

	Product code	Description	Size	Qty
A	CIL0705	trinket box	90 x 60 x 25mm	1 pc
	CIL0706	trinket box, perforated lid	105 x 70 x 25mm	1 pc
	CIL0707	trinket box	120 x 80 x 25mm	1 pc
B	CIL76-0175	trinket box, perforated sides	70 x 30 x 15mm	1 pc
	CIL76-0185	trinket box, wire mesh lid	90 x 60 x 25mm	1 pc
	CIL76-0186	trinket box, wire mesh lid	105 x 70 x 25mm	1 pc
	CIL76-0187	trinket box, wire mesh lid	120 x 80 x 25mm	1 pc
C	CIL76-0189	trinket box, perforated sides and base	150 x 100 x 25mm	1 pc
	CIL0406	fine mesh basket	80 x 80 x 35mm	1 pc
	CIL76-0204	fine mesh basket	105 x 70 x 25mm	1 pc
D	SJ/165-009	fine mesh basket, with lid	80 x 80 x 60mm	1 pc

Sterilisation Trays and Accessories - Nic-Nac Baskets

NIC-NAC BASKETS

A range of small Nic-Nac Baskets manufactured from fine perforated stainless steel with scalloped hinged lids and catch. Designed to accommodate those small items that get lost in a tray and for transportation through the washing process. Compatible with hydrogen peroxide, EtO and steam sterilisation up to 134°C.

	Product code	Size	Qty
A	NIC001R	80 x 30mm	1 pc
B	NIC002	178 x 127 x 25mm	1 pc
C	NIC003	229 x 178 x 25mm	1 pc
C	NIC004DEEP	180 x 180 x 50mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays



SIEVE FOR SMALL PARTS

Small part sieves and accessories to help you to organise your work processes.
Can be used for reprocessing, transportation, and storage for the protection of instruments or small parts. Compatible with steam sterilisation up to 134°C.

	Product code	Description	Size	Qty
A	08405	without lid and clips	236 x 121 x 46mm	1 pc
	08947	lid for ref 08405, with clips	236 x 121 x 8mm	1 pc
B	08406	set with lid and clips	121 x 121 x 29mm	1 pc
C	08408	set with lid and clips	121 x 73 x 29mm	1 pc
D	08409	without lid and clips	121 x 73 x 29mm	1 pc
E	08430	lid for ref 08409, with clips	121 x 73 x 7mm	1 pc
	08448	set with lid and clips	121 x 121 x 29mm	1 pc
F	08449	set with lid and clips	121 x 121 x 47mm	1 pc
G	08450	without lid and clips	121 x 121 x 29mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays

SIEVE FOR SMALL PARTS

Small part sieves and accessories to help you to organise your work processes.
Can use for reprocessing, transportation, and storage for the protection of instruments or small parts.
Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
08452	lid for ref 08450, with clips	121 x 121 x 7mm	1 pc
08451	without lid and clips	121 x 121 x 47mm	1 pc
08454	with lid and clips	236 x 73.5 x 28mm	1 pc
A 08457	without lid and clips	236 x 73,5 x 28mm	1 pc
08461	lid for ref 08457, with 2 clips	236 x 73.5 x 7mm	1 pc
B 08463	set with lid and clips	236 x 121 x 28mm	1 pc
C 08464	set with lid and clips	236 x 121 x 28mm	1 pc
08465	without lid and clips	236 x 121 x 28mm	1 pc
08467	lid ref 08465, with 2 clips	236 x 121 x 7mm	1 pc
08466	without lid and clips	236 x 121 x 28mm	1 pc
08468	lid for ref 08466, with 2 clips	236 x 121 x 7mm	1 pc
D 08442	set with lid and clips	63 x 63 x 23.4mm	1 pc
E 08434	with lid	80 x 80 x 40mm	1 pc
F 08435	with lid	80 x 40 x 30mm	1 pc
G 08896	with lid	40 x 40 x 25mm	1 pc





Sterilisation Trays and Accessories - Sieve Trays

SIEVE FOR SMALL PARTS

Small part sieves and accessories to help you to organise your work processes. Can be used for reprocessing, transportation, and storage for the protection of instruments or small parts. Compatible with steam sterilisation up to 134°C.

	Product code	Description	Size	Qty
A	08898	with lid	160 x 40 x 20mm	1 pc
B	08899	with lid	200 x 40 x 40mm	1 pc
C	08900	with lid	80 x 40 x 40mm	1 pc
D	08901	with lid	80 x 80 x 40mm	1 pc
E	08906	with lid	250 x 200 x 25mm	1 pc
	08907	with lid	300 x 180 x 25mm	1 pc
F	08772	with lid	210 x 80 x 40mm	1 pc
	08948	with lid	200 x 60 x 40mm	1 pc
	08875	with lid	300 x 150 x 70mm	1 pc
G	08950	with lid	300 x 150 x 70mm	1 pc
H	08762	without lid	200 x 50 x 50mm	1 pc





Sterilisation Trays and Accessories - Sieve Trays

SIEVE FOR SMALL PARTS WITH LENGTH DIVIDERS


Small part sieves and accessories to help you to organise your work processes. Can be used for reprocessing, transportation, and storage for the protection of instruments or small parts. Lengthwise dividers for small parts sieves. For individual division. Includes two clips.

	Product code	Description	Size	Qty
A	08410	lengthwise divider 1/1	224 x 24mm	1 pc
	08412	lengthwise divider 1/1	224 x 42mm	1 pc
B	08482	lengthwise divider 2/3	300 x 43mm	1 pc
C	08436	lengthwise divider 1/1	465 x 43mm	1 pc
D	08414	crosswise divider 1/1	109 x 24mm	1 pc
E	08437	crosswise divider 1/3	71 x 43mm	1 pc
F	08481	crosswise divider 1/1	61.5 x 24mm	1 pc
	08483	crosswise divider 1/1	232 x 43mm	1 pc
	08484	crosswise divider 2/3	152 x 43mm	1 pc
	08485	crosswise divider 2/3	152 x 43mm	1 pc

COMPARTMENTAL SIEVE FOR SMALL PARTS

A Compartmental Sieve Tray for small parts with 65 inner compartments. Made from stainless steel.

	Product code	Description	Size	Qty
G	08880	65 inner compartments, with lid, hinge, spring cap and side walls	40 x 15mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays



SIEVE TRAYS

Our range of stainless steel Sieve Trays, with and without lids, helps you guarantee optimal and efficiency in your work process. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
A 08808	without lid	480 x 250 x 100mm	1 pc
08809	without lid	560 x 250 x 50mm	1 pc
08810	without lid	560 x 250 x 70mm	1 pc
B 08415	without lid	240 x 250 x 50mm	1 pc
08416	without lid	480 x 250 x 50mm	1 pc
08417	without lid	480 x 250 x 70mm	1 pc
08580	with lid	540 x 250 x 50mm	1 pc
C 08761	with lid	400 x 280 x 70mm	1 pc
08581	with lid	540 x 250 x 100mm	1 pc
D 08801	rochard sieve, with lid	725 x 275 x 113mm	1 pc




Sterilisation Trays and Accessories - Sieve Trays

SIEVE TRAYS AND MATCHING LIDS

Sieve Trays with matching lids, made of Chrome Nickel Steel including accessories as well. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
A 08455	set 1/6 XL, with lid, set comprising: 1 x system small parts sieve 1/6 XL with lid and clips	396 x 73.5 x 47mm	1 pc
B 08456	1/6, without lid	396 x 73.5 x 47mm	1 pc
C 08459	1/6 XL, without lid	396 x 73.5 x 47mm	1 pc
D 08453	set 1/6, with lid, set comprising: 1 x system small parts sieve 1/6 with lid and clips	396 x 73.5 x 47mm	1 pc
E 08460	lid for ref 08456 with 2 clips	236 x 73.5 x 7mm	1 pc
08462	lid for ref 08459 with 2 clips	396 x 73.5 x 7mm	1 pc
08470	set 1/3L, with lid, set comprising: 1 x system small parts sieve 1/3L with lid and clips	466 x 73.5 x 47mm	1 pc



Sterilisation Trays and Accessories - Sieve Trays

SIEVE TRAY SETS

Sieve Trays Sets, made of Chrome Nickel Steel including accessories as well. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
A 08469	1 x system small parts sieve with lid and clips	466 x 73.5 x 28mm	1 pc
B 08471	1 x system small parts sieve with lid and clips	236 x 236 x 28mm	1 pc
C 08472	1 x system small parts sieve with lid and clips	236 x 236 x 46mm	1 pc

SIEVE TRAYS WITH DROP HANDLES

Sieve Trays with drop handles are made of stainless steel and without lids. Compatible with steam sterilisation up to 134°C.

Product code	Size	Qty
D 08908	400 x 250 x 70mm	1 pc
E 08909	250 x 240 x 50mm	1 pc



Sterilisation Trays and Accessories - Flat Base Baskets



FLAT BASE BASKETS

Our Flat Base Baskets range are stainless steel for instrument washing, disinfection, sterilisation, storage, and handling. Compatible with steam sterilisation processes up to 137°C.

Product code	Description	Size	Qty
CIL76-0206	flat base basket, perforated sides	150 x 50 x 35mm	1 pc
A CIL76-0049	flat base basket, perforated sides	480 x 250 x 70mm	1 pc
CIL76-0039	flat base basket, perforated sides	380 x 240 x 50mm	1 pc
CIL76-0047	flat base basket, perforated sides	480 x 250 x 50mm	1 pc
CIL76-0048	flat base basket, perforated sides	480 x 250 x 60mm	1 pc
B CIL76-0030	flat base basket, perforated sides	250 x 240 x 50mm	1 pc
CIL76-0057	flat base basket, perforated sides	480 x 360 x 60mm	1 pc
CIL76-0058	flat base basket, perforated sides	480 x 480 x 50mm	1 pc
CIL76-0060	flat base basket, perforated sides	480 x 480 x 70mm	1 pc
CIL76-0066	flat base basket, perforated sides	520 x 250 x 50mm	1 pc
CIL76-0076	flat base basket, perforated sides	550 x 250 x 60mm	1 pc
CIL76-0079	flat base basket, perforated sides	550 x 430 x 50mm	1 pc
CIL76-0086	flat base basket, perforated sides	600 x 300 x 60mm	1 pc
CIL76-0087	flat base basket, perforated sides	600 x 350 x 60mm	1 pc
CIL76-0088	flat base basket, perforated sides	670 x 300 x 70mm	1 pc
CIL76-0089	flat base basket, perforated sides	700 x 250 x 60mm	1 pc
CIL76-0090	flat base basket, perforated sides	700 x 300 x 90mm	1 pc
CIL76-0095	lid for flat base basket	250 x 240mm	1 pc
CIL76-0096	lid for flat base basket	480 x 250mm	1 pc
CIL76-0097	lid for flat base basket	480 x 480mm	1 pc
CIL0202LID	lid for wire mesh basket	240 x 250mm	1 pc
CIL0210LID	lid for wire mesh basket	480 x 250mm	1 pc
C CIL0202	wire mesh basket	250 x 240 x 60mm	1 pc
D CIL0210	wire mesh basket	480 x 250 x 60mm	1 pc

FLAT BASE DIN BASKETS

Product code	Description	Size	Qty
SJ-131-020	din basket, with fitted lid on top	650 x 250 x 70mm	1 pc
SJ-275-275	stainless steel din basket perforated sides, flat base	275 x 275 x 50mm	1 pc





Sterilisation Trays and Accessories - Mesh Baskets

MESH BASKETS AND LIDS


Stainless steel Fine Mesh Baskets with lids for instrument washing, disinfection, sterilisation, storage, and handling. Compatible with steam sterilisation up to 137°C, and plasma sterilisation processes.

Product code	Description	Size	Qty
CIL76-0143	all wire mesh basket	480 x 250 x 50mm	1 pc
CIL76-0200	all wire mesh basket	80 x 40 x 30mm	1 pc
A CIL76-0202	fine mesh basket	80 x 80 x 55mm	1 pc
CIL76-0210	all wire mesh basket	190 x 100 x 85mm	1 pc
CIL76-0215	fine mesh basket	200 x 100 x 35mm	1 pc
CIL76-0222	fine mesh basket	240 x 90 x 40mm	1 pc
CIL76-0250	fine mesh basket	120 x 155 x 40mm	1 pc
CIL76-0253	fine mesh lid	120 x 155mm	1 pc
CIL76-0256	fine mesh basket, take off lid, hinges and locking clip	120 x 155 x 50mm	1 pc
CIL76-0252	fine mesh basket, with lid	235 x 155 x 40mm	1 pc
CIL76-0475	all wire mesh basket	250 x 170 x 70mm	1 pc

PERFORATED PLATE MESH TRAYS

Stainless steel Perforated Plate Mesh Trays. Suited for soft packaging. Compatible with steam sterilisation up to 134°C.

Product code	Description	Size	Qty
B 08432	with drop handles, without lid	480 x 250 x 105mm	1 pc
08473	with drop handles, without lid	480 x 250 x 60mm	1 pc
08474	with drop handles, without lid	480 x 250 x 65mm	1 pc
08475	with drop handles, without lid	480 x 250 x 100mm	1 pc
C 08476	with drop handles, without lid	240 x 250 x 60mm	1 pc
08478	with drop handles, without lid	240 x 250 x 65mm	1 pc
08579	with turn-lock and lid	478 x 247mm	1 pc
08583	without lid	480 x 250 x 100mm	1 pc
D 08934	flushing mesh basket, with lid	240 x 250 x 60mm	1 pc
E 08942	flushing mesh basket, with lid	480 x 250 x 60mm	1 pc
08943	flushing mesh basket, with lid	240 x 250 x 60mm	1 pc
08945	flushing mesh basket, with lid	220 x 140 x 54mm	1 pc
08975	without lid	120 x 250 x 60mm	1 pc
08976	without lid	240 x 250 x 66mm	1 pc



Sterilisation Trays and Accessories - Micro Mesh Trays

MICRO MESH TRAYS

Micro Mesh Trays are designed to hold smaller items during washing, sterilisation, storage, and transportation. Micro mesh allows you to keep larger items that have small loose parts together. Micro mesh size is approximately 2mm x 2mm. Compatible with steam and gas plasma sterilisation.

	Product code	Description	Size	Qty
A	MT-5500	with lids	40 x 40 x 20mm	1 pc
B	MT-5600	with lids	80 x 40 x 20mm	1 pc
C	MT-5650	with lids	160 x 40 x 20mm	1 pc
D	MT-5700	with lids	70 x 70 x 30mm	1 pc
E	MT-5800	with lids	80 x 80 x 40mm	1 pc
F	MT-5900	with lids	90 x 90 x 50mm	1 pc
G	MT-6020	with lids	180 x 130 x 20mm	1 pc
H	MT-7100	with lids	225 x 110 x 50mm	1 pc
I	MT-7200	with lids	270 x 135 x 40mm	1 pc
J	MT-7300	with lids	400 x 250 x 50mm	1 pc
K	MT-9001	with drop handles	480 x 250 x 50mm	1 pc
L	MT-9101	with drop handles	250 x 240 x 50mm	1 pc
M	MT-9521	with drop handles	700 x 210 x 50mm	1 pc



Sterilisation Trays and Accessories - Stainless Steel Baskets

PREMIUM STAINLESS STEEL BASKETS

Stainless Steel Baskets for holding instruments during sterilisation. Compatible with steam sterilisation up to 137°C.

Product code	Description	Size	Qty
A BSK316	stainless steel basket	480 x 250 x 50mm	1 pc
B BSK317	stainless steel basket	250 x 250 x 50mm	1 pc
C BSK318	stainless steel basket	480 x 480 x 50mm	1 pc
LID316	stainless steel perforated lid to fit BSK316	480 x 250mm	1 pc
LID317	stainless steel perforated lid to fit BSK317	250 x 250mm	1 pc
LID318	stainless steel perforated lid to fit BSK318	480 x 480mm	1 pc
BSK316H	stainless steel basket, with handle	480 x 250 x 50mm	1 pc
BSK317H	stainless steel basket, with handle	250 x 250 x 50mm	1 pc
BSK318H	stainless steel basket, with handle	480 x 480 x 50mm	1 pc
BELBSK116	stainless steel basket, mesh base with handles	480 x 250 x 60mm	1 pc
BELBSK117	stainless steel basket, with handle	250 x 240 x 60mm	1 pc



Sterilisation Trays and Accessories - Stainless Steel Baskets

STRINDBERG TRAYS


A series of stainless steel trays designed to protect delicate instruments during washing and sterilisation. Instruments are held in place with a firm clip protected with a silicone shield.

Product code	Description	Size	Qty
A RON005	full tray, 20 instrument slots	275 x 180 x 25mm	1 pc
RON006	half tray, 9 instrument slots	131 x 180 x 50mm	1 pc
RON007	quarter tray, 5 instrument slots	73 x 180 x 25mm	1 pc


TRAYS AND CASSETTES

Stainless steel tray for the storage of instruments.

Product code	Description	Size	Qty
B 08940	storage tray for instruments, with silicone holder 5-fold with fixing bracket	180 x 85 x 23mm	1 pc



Sterilisation Trays and Accessories - CBM Baskets



MINI AND MICRO BASKETS

Small, stainless steel baskets to hold find items such as screws and bolts. Available with lids and designed to transport small items through the washing and autoclave process.

	Product code	Description	Size	Qty
A	CBM63J	mini basket with lid, without handles	200 x 80 x 35mm	1 pc
B	CBM67F2	micro basket with lid, without handles	170 x 135 x 35mm	1 pc
C	CBM63M	mini basket with lid, without handles	220 x 60 x 35mm	1 pc
D	CBM76D	mini basket with lid and handles	267 x 125 x 49mm	1 pc

STERILISATION BASKETS SYSTEM


CBM Sterilisation Baskets Systems are manufactured with stainless steel and can be used within sterilisation containers or with sterilisation wrap. Due to electro polishing, they are resistant to colour changes and have a surface treatment that increases their water drainage capacity.

	Product code	Description	Size	Qty
E	CBM71/1	stainless steel basket, without lid	240 x 240 x 50mm	1 pc
	CBM71LID	stainless steel basket, with lid	240 x 240mm	1 pc
	CBM74/1	stainless steel basket, without lid	480 x 240 x 50mm	1 pc
	CBM74LID	stainless steel basket, with lid	480 x 240mm	1 pc
	CBM75/1	stainless steel basket, without lid	540 x 240 x 50mm	1 pc
	CBM75LID	stainless steel basket, with lid	540 x 240mm	1 pc
	CBM79/2	stainless steel basket, without lid	650 x 240 x 70mm	1 pc
	CBM79/LID	stainless steel basket, with lid	650 x 240mm	1 pc

SCOPE BASKETS SYSTEM

A range of baskets for the reprocessing, storage, and transportation of arthroscopes and endoscopes, available with integrated silicone instrument holders.

	Product code	Description	Size	Qty
F	CBM48/S	scope basket for 1 scope, with lid	480 x 150 x 80mm	1 pc
G	CBM66/S	scope basket for 1 scope, with lid	450 x 70 x 60mm	1 pc
H	CBM64ECV	empty scope basket with lid and handles	650 x 150 x 70mm	1 pc
I	CBM66AV	empty scope basket with lid and handles	580 x 150 x 80mm	1 pc
	CBM66/FEN	flexible scope basket with lid and 4 silicone supports	580 x 150 x 80mm	1 pc



Sterilisation Trays and Accessories - Plastic Sterilisation Trays

STERILISATION TRAYS

Sterilisation trays made of plastic and stainless steel. The plastic trays and cassettes are used for sterilisation and instrument storage.

Product code	Description	Colour	Size	Qty
A 08422	lid black-transparent with 4mm holes	beige	150 x 65 x 26mm	1 pc
B 08424	lid black-transparent with 6mm holes	beige	230 x 175 x 42mm	1 pc
C 08425	lid black-transparent with 6mm holes	beige	300 x 200 x 50mm	1 pc

WASH TRAYS

Wash Trays made for Da Vinci® systems.

Product code	Description	Size	Qty
VAS4XI	washing tray for Endo Wrist®	660 x 212 x 50mm	1 pc
08798	sieve tray	640 x 340 x 100mm	1 pc
08132	mesh basket with lid, made of stainless steel	670 x 100 x 70mm	1 pc




Sterilisation Trays and Accessories - Stainless Steel Scope Baskets

STAINLESS STEEL SCOPE BASKETS

Stainless steel wire mesh baskets for one or two endoscopes or arthroscopes with integral silicone instrument holders. This lockable range protects scope instruments whilst washing and sterilising.

Product code	Description	Size	Qty
A SB-1210	with silicone inserts	180 x 80 x 52mm	1 pc
B SB-4210	with silicone inserts	460 x 80 x 52mm	1 pc
C SB-4610	with silicone inserts	460 x 160 x 52mm	1 pc
D SB-6520	with silicone inserts, includes small lens basket with fixed bracket	640 x 150 x 77mm (80 x 80 x 40mm)	1 pc
CIL2728	scope basket	290 x 80 x 52mm	1 pc
CIL2727	scope basket	460 x 80 x 52mm	1 pc
CIL2726	scope basket	670 x 80 x 52mm	1 pc
CIL2701	scope basket, with small mesh basket	640 x 150 x 77mm	1 pc
CIL2729	scope basket for 2 endoscopes	460 x 80 x 52mm	1 pc
CIL2730	double scope basket	290 x 80 x 52mm	1 pc
CIL2751	double scope basket	290 x 160 x 52mm	1 pc
CIL2753	double scope basket	460 x 160 x 52mm	1 pc
CIL76-0001	endoscope / arthroscope basket	640 x 150 x 77mm	1 pc
CIL76-0010	universal scope basket	580 x 160 x 52mm	1 pc
CIL76-0011	universal scope basket	290 x 160 x 52mm	1 pc
CIL76-0012	universal scope basket	460 x 160 x 52mm	1 pc
SJ-102-001P	stainless steel, with 1 pin	460 x 80 x 52mm	1 pc
SJ-102-002P	stainless steel, with 2 pin	460 x 80 x 52mm	1 pc
SJ-103-002	with small mesh basket	640 x 150 x 77mm	1 pc





Sterilisation Trays and Accessories - Silicone Inserts


SILICONE INSERTS

For custom configuration within a tray to accommodate various sizes of instruments.

	Product code	Description	Opening diameter	Size	Qty
A	A237593	straight opening, 3 slots	N/A	57 x 15mm	1 pc
B	A237596	straight opening, 6 slots	N/A	57 x 15mm	1 pc
	A387596	straight opening, 6 slots	N/A	95 x 19mm	1 pc
	A3875912	straight opening, 12 slots	N/A	95 x 19mm	1 pc
C	A237591K	round opening, 1 slot	6mm	57 x 15mm	1 pc
D	A237592K	round opening, 2 slots	6mm	57 x 15mm	1 pc
E	A387592	round opening, 2 slots	6mm	95 x 19mm	1 pc
	A2616902	round opening, 2 slots	10mm	63 x 25mm	1 pc
	A2616901	round opening, 2 slots	10 & 13mm	53 x 25mm	1 pc
	A387593	round opening, 3 slots	6mm	95 x 19mm	1 pc
F	A237593K	round opening, 3 slots	6mm	57 x 15mm	1 pc
	A387594	round opening, 4 slots	6mm	95 x 19mm	1 pc
	A387595	round opening, 5 slots	6mm	95 x 19mm	1 pc
G	A32759P2	long, round opening, 1 slot	7mm	89 x 32mm	1 pc
H	A32759P1	long, round opening, 1 slot	16mm	89 x 32mm	1 pc
I	A327593	long, round opening, 3 slots	7mm	89 x 32mm	1 pc
J	A327593L	long, round opening, 3 slots	16 & 7 & 13mm	89 x 32mm	1 pc



A



B




C



D



E



Sterilisation Trays and Accessories - Corner Protectors

SURGICAL TRAY CORNER PROTECTORS

Tray Corner Protectors are made from medical-grade silicone and are intended for use in the organization and storage of surgical trays. Compatible with steam sterilisation. Single use.

Product code	Description	Colour	Size	Qty
A CS-9125C	foam, 4 pre-cut corners	blue	152 x 152 x 32mm	80 sheets

CORNER PROTECTORS

Tray Corner Protectors are made from medical-grade silicone and are intended for use in the organization and storage of surgical trays. Heat-resistant up to 134°C.

Product code	Material	Colour	Size	Qty
B 42006	silicone	pink	73.1 x 73.1 x 18.5mm	80 pcs
42007	card	white	100 x 100 x 50mm	600 pcs
C 42008	foam	blue	75 x 75 x 32mm	288 pcs
42010	card	white	100 x 100 x 50mm	600 pcs
42011	silicone	yellow	76 x 54 x 33mm	4 pcs
D 50062	silicone	violet	76 x 54 x 33mm	4 pcs

CARD AND CREPE CORNER PROTECTORS

Tray Corner Protectors are intended for use in the organization and storage of surgical trays. Suitable for steam sterilisation, they are available in crepe or card.

Product code	Material	GSM	Qty
CC001	card	171gsm	200 pcs
E CC002	crepe	171gsm	100 pcs

A



K

Sterilisation Trays and Accessories - Soaker Sheets

SOAKER SHEETS

140gsm Soaker Sheets are for use under trays to protect inner wraps from sharp corners and feet on trays and to absorb condensate during sterilisation. Soaker Sheets can also be used as a shelf liner for added protection. Single use.



Product code	Size	Qty
A SOK001SS	750 x 500mm	50 sheets
SOK002SS	500 x 375mm	100 sheets
SOK003SS	610 x 305mm	100 sheets
SOK004SS	500 x 200mm	100 sheets
SOK005SS	500 x 500mm	100 sheets
SOK006SS	250 x 250mm	100 sheets
SOK007SS	500 x 280mm	100 sheets
SOK008SS	1219 x 560mm	50 sheets
SOK009SS	1524 x 560mm	50 sheets
SOK010SS	914 x 508mm	50 sheets
SOK011SS	914 x 588mm	50 sheets
SOK012SS	1000 x 600mm	50 sheets
SOK013SS	1295 x 350mm	50 sheets
SOK016SS	200 x 125mm	1000 sheets



Sterilisation Trays and Accessories - Protective Mats

ABSORBENT PROTECTIVE MATS

Absorbent Mats are intended for placement on the bottom of a sterilisation tray to absorb wetness and humidity during processing. The 'dry' version, with an added absorbent layer, is also available to provide rapid moisture dispersion during steam sterilisation. Compatible with steam sterilisation. Single use.

PURPLE

Product code	Description	Size	Qty
A MAT001	double din	750 x 500mm	50 sheets
MAT002	double din	500 x 375mm	100 sheets
MAT003	double din	610 x 305mm	100 sheets
MAT004	din	500 x 250mm	100 sheets
MAT005	half din	250 x 250mm	100 sheets
MAT007	double din	500 x 500mm	100 sheets

PURPLE DRY

Product code	Description	Size	Qty
A DRYMAT001	double din	750 x 500mm	50 sheets
DRYMAT002	double din	500 x 375mm	100 sheets
DRYMAT003	double din	610 x 305mm	100 sheets
DRYMAT004	din	500 x 250mm	100 sheets
DRYMAT005	half din	250 x 250mm	100 sheets
DRYMAT006	din	500 x 280mm	100 sheets
DRYMAT007	double din	500 x 500mm	100 sheets
DRYMAT008	double din	600 x 600mm	100 sheets
DRYMAT009	double din	750 x 750mm	100 sheets
DRYMAT010	double din	900 x 900mm	100 sheets




GREEN DRY

Product code	Description	Size	Qty
B DRYMATGRN001	double din	750 x 500mm	50 sheets
DRYMATGRN002	double din	500 x 375mm	100 sheets
DRYMATGRN003	double din	610 x 305mm	100 sheets

A

B

K



Sterilisation Trays and Accessories - Protective Mats and Labels

GRID INSERTS AND MATS

Grid inserts and mats made of silicone. Heat resistant. Suitable for washer-disinfectors. Flexible.

Product code	Description	Colour	Heat-resistant	Size	Qty
A 08601	honeycomb structure, for blue perforated mesh trays	blue	up to +230°C	480 x 250mm	1 pc
B 08609	honeycomb structure, for blue perforated mesh trays	blue	up to +230°C	480 x 250mm	1 pc
C 08605	lattice structure	blue	-40°C up to +300°C	600 x 400 x 2mm	1 pc

SILICONE NUB MATS

Silicone nub mats. Suitable for washer-disinfectors. Flexible.

Product code	Description	Colour	Size	Qty
08426	compatible with ref 08422	blue	150 x 55mm	1 pc
08428	compatible with ref 08424	blue	200 x 135mm	1 pc
08429	compatible with ref 08425	blue	250 x 155mm	1 pc

SILICONE MESH COVER

Silicone mesh cover for system tray. Heat resistant up to 134°C. Flexible.

Product code	Description	Size	Qty
D 08881	for sieve trays 1/1	430 x 230mm/ 415 x 210mm	1 pc

STAINLESS STEEL SIEVE TRAY LABELS

Sieve tray labels are vital since hospitals are constantly in search of ever more sophisticated means of tracking sieve trays. Sieve tray labels are an ideal solution to this problem in terms of both colour coding and providing information. Heat resistant up to 134°C. For permanent use.

Product code	Description	Qty
E 01590	clip label, suitable for double frame	25 pcs
F 01445	label	50 pcs
G 01588	tie on label with hole (6 x 10mm)	5 pcs
H 01452	label with hole (4 x 9mm)	50 pcs





LABEL HOLDER

Stainless Steel clip label holder for permanent use and heat resistant.

	Product code	Description	Size	Qty
A	01497	label holder, compatible with ref 01550	50 x 20mm	25 pcs
B	01498	label holder, compatible with ref 01948	60 x 40mm	25 pcs

INSERTION LABEL


Compatible with clip label holder ref 01497 and dutiable label holders on containers. Heat resistant up to 134°C. For permanent use. Suitable for laser marking.

	Product code	Colour	Size	Qty
C	01550	white	56 x 18mm	10 pcs
D	01551	red	56 x 18mm	10 pcs

ALUMINIUM LABEL

Suitable for slim label slots on containers and in ref 01497. Suitable for laser marking.

	Product code	Colour	Size	Qty
E	28500	silver	18 x 50mm	25 pcs



CONTAINERS AND ACCESSORIES

CBM Containers are reusable and are installed with aluminium filters to protect and maintain the sterility of surgical instruments throughout the sterilisation process conforming to EN285 and are fully compliant with EN969-8 and ISO 11607-1.

Container Filters for steam sterilisation are available as cotton or synthetic filters with steam indicators in a rectangular or round shape.

Easy to handle paper seals are an easy way to seal containers quickly, simply, and safely.

We have a variety of container labels that come in different designs that offers a simple solution to supplement your documentation.

[CBM Containers](#)

[Container Filters](#)

[Paper Seals](#)

[Container Labels](#)

Please note: Unless specified, all products are not made with natural rubber latex.



A



B



C



D



Containers and Accessories – CBM Containers and Container Filters

CBM CONTAINERS

The CBM Container is installed with permaside aluminium filters, used in steam sterilisers conforming to EN285, and fully compliant with EN969-8 and ISO 11607-1. Anodised aluminium increases corrosion resistance, particularly against modern alkaline washing cycles. The container lid is manufactured from stainless steel for improved scratch resistance. CBM container also employ permanent PTFE filter cartridges which are certified for more than 2000 cycles.

	Product code	Size	Qty
A	CBM2072PS	600 x 300 x 160mm	1 pcs

CONTAINER FILTERS

High quality white paper filters made of non-woven material with reduced particle residues. Steam indicators in a range of shapes and sizes. Single use.

	Product code	Description	Size	Qty
B	01140	round paper filter, with steam indicator	190mm (hole 6mm)	500 pcs
	01152	round paper filter, with steam indicator	190mm (hole 6mm)	500 pcs
C	01148	small paper set filter, round with steam indicator	118 x 235mm	200 pcs
D	01315	paper filter, with steam indicator	240 x 240mm	250 pcs

K

Containers and Accessories


A

STER. DATE _____

EXP. DATE _____

STER. NO. _____

NAME _____

REF 01008  STEAM BROWN


B

STER. DATE _____

EXP. DATE _____

STER. NO. _____

NAME _____


REF01006  STEAM BROWN

C

Datum: _____

Exp. dat.: _____

Name: _____

REF 01004  Dampf = Braun

D

STER. DAT. _____

STER. DATE _____


VERFALL DAT. _____

EXP. DATE _____

STER. NR. _____

STER. NO. _____

NAME _____

REF 01135  DAMPF BRAUN STEAM BROWN

E

STER. DAT. _____

STER. DATE _____


VERFALL-DAT. _____

EXP. DATE _____

STER. NR. _____

STER. NO. _____

NAME _____

DAMPF BRAUN  STEAM BROWN REF01136


F

Date: _____

Exp. Date: _____

Autoclave #: _____ Lot: _____

Name: _____

REF 01169  STEAM = BROWN


G

STER. DATE _____

EXP. DATE _____

STER. NO. _____

NAME _____

REF 01400  STEAM BROWN

Containers and Accessories – Paper Seals and Container Labels

PAPER SEALS

An environmentally alternative to plastic seals. Our white paper seals with protocol label can be attached to all standard containers. With steam indicator, pink to brown. Single use.

Product code	Description	Size	Qty
A 01008	seal with protocol label, with steam indicator, self-adhesive heat resistant up to 134°C	74 x 35mm	1000 pcs
B 01006	seal with protocol, with steam indicator	74 x 35mm	1000 pcs
C 01004	seal with protocol, with steam indicator	117 x 38mm	500 pcs

CONTAINER LABELS

To support the identification of containers throughout the reprocessing cycle, you can insert container labels into the slot provided on the container. With steam indicator and documentation fields. Single use.

Product code	Description	Size	Qty
D 01135	container label, non-adhesive	39.5 x 68mm	1000 pcs
E 01136	container label	39.5 x 67mm	1000 pcs
F 01169	container label	40 x 59 x 15mm	1000 pcs
G 01400	container label, non-adhesive	35 x 80mm	80 pcs

K

Containers and Accessories

INSTRUMENT IDENTIFICATION

The KeyDot is a small, laser-etched, 2D Data Matrix barcode label that can be applied directly to an instrument's surface. Tested for sterilisation efficacy and biocompatibility for use on stainless steel surgical instruments, it has been conveniently incorporated into most of the instrument tracking software applications available today. With the KeyDot there is no specialised training or intricate equipment needed so application can be done on-site, eliminating the hassle of sending your instrument sets out for this tracking solution.

Our Identification Tape is designed to help you organise and identify instruments as they move between sterile processing and various departments. Using tape can help save you time, effort, and money by offering effective organisation and identification of your surgical instrumentation. It is easy to apply and is compatible with steam sterilisation. We have the broadest selection of tape available and it is colourful, easy to use, and economical.

Adhesive Remover is designed to dissolve tape residue on instruments, trays, basins, and other surfaces. Adhesive Remover is the plant based solution you need to remove residue without causing any damage.

KeyDot

Roll Tape

Adhesive Remover

Instrument Marking Tape

Colour Coding

Please note: Unless specified, all products are not made with natural rubber latex.



1 clean application area




2 lift KeyDot from sheet



3 place KeyDot on instrument




4 apply firm pressure and sterilise to secure adhesive



KEYDOT

With the KeyDot there is no specialised training or intricate equipment required. Applying a KeyDot to a surgical instrument is quick and easy: simply clean, lift, place, and press. The instrument then needs to be sterilised to cure the adhesive and assure firm adhesion.

	Product code	Diameter	Qty
A	ID-405	3.2mm	1 pc
B	ID-400	4.8mm	1 pc
C	ID-410	6.4mm	1 pc
D	ID-420	9.5mm	1 pc





Instrument Identification – Roll Tape and Adhesive Remover

ROLL TAPE

Compatible with steam sterilisation. Total length 6350mm. Available in three widths.

Colour	3mm	6mm	9mm
A white	T-01-01	T-02-01	T-03-01
B blue	T-01-02	T-02-02	T-03-02
C green	T-01-03	T-02-03	T-03-03
D red	T-01-04	T-02-04	T-03-04
E yellow	T-01-05	T-02-05	T-03-05
F violet	T-01-06	T-02-06	T-03-06
G orange	T-01-07	T-02-07	T-03-07
H brown	T-01-08	T-02-08	T-03-08
I black	T-01-09	T-02-09	T-03-09

ADHESIVE REMOVE

Our citrus based adhesive remover is designed to remove adhesive residue from instruments, trays, basins and other surfaces without causing any damage.

Product code	Size	Qty
J BC-200-1	250ml	1 bottle
BC-200-6	250ml	6 bottles




K

Instrument Identification


SOLID COLOUR INSTRUMENT MARKING TAPE TAPE 3MM WIDTH

Compatible with steam sterilisation. Customisable.

Each 216 x 279mm sheet includes:
462 pieces of 3 x 25.4mm
88 pieces of 3 x 38.1mm
22 pieces of 3 x 50.8mm
22 pieces of 3 x 127mm
(total of 19m of tape)




	Product code	Colour
A	T-110-02	blue
B	T-110-03	green
C	T-110-04	red
D	T-110-05	yellow
E	T-110-06	purple
F	T-110-12	grey



SOLID COLOUR INSTRUMENT MARKING TAPE TAPE 6.35MM WIDTH

Compatible with steam sterilisation. Customisable.

Each 216 x 279mm sheet includes:
231 pieces of 6.35 x 25.4mm
44 pieces of 6.35 x 38.1mm
11 pieces of 6.35 x 50.8mm
11 pieces of 6.35 x 127mm
(total of 9.50m of tape)



	Product code	Colour
A	T-120-02	blue
B	T-120-03	green
C	T-120-04	red
D	T-120-05	yellow
E	T-120-07	orange
F	T-120-11	pink



Instrument Identification – Instrument Marking Tape




PATTERNED INSTRUMENT MARKING TAPE – 6.35MM WIDTH

Compatible with steam sterilisation and customisable:
One 216 x 279mm sheet pre-cut into 25, 38, 50.8, and 127mm lengths

DID YOU KNOW?

Key Surgical has a specialised Clinical Team which can bring added value to your product experience. We offer free training on how to use CSSD products as well as several training modules. Some of these include: Central Decontamination Unit Reprocessing, Managing Risks in Decontamination, Aseptic Presentation training and more. Additional services from our team also include 180 and 360 Central Decontamination Unit Audits which can cover just areas such as sterile store or the entire decontamination room from clean room to sterile store and transportation.



Instrument Identification – Colour Coding

COLOURCURLS®

Made of polypropylene for identifying instruments. Single use.

	Product code	Description	Colour	Diameter	Qty
A	09620	ColourCurls® mini	light blue	3.5mm	500 pcs
B	09621	ColourCurls® mini	dark blue	3.5mm	500 pcs
C	09622	ColourCurls® mini	light green	3.5mm	500 pcs
D	09623	ColourCurls® mini	dark green	3.5mm	500 pcs
E	09625	ColourCurls® mini	pink	3.5mm	500 pcs
F	09626	ColourCurls® mini	dark red	3.5mm	500 pcs
G	09627	ColourCurls® mini	light yellow	3.5mm	500 pcs
H	09628	ColourCurls® mini	orange	3.5mm	500 pcs
I	09629	ColourCurls® mini	brown	3.5mm	500 pcs
J	09630	ColourCurls® mini	black	3.5mm	500 pcs
K	09631	ColourCurls® mini	light violet	3.5mm	500 pcs
L	09632	ColourCurls® mini	turquoise	3.5mm	500 pcs
M	09639	ColourCurls® maxi	light blue	6mm	500 pcs
N	09640	ColourCurls® maxi	dark blue	6mm	500 pcs
O	09641	ColourCurls® maxi	light green	6mm	500 pcs
P	09642	ColourCurls® maxi	dark green	6mm	500 pcs
Q	09643	ColourCurls® maxi	white	6mm	500 pcs
R	09644	ColourCurls® maxi	pink	6mm	500 pcs
S	09645	ColourCurls® maxi	dark red	6 mm	500 pcs
T	09646	ColourCurls® maxi	light yellow	6mm	500 pcs
U	09647	ColourCurls® maxi	orange	6mm	500 pcs
V	09648	ColourCurls® maxi	brown	6mm	500 pcs
W	09649	ColourCurls® maxi	black	6mm	500 pcs
X	09650	ColourCurls® maxi	light violet	6mm	500 pcs
Y	09651	ColourCurls® maxi	turquoise	6mm	500 pcs
Z	09695	ColourCurls® maxi	grey	6mm	500 pcs
AA	09664	ColourCurls® supermaxi	dark blue	7mm	100 pcs
BB	09665	ColourCurls® supermaxi	light green	7mm	100 pcs
CC	09666	ColourCurls® supermaxi	dark green	7mm	100 pcs
DD	09667	ColourCurls® supermaxi	white	7mm	100 pcs
EE	09668	ColourCurls® supermaxi	pink	7mm	100 pcs
FF	09669	ColourCurls® supermaxi	dark red	7mm	100 pcs
GG	09670	ColourCurls® supermaxi	light yellow	7mm	100 pcs
HH	09671	ColourCurls® supermaxi	orange	7mm	100 pcs
II	09672	ColourCurls® supermaxi	brown	7mm	100 pcs
JJ	09673	ColourCurls® supermaxi	black	7mm	100 pcs
KK	09675	ColourCurls® supermaxi	turquoise	7mm	100 pcs

Instrument Identification – Colour Coding



FIXING RINGS

Made of plastic for attaching tray labels. Heat-resistant up to 134°C. Single use.

	Product code	Colour	Diameter	Qty
A	09423	orange	10mm	100 pcs
B	09424	violet	10mm	100 pcs
C	09425	black	10mm	100 pcs
D	09426	transparent	10mm	100 pcs
E	09596	dark blue	10mm	100 pcs
F	09598	ochre yellow	10mm	100 pcs
G	09585	pink	10mm	100 pcs
H	09595	light green	10mm	100 pcs
I	09597	brown	10mm	100 pcs
J	09418	blue	10mm	100 pcs
K	09422	yellow	10mm	100 pcs

K

Instrument Identification

INSTRUMENT PROTECTION

Instrument Tip Caps are designed to extend the life of your instruments. Our tip caps protect your valuable instruments from unnecessary wear or damage, helping reduce repair costs while giving you the peace of mind that your instrument investment is being protected. In addition to protecting and extending the life of your valuable instruments, tip caps also help to protect sterilisation packs from puncture and protect employees from sharp instruments.

Available in a variety of shapes and sizes to fit all your instrument protection needs, we have the perfect tip cap for each of your various instruments, from flat caps to twin caps to instrument holders, and they are available with or without vents, depending on your preference. Our tip caps are radiopaque, single use, and compatible with steam sterilisation.

Our Instrument Holders are designed to protect various delicate instruments using an easy-to-fold design. Instrument Sleeves are single use, and compatible with steam and gas plasma sterilisation.

Flat Claps

Round Caps

Osteotome Caps

Retractor Caps

Twin Caps

MIS Caps

Skin Hook Cap

Tapered Caps

Instrument Holders

Tube and Cable Binders

Protection Boards

Please note: Unless specified, all products are not made with natural rubber latex.





FLAT CAPS

Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Description	Colour	Size	Qty
A	3-01-07	non-vented	orange	9.5 x 25.4mm	100 pcs
B	3-02-07	non-vented	orange tint	9.5 x 25.4mm	100 pcs
C	3-03-07	non-vented	clear	9.5 x 25.4mm	100 pcs
D	3-04-07	with vents	clear	9.5 x 25.4mm	100 pcs
E	3-08-07	with vents	orange tint	9.5 x 25.4mm	100 pcs
F	3-15-07	with vents	orange	9.5 x 25.4mm	100 pcs
G	3-01-08	non-vented	brown	15.9 x 25.4mm	100 pcs
H	3-02-08	non-vented	brown tint	15.9 x 25.4mm	100 pcs
I	3-03-08	non-vented	clear	15.9 x 25.4mm	100 pcs
J	3-04-08	with vents	clear	15.9 x 25.4mm	100 pcs
K	3-08-08	with vents	brown tint	15.9 x 25.4mm	100 pcs
L	3-15-08	with vents	brown	15.9 x 25.4mm	100 pcs
M	3-03-09	non-vented	clear	25.4 x 25.4mm	100 pcs
N	3-01-09	non-vented	black	25.4 x 25.4mm	100 pcs
O	3-02-09	non-vented	black tint	25.4 x 25.4mm	100 pcs
P	3-04-09	with vents	clear	25.4 x 25.4mm	100 pcs
Q	3-15-09	with vents	black	25.4 x 25.4mm	100 pcs
R	3-08-09	with vents	black tint	25.4 x 25.4mm	100 pcs

PRODUCTS SHOWN ACTUAL SIZE

Instrument Protection – Flat Caps



PRODUCTS SHOWN ACTUAL SIZE

PRECISION FLAT CAPS


For a more precise fit for enhanced instrument protection. With inner nubs. Vented. Single use.

Product code	Description	Colour	Size	Qty
A 30000	with vents, with raised inner nubs	orange tint	1.17 x 9.56 x 25.4mm	100 pcs
B 30001	with vents, with raised inner nubs	clear	3.18 x 25.4 x 25.4mm	100 pcs
C 30002	with vents, with raised inner nubs	black	3.18 x 25.4 x 25.4mm	100 pcs
D 30005	with vents	clear	1.57 x 9.53 x 25.4mm	100 pcs
E 30016	with vents	clear	1.98 x 17.8 x 25.4mm	100 pcs

FLAT CAP ASSORTMENT

Product code	Includes one of each of the following flat caps
F 3-08-100	3-08-01, 3-08-02, 3-08-03, 3-08-04, 3-08-05, 3-08-06, 3-08-07, 3-08-08, 3-08-09, 3-08-11, 3-08-12, 3-08-13, 3-08-15, 3-08-17, 3-08-18, 3-08-19

Instrument Protection – Round Caps



PRODUCTS SHOWN ACTUAL SIZE

ROUND CAPS

Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Description	Colour	Size	Qty
A	3-01-01	non-vented	white	1.6 x 19.1mm	100 pcs
	3-01-01S	N/A	black	1.6 x 20.3mm	100 pcs
B	3-02-01	non-vented	white tint	1.6 x 19.1mm	100 pcs
C	3-03-01	non-vented	clear	1.6 x 19.1mm	100 pcs
D	3-04-01	with vents	clear	1.6 x 19.1mm	100 pcs
E	3-08-01	with vents	white tint	1.6 x 19.1mm	100 pcs
F	3-15-01	with vents	white	1.6 x 19.1mm	100 pcs
	3-48-01	with vents	red tint	1.6 x 19.1mm	100 pcs
	3-60-01	pvc plastisol (Ia-6120-80)	black	2.0 x 5.92mm	100 pcs
G	3-01-02	non-vented	blue	2.0 x 19.1mm	100 pcs
H	3-02-02	non-vented	blue tint	2.0 x 19.1mm	100 pcs
I	3-03-02	non-vented	clear	2.0 x 19.1mm	100 pcs
J	3-04-02	with vents	clear	2.0 x 19.1mm	100 pcs
K	3-08-02	with vents	blue tint	2.0 x 19.1mm	100 pcs
L	3-15-02	with vents	blue	2.0 x 19.1mm	100 pcs
M	3-01-03	non-vented	green	2.8 x 19.1mm	100 pcs
N	3-02-03	non-vented	green tint	2.8 x 19.1mm	100 pcs
O	3-03-03	non-vented	clear	2.8 x 19.1mm	100 pcs
P	3-04-03	with vents	clear	2.8 x 19.1mm	100 pcs
Q	3-08-03	with vents	green tint	2.8 x 19.1mm	100 pcs
R	3-15-03	with vents	green	2.8 x 19.1mm	100 pcs
S	3-01-04	non-vented	red	3.2 x 25.4mm	100 pcs
T	3-02-04	non-vented	red tint	3.2 x 25.4mm	100 pcs
U	3-03-04	non-vented	clear	3.2 x 25.4mm	100 pcs
V	3-04-04	with vents	clear	3.2 x 25.4mm	100 pcs
W	3-08-04	with vents	tinted red	3.2 x 25.4mm	100 pcs
X	3-15-04	with vents	red	3.2 x 25.4mm	100 pcs
Y	3-01-05	non-vented	yellow	4.7 x 25.4mm	100 pcs
Z	3-02-05	non-vented	yellow tint	4.7 x 25.4mm	100 pcs

















K

Instrument Protection

Instrument Protection – Round Caps








ROUND CAPS

Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Description	Colour	Size	Qty
	3-03-05	non-vented	clear	4.7 x 25.4mm	100 pcs
	3-04-05	with vents	clear	4.7 x 25.4mm	100 pcs
	3-08-05	with vents	yellow tint	4.7 x 25.4mm	100 pcs
	3-15-05	with vents	yellow	4.7 x 25.4mm	100 pcs
	3-01-06	non-vented	purple	10.3 x 19.1mm	100 pcs
	3-02-06	non-vented	purple tint	10.3 x 19.1mm	100 pcs
	3-03-06	N/A	clear	10.3 x 19.1mm	100 pcs
	3-04-06	with vents	clear	10.3 x 19.1mm	100 pcs
	3-08-06	with vents	purple tint	10.3 x 19.1mm	100 pcs
	3-15-06	with vents	purple	10.3 x 19.1mm	100 pcs
	3-05-27	non-vented	red	7.9 x 38.1mm	100 pcs
	3-05-28	non-vented	navy	11.2 x 25.4mm	100 pcs
	3-05-29	non-vented	navy	6.4 x 31.8mm	100 pcs
	3-05-30	non-vented	navy	4.0 x 25.4mm	100 pcs
	3-04-98	with vents	clear	3.9 x 18.7mm	100 pcs
	3-04-99	with vents	clear	4.7 x 38.1mm	100 pcs

PRECISION ROUND CAPS

For a more precise fit for enhanced instrument protection. With inner nubs. Vented. Single use.

	Product code	Description	Colour	Size	Qty
	30006	with vents	tinted white	1.6 x 19.1mm	100 pcs
	30007	with vents	blue	2 x 19.1mm	100 pcs
	30008	with vents	green	2.8 x 19.1mm	100 pcs
	30009	with vents	red	3.2 x 25.4mm	100 pcs
	30010	with vents	yellow	4.8 x 25.4mm	100 pcs
	30012	with vents	clear	2 x 19.1mm	100 pcs
	30013	with vents	clear	2 x 19.1mm	100 pcs

PRODUCTS SHOWN ACTUAL SIZE

K

Instrument Protection

Instrument Protection – Osteotome Caps



PRODUCTS SHOWN ACTUAL SIZE

OSTEOTOME CAPS

Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Description	Colour	Size	Qty
A	3-03-11	non-vented	clear	6.4 x 19.1mm	100 pcs
B	3-04-11	with vents	clear	6.4 x 19.1mm	100 pcs
C	3-08-11	with vents	white tint	6.4 x 19.1mm	100 pcs
D	3-03-12	non-vented	clear	9.5 x 25.4mm	100 pcs
E	3-04-12	with vents	clear	9.5 x 25.4mm	100 pcs
F	3-08-12	with vents	blue tint	17.8 x 25.4mm	100 pcs
G	3-03-13	non-vented	clear	17.8 x 25.4mm	100 pcs
H	3-04-13	with vents	clear	17.8 x 25.4mm	100 pcs
I	3-08-13	with vents	green tint	17.8 x 25.4mm	100 pcs
J	3-03-15	non-vented	clear	17.8 x 25.4mm	100 pcs
K	3-04-15	with vents	clear	17.8 x 25.4mm	100 pcs
L	3-08-15	with vents	yellow tint	17.8 x 25.4mm	100 pcs
M	3-03-09	non-vented	clear	25.4 x 25.4mm	100 pcs
N	3-04-09	with vents	clear	25.4 x 25.4mm	100 pcs
O	3-08-09	with vents	black tint	25.4 x 25.4mm	100 pcs
P	3-04-17	with vents	clear	31.8 x 31.8mm	100 pcs
Q	3-08-17	with vents	purple tint	31.8 x 31.8mm	100 pcs
R	3-04-18	with vents	clear	38.1 x 38.1mm	100 pcs
S	3-08-18	with vents	red tint	38.1 x 38.1mm	100 pcs
T	3-04-19	with vents	clear	44.5 x 31.8mm	100 pcs
U	3-08-19	with vents	orange tint	44.5 x 31.8mm	100 pcs

PRECISION OSTEOTOME CAPS

For a more precise fit for enhanced instrument protection. With inner nubs. Single use.

	Product code	Description	Colour	Size	Qty
V	30004	with vents	clear	6.4 x 19.1mm	50 pcs

V



Instrument Protection – Retractor Caps



PRODUCTS SHOWN ACTUAL SIZE

SELF-RETAINING RETRACTOR CAPS


Protect a variety of self-retaining retractors including Weitlaner, D’Errico-Adson, Beckman, and Adson retractors. Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Colour	Size	Qty
A	01472	purple tint	25 x 13mm	50 pcs
B	01473	clear	25 x 13mm	50 pcs
C	01474	green tint	32 x 13mm	50 pcs
D	01475	clear	32 x 13mm	50 pcs
E	3-09-07	blue tint	21.6 x 20.3mm	100 pcs
F	3-09-08	purple tint	25.4 x 25.4mm	100 pcs
G	3-09-09	orange tint	25.4 x 31.8mm	100 pcs
H	3-09-10	green tint	25.4 x 44.5mm	100 pcs

GELPI RETRACTOR CAPS

Protect a variety of Gelpi retractors. Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Colour	Size	Qty
I	3-09-03	clear	7.9 x 38.1mm	100 pcs
J	3-09-04	clear	7.9 x 50.8mm	100 pcs
K	3-09-05	clear	7.9 x 63.5mm	100 pcs





PRODUCTS SHOWN ACTUAL SIZE

Instrument Protection – Twin Caps, MIS Caps and Skin Hook Cap

TWIN CAPS

Protect instruments while keeping them in an open position for cleaning. Perfect for scissors and forceps. Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Description	Colour	Size	Qty
A	3-06-01	with vents	white tint	1.6 x 20.3mm	50 pcs
B	3-06-02	with vents	blue tint	2.0 x 21.1mm	50 pcs
C	3-06-03	with vents	green tint	2.8 x 21.7mm	50 pcs
D	3-06-04	with vents	red tint	3.2 x 27.8mm	50 pcs
E	3-06-05	with vents	yellow tint	4.7 x 28.2mm	50 pcs
F	3-06-07	with vents	orange tint	27.6 x 19.1mm	50 pcs

TWIN CAP ASSORTMENT

Product code Includes one of each of the following twin caps

3-06-50 3-06-01, 3-06-02, 3-06-03, 3-06-04, 3-06-05, 3-06-07

MINIMALLY INVASIVE SURGICAL INSTRUMENT (MIS) CAPS

Instrument Caps are designed to protect endoscopes, probes, veress needles, and trocars during sterilisation. Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Description	Colour	Size	Qty
G	3-07-03	with vents	clear	4.3 x 28.7mm	50 pcs
H	3-07-05	with vents	clear	4.3 x 50.8mm	50 pcs
I	3-07-09	with vents	clear	9.1 x 66.0mm	50 pcs

SKIN HOOK CAP

Protect delicate single and double prong skin hook instrumentation. Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Description	Colour	Size	Qty
J	3-10-02	single prong, 5mm hook	clear	6.3 x 19.1mm	100 pcs



Instrument Protection – Tapered Caps



TAPERED CAPS

Radiopaque. Compatible with steam sterilisation. Single use.

	Product code	Description	Colour	Diameter	Length	Qty
A	3-11-01	with vents	white tint	1.5mm to 3.3mm	38.1mm	50 pcs
B	3-11-02	with vents	blue tint	2.5mm to 4.3mm	38.1mm	50 pcs
C	3-11-03	with vents	green tint	3.5mm to 5.3mm	38.1mm	50 pcs
D	3-11-04	with vents	red tint	4.5mm to 6.2mm	38.1mm	50 pcs
E	3-11-05	with vents	yellow tint	5.5mm to 7.2mm	38.1mm	50 pcs
F	3-11-06	with vents	purple tint	6.5mm to 8.2mm	38.1mm	50 pcs
G	3-11-07	with vents	orange tint	7.5mm to 9.2mm	38.1mm	50 pcs
H	3-11-08	with vents	brown tint	8.5mm to 10.2mm	38.1mm	50 pcs
I	3-11-09	with vents	black tint	9.4mm to 11.2mm	38.1mm	50 pcs
J	3-12-01	with vents	white tint	2mm to 3.1mm	52.1mm	50 pcs
K	3-12-02	with vents	blue tint	2.5mm to 4.9mm	52.3mm	50 pcs
L	3-12-03	with vents	green tint	3.5mm to 5.8mm	53.3mm	50 pcs
M	3-12-04	with vents	red tint	4.5mm to 6.9mm	53.3mm	50 pcs
N	3-12-05	with vents	yellow tint	5.5mm to 7.8mm	53.3mm	50 pcs
O	3-12-06	with vents	purple tint	6.5mm to 8.9mm	50.8mm	50 pcs
P	3-12-07	with vents	orange tint	7.5mm to 9.8mm	50.8mm	50 pcs
Q	3-12-08	with vents	brown tint	8.5mm to 10.8mm	50.8mm	50 pcs
R	3-12-09	with vents	black tint	9.4mm to 11.8mm	50.8mm	50 pcs
S	3-13-03	with vents	green tint	13.4mm to 14.7mm	38.1mm	50 pcs
T	3-13-05	with vents	yellow tint	15.4mm to 16.8mm	38.1mm	50 pcs
U	3-13-07	with vents	orange tint	17.4mm to 18.8mm	38.1mm	50 pcs
V	3-14-03	with vents	green tint	13.4mm to 15.2mm	50.8mm	50 pcs
W	3-14-05	with vents	yellow tint	15.4mm to 17.2mm	50.8mm	50 pcs
X	3-14-07	with vents	orange tint	17.4mm to 19.3mm	50.8mm	50 pcs
Y	3-05-02	non-vented	blue tint	6.35mm	38.1mm	50 pcs
Z	3-05-06	non-vented	purple tint	10.16mm	48.3mm	50 pcs
AA	3-05-07	non-vented	orange tint	11.43mm	44.5mm	50 pcs

K

Instrument Protection



PRODUCTS SHOWN ACTUAL SIZE

Instrument Protection – Instrument Holders


INSTRUMENT SLEEVES

Easy-to-fold protective sleeves conveniently protect your various delicate instruments during sterilisation. Available in three versatile sizes. Compatible with steam and gas plasma sterilisation. Single use.

Product code	Description	Size	Qty
A CS-1010W	white paper	25.4mm	1000 pcs
B CS-2025W	white paper	50.8mm	1000 pcs
C CS-3050W	white paper	76.2mm	1000 pcs



Instrument Protection – Instrument Holders



WEBMESH GUARD

Instrument protector for low heat and plasma sterilisation. Made of Polyethylene.

	Product code	Colour	Size	Qty
A	01851	green	15 x 76mm	24 pcs
B	01852	blue	20 x 76mm	24 pcs
C	01853	red	40 x 76mm	24 pcs

POINT PROTECTORS

Easy-to-fold protective sleeves conveniently protect your various delicate instruments during sterilisation. Available in three versatile sizes. Compatible with steam sterilisation. Single use.

	Product code	Description	Size	Qty
D	PPT001	white card	160 x 50mm	2000/box
	PPT002	white card	160 x 50mm	1000/box
	PPT003	white card	160 x 50mm	2000/box
	PPT004	white card	160 x 50mm	1000/box





Instrument Protection – Instrument Holders

TIP PROTECTORS

To protect high-quality instruments from damage during transportation, handling and sterilisation. However, the shape of the protectors is not adapted to the specific instrument. Heat-resistant up to 134°C. Single use.

	Product code	Material	Size	Qty
	01469	cardboard	105 x 25mm	1000 pcs
A	01470	cardboard	60 x 10mm	1000 pcs
B	01471	transparent film	105 x 25mm	1000 pcs
C	01476	transparent film	60 x 10mm	1000 pcs
D	01477	transparent film	140 x 50mm	100 pcs
	01492	carton	140 x 50mm	100 pcs
	01499	carton	105 x 20mm	100 pcs
E	01478	knife cover, with 16 holes	246 x 33mm	100 pcs
F	01479	knife cover, with 8 holes	123 x 33mm	100 pcs

NEEDLE PROTECTORS

Simple-to-assemble paper based product which are designed to protect chisel and other sharp instruments during sterilisation, transport and storage. Single use.

	Product code	Description	Size	Qty
	PPT005	white card	160 x 50mm	2000/box
G	PPT006	white card	160 x 50mm	1000/box

CLEAR POINT PROTECTORS

Simple-to-assemble plastic based product which are designed to protect chisels and other sharp instruments during sterilisation, transport and storage. Ideal for use with products that require identification. Single use.

	Product code	Description	Size	Qty
H	PPT007C	clear plastic	160 x 50mm	200/box
	PPT008C	clear plastic	160 x 50mm	200/box

Instrument Protection – Tube and Cable Binders


SILICONE TUBES

Silicone tubes can be used to bundle different instruments. Length 25m. Single use.

Product code	Colour	Size	Qty
09707	clear	1.0 x 2.0mm	25 pcs
09024	clear	2.0 x 1.0mm	25 pcs
09706	clear	2.5 x 0.5mm	25 pcs
09183	clear	2.5 x 1.2mm	25 pcs
09149	clear	2.5 x 1.5mm	25 pcs
09233	clear	4.0 x 0.5mm	25 pcs
09023	clear	4.0 x 1.0mm	25 pcs
09150	clear	4.0 x 1.2mm	25 pcs
09115	clear	4.0 x 1.5mm	25 pcs
09151	clear	6.0 x 1.2mm	25 pcs
09193	clear	6.0 x 2.5mm	25 pcs
A 09018	clear	7.0 x 2.0mm	25 pcs
09116	clear	7.0 x 3.0mm	25 pcs
09152	clear	8.0 x 1.2mm	25 pcs
09701	clear	8.0 x 2.5mm	25 pcs
09700	clear	10.0 x 2.5mm	25 pcs
09142	clear	10.6 x 1.2mm	25 pcs
09153	clear	12.0 x 1.2mm	25 pcs
09705	clear	20.0 x 1.0mm	25 pcs
B 09119	red	5.0 x 1.5mm	25 pcs
09165	green	2.0 x 1.0mm	25 pcs
C 09710	green	2.5 x 0.5mm	25 pcs
09711	green	4.0 x 0.5mm	25 pcs
09712	green	6.0 x 1.2mm	25 pcs
09713	green	8.0 x 1.2mm	25 pcs
D 09727	white	2.0 x 1.0mm	25 pcs
09166	blue	2.0 x 1.0mm	25 pcs
09154	blue	3.0 x 2.0mm	25 pcs
09127	blue	5.0 x 1.5mm	25 pcs
09136	orange	5.0 x 1.5mm	25 pcs



Instrument Protection – Tube and Cable Binders



SILICONE CABLES

For bundling instruments, heat-resistant up to 134°C. Single use.

	Product code	Colour	Size	Qty
A	09235	red	110mm	100 pcs
B	09236	yellow	110mm	100 pcs
C	09237	green	110mm	100 pcs
	09240	yellow	210mm	100 pcs
	09241	green	210mm	100 pcs
	09242	blue	210mm	100 pcs
	09252	green	60mm	100 pcs
D	09238	blue	110mm	100 pcs
E	09285	blue	190mm	100 pcs
F	09586	yellow	210mm	100 pcs

TUBE AND CABLE BINDERS

Tube and cable binders for bundling instruments. Single use.

	Product code	Colour	Size	Qty
G	09158	white	25mm	100 pcs
H	09159	black	25mm	100 pcs

SILICONE STRING

String for bundling instrument. Heat-resistant up to 134°C. Single use.

	Product code	Colour	Length	Qty
	09134	clear	25m	1 pc
I	09161	yellow	25m	1 pc





Instrument Protection – Protection Boards

SHEETS

Silicone sheets for packing tables protect instruments from knocks and also aid handling thanks to their non-slip properties.

	Product code	Colour	Thickness	Qty
A	08700	transparent	2mm	1 pc
B	08705	red-brown	2mm	1 pc
C	08703	blue	2mm	1 pc
D	08701	leaf-green	2mm	1 pc

SILICONE BOARDS

The silicone boards protect the instruments, make handling easier thanks to their anti-slip properties and take away the unpleasant cold, glare and noise effects of metal surfaces from the staff. Uncoloured.

	Product code	Size	Qty
	09189	300 x 170 x 10mm	1 pc
E	09190	170 x 150 x 10mm	1 pc
F	09192	210 x 150 x 10mm	1 pc
G	09191	300 x 210 x 10mm	1 pc
	09194	50 x 50 x 5mm	1 pc
	09195	1000 x 2300 x 5mm	1 pc
	09198	450 x 530 x 5mm	1 pc



TAGS

Locking Tags are designed to help you identify your instruments, trays, and trolleys as they move from sterilisation to staging areas to procedures. Our tags allow you to indicate special handling instructions, secure trays and carts from tampering and quickly identify priority instruments for repair and sterilisation. Tags are available in a wide range of styles, colours, and locking configurations. Printed tags with special handling instructions, as well as uniquely numbered tags, are available as a standard product. Need a customised tag? We can do that for you!

Breakaway Tags are also available in a range of colours and can be written on with a permanent lab marker. They must be cut or broken to open, providing a secure lock.

Write-On Tags allow you to write specific handling instructions on each tag when more detailed communication is necessary. These tags attach to instruments and trays with ties and are larger and more visible than our other tags, ensuring that detailed instructions will be seen by all.

Metal Tray Tags easily clip onto the edge of a tray and provide ample space for identification, tracking and barcoding needs. Stainless steel for durability.

Plastic Tray Tags

Breakaway Tags

Locking Tags

One-Loop Tags

Self-Tie Tags

Identification Tags

Washer Tags


Write-On Tags

Metal Tray Tags

Please note: Unless specified, all products are not made with natural rubber latex.



Tags – Plastic Tray Tags



WHITE TAGS SHOWN ACTUAL SIZE

PLASTIC TRAY TAGS

Durable, plastic tags designed for use on instrument trays. Tray tags provide a flat surface for barcode placement. Great for colour coding. Compatible with steam and gas plasma sterilisation. Tags include a built-in clip on the back.

	Product code	Colour	Size	Qty
A	PT-01-01	white	30 x 30mm	10 pcs
B	PT-01-02	blue	30 x 30mm	10 pcs
C	PT-01-03	green	30 x 30mm	10 pcs
D	PT-01-04	red	30 x 30mm	10 pcs
E	PT-01-05	yellow	30 x 30mm	10 pcs
F	PT-01-06	purple	30 x 30mm	10 pcs
G	PT-01-07	orange	30 x 30mm	10 pcs
H	PT-01-09	black	30 x 30mm	10 pcs
I	PT-01-14	light blue	30 x 30mm	10 pcs
J	PT-02-01	white	55 x 30mm	10 pcs
K	PT-02-02	blue	55 x 30mm	10 pcs
L	PT-02-03	green	55 x 30mm	10 pcs
M	PT-02-04	red	55 x 30mm	10 pcs
N	PT-02-05	yellow	55 x 30mm	10 pcs
O	PT-02-06	purple	55 x 30mm	10 pcs
P	PT-02-07	orange	55 x 30mm	10 pcs
Q	PT-02-09	black	55 x 30mm	10 pcs
R	PT-02-14	light blue	55 x 30mm	10 pcs



WHITE TAGS SHOWN ACTUAL SIZE

Tags – Plastic Tray Tags

PLASTIC TRAY TAGS

Durable, plastic tags designed for use on instrument trays. Tray tags provide a flat surface for barcode placement. Great for colour coding. Compatible with steam and gas plasma sterilisation. Tags include a built-in clip on the back.

	Product code	Colour	Size	Qty
A	PT-04-01	white	70 x 30mm	10 pcs
B	PT-04-02	blue	70 x 30mm	10 pcs
C	PT-04-03	green	70 x 30mm	10 pcs
D	PT-04-04	red	70 x 30mm	10 pcs
E	PT-04-05	yellow	70 x 30mm	10 pcs
F	PT-04-06	purple	70 x 30mm	10 pcs
G	PT-04-07	orange	70 x 30mm	10 pcs
H	PT-04-09	black	70 x 30mm	10 pcs
I	PT-04-14	light blue	70 x 30mm	10 pcs
J	PT-05-01	white	90 x 30mm	10 pcs
K	PT-05-02	blue	90 x 30mm	10 pcs
L	PT-05-03	green	90 x 30mm	10 pcs
M	PT-05-04	red	90 x 30mm	10 pcs
N	PT-05-05	yellow	90 x 30mm	10 pcs
O	PT-05-06	purple	90 x 30mm	10 pcs
P	PT-05-07	orange	90 x 30mm	10 pcs
Q	PT-05-09	black	90 x 30mm	10 pcs
R	PT-05-14	light blue	90 x 30mm	10 pcs

K

Tags



TAGS SHOWN ACTUAL SIZE

PLASTIC TRAY TAGS

Durable, plastic tags designed for use on instrument trays. Tray tags provide a flat surface for barcode placement. Great for colour coding. Compatible with steam and gas plasma sterilisation. Tags designed to be attached to a tray using a fastener such as a ring or clip.

Product code	Colour	Size	Qty
A PT-03-01	white	60 x 40mm	10 pcs
B PT-03-02	blue	60 x 40mm	10 pcs
C PT-03-03	green	60 x 40mm	10 pcs
D PT-03-04	red	60 x 40mm	10 pcs
E PT-03-05	yellow	60 x 40mm	10 pcs
F PT-03-06	purple	60 x 40mm	10 pcs
G PT-03-07	orange	60 x 40mm	10 pcs
H PT-03-09	black	60 x 40mm	10 pcs
I PT-03-14	light blue	60 x 40mm	10 pcs
J PT-06-01	white	90 x 33mm	10 pcs
K PT-06-02	blue	90 x 33mm	10 pcs
L PT-06-03	green	90 x 33mm	10 pcs
M PT-06-04	red	90 x 33mm	10 pcs
N PT-06-05	yellow	90 x 33mm	10 pcs
O PT-06-06	purple	90 x 33mm	10 pcs
P PT-06-07	orange	90 x 33mm	10 pcs
Q PT-06-09	black	90 x 33mm	10 pcs
R PT-06-14	light blue	90 x 33mm	10 pcs

TRAY TAG RING AND CLIP

Durable, stainless steel ring and clip are designed to be used with Plastic Tray Tags to attach the tag to the tray. The clip specifically fits Plastic Tray Tags PT-06-01 through to PT-06-14.

Product code	Decription	Qty
S TG-100	ring, 20mm diameter	100 pcs
T TG-200	clip	100 pcs



PLASTIC TRAY TAGS



Tray tags simplify the identification of sets or instruments and provide a useful means of labelling items for tracking and tracing. All tags are produced to withstand the rigours of the decontamination cycle. Available in 18 colours and also in stainless steel. Tray tags can be attached to trays using either a split ring or a cable tie.

Product code	Colour	Size	Qty
A TAG020	black	90 x 38mm	100 pcs
B TAG018	grey	90 x 38mm	100 pcs
C TAG027	sky blue	90 x 38mm	100 pcs
D TAG026	royal blue	90 x 38mm	100 pcs
E TAG023	turquoise	90 x 38mm	100 pcs
F TAG017	green	90 x 38mm	100 pcs
G TAG021	lime	90 x 38mm	100 pcs
H TAG013	yellow	90 x 38mm	100 pcs
I TAG014	orange	90 x 38mm	100 pcs
J TAG015	red	90 x 38mm	100 pcs
K TAG012	tan	90 x 38mm	100 pcs
L TAG011	white	90 x 38mm	100 pcs
TAG016	blue	90 x 38mm	100 pcs
TAG019	purple	90 x 38mm	100 pcs
TAG022	pink	90 x 38mm	100 pcs
TAG024	mustard	90 x 38mm	100 pcs
TAG025	terracotta	90 x 38mm	100 pcs
TAG55WH	white	90 x 55mm	100 pcs



WHITE TAG SHOWN ACTUAL SIZE








PRODUCTS SHOWN 60% ACTUAL SIZE

Tags – Breakaway Tags

BREAKAWAY TAGS

Versatile, durable, and economical. Tamper-resistant. Tags cinch closed until tight and must be cut or broken to open, providing a secure lock. Select styles are compatible with steam and gas plasma sterilisation. Single use.

	Product code	Description	Colour	Size	Qty
A	LR-102	sterilisable	blue	89mm	100 pcs
B	LR-103	sterilisable	green	89mm	100 pcs
C	LR-104	sterilisable	red	89mm	100 pcs
D	LR-105	sterilisable	yellow	89mm	100 pcs
E	LR-103N	sterilisable, uniquely numbered	green	89mm	100 pcs
F	LR-105N	sterilisable, uniquely numbered	yellow	89mm	100 pcs
G	LR-312	sterilisable	blue	139.7mm	100 pcs
H	LR-313	sterilisable	green	139.7mm	100 pcs
I	LR-314	sterilisable	red	139.7mm	100 pcs
J	LR-315	sterilisable	yellow	139.7mm	100 pcs
K	LR-322	non-sterilisable, uniquely numbered	blue	152.4mm	100 pcs
L	LR-323	non-sterilisable, uniquely numbered	green	152.4mm	100 pcs
M	LR-324	non-sterilisable, uniquely numbered	red	152.4mm	100 pcs
N	LR-325	non-sterilisable, uniquely numbered	yellow	152.4mm	100 pcs
O	LR-327	non-sterilisable, uniquely numbered	orange	152.4mm	100 pcs
P	LR-329	non-sterilisable, uniquely numbered	black	152.4mm	100 pcs



PRODUCTS SHOWN ACTUAL SIZE

LOCKING TAGS

Versatile, durable, and economical. Tamper-resistant. Compatible with steam and gas plasma sterilisation. Customisable. 25.4mm length. Single use. Quantities available in 100 pcs and 1000 pcs.



	Product code	Description	Colour	Qty
A	LT-150	clean	white	100 pcs
B	LT-210	quick turnover	blue	100 pcs
C	LT-310	loaner	green	100 pcs
D	LT-450	repair	red	100 pcs
E	LT-550	sharpen	yellow	100 pcs
F	LT-720	broken	orange	100 pcs
G	LOCK994G	good	green	100 pcs
H	LT-100N	uniquely numbered	white	100 pcs
I	LT-300N	uniquely numbered	green	100 pcs
J	LT-400N	uniquely numbered	red	100 pcs
K	LT-500N	uniquely numbered	yellow	100 pcs
L	LT-900N	uniquely numbered	black	100 pcs

SNAPLOCK LOCKING TAGS

Our range of non-repetitive numbered locks for full product tracking and traceability. Compatible with steam sterilisation systems. Single use.

	Product code	Description	Colour	Qty
	SNAPLOCKG001	uniquely numbered	green	100 pcs
M	SNAPLOCKG002	no numbers	green	100 pcs
	SNAPLOCKR002	uniquely numbered	red	100 pcs
	SNAPLOCKY003	uniquely numbered	yellow	100 pcs





PRODUCTS SHOWN ACTUAL SIZE

LOCKING TAGS

Versatile, durable, and economical. Tamper-resistant. Compatible with steam and gas plasma sterilisation. Customisable. 25.4mm length. Single use. Quantities available in 100 pcs and 1000 pcs.

Product code	Colour	Qty
A LT-100	white	100 pcs
B LT-110	pink	100 pcs
C LT-112	gold	100 pcs
D LT-116	purple	100 pcs
E LT-200	blue	100 pcs
F LT-300	green	100 pcs
G LT-400	red	100 pcs
H LT-500	yellow	100 pcs
I LT-600	grey	100 pcs
J LT-700	orange	100 pcs
K LT-800	brown	100 pcs
L LT-900	black	100 pcs

LOCKING TAGS

Plastic seals for Aesculap containers. Heat resistant up to 134°C. Single use.

Product code	Colour	Qty
M 01020	blue	1000 pcs
N 01023	red	1000 pcs
O 01024	yellow	1000 pcs
P 02623	light blue	1000 pcs
Q 01241	green	1000 pcs
R 01243	red	1000 pcs
S 01254	yellow	1000 pcs





CASE LOCKING TAGS

A versatile, durable, and economical option for identifying and securing cases and containers. Tamper-resistant. Compatible with steam and gas plasma sterilisation. 22mm operational length. Smallest size slot tag will fit is 1.6mm x 10mm. Single use. Quantities available in 100 pcs and 1000 pcs.

	Product code	Colour	Qty
A	LT-001	white	100 pcs
B	LT-002	blue	100 pcs
C	LT-003	green	100 pcs
D	LT-004	red	100 pcs
E	LT-007	orange	100 pcs



INDICATOR DOT LOCKING TAGS

Indicator dots on tags change colour to indicate exposure to steam sterilisation (blue dot to dark). Versatile, durable, and economical. Tamper-resistant. Compatible with steam sterilisation. 25.4mm length. Single use. Quantities available in 100 pcs and 1000 pcs.


	Product code	Colour	Qty
F	LS-007	orange	100 pcs
G	LS-100	white	100 pcs
H	LS-400	red	100 pcs




INDICATOR DOT LOCKING TAGS

Rectangular bars with universal seals. Compatible with steam sterilisation. Single use.

	Product code	Colour	Qty
I	01013	blue	1000 pcs
J	01212	red	1000 pcs
K	01213	green	1000 pcs
L	01214	yellow	1000 pcs
M	01215	black	1000 pcs





INDICATOR DOT LOCKING TAGS


Suitable for Aesculap and Wagner containers amongst others. Rectangular bar with '2011' universal seals. Compatible with steam sterilisation. Single use.

	Product code	Colour	Qty
A	01220	dark blue	1000 pcs
B	01221	red	1000 pcs
C	01233	green	1000 pcs
D	01238	yellow	1000 pcs
E	01239	white	1000 pcs
F	02640	orange	1000 pcs
G	02507	blue	1000 pcs
H	02510	red	1000 pcs
I	02511	yellow	1000 pcs

INDICATOR DOT LOCKING TAGS

For Aesculap containers. Rectangular bar with seals. Compatible with steam sterilisation. Single use.

	Product code	Colour	Qty
J	02670	light blue	1000 pcs
K	01228	blue	1000 pcs
L	01230	white	1000 pcs
M	01231	red	1000 pcs
N	01237	violet	1000 pcs





Tags – Locking Tags

INDICATOR DOT LOCKING TAGS

For Aesculap and Wagner containers. Compatible with steam sterilisation. *For Martin containers. Single use.

	Product code	Colour	Qty
A	02632	red	1000 pcs
B	02633	blue	1000 pcs
C	02637*	red	500 pcs

INDICATOR DOT LOCKING TAGS

Rectangular bar with 1-2-3 universal seals. Compatible with steam sterilisation. Single use.

	Product code	Colour	Qty
D	02601	natural colour	1000 pcs
E	02603	white	1000 pcs
F	02605	yellow	1000 pcs
G	02607	red	1000 pcs
H	02609	blue	1000 pcs
I	02611	green	1000 pcs

INDICATOR DOT LOCKING TAGS


For Aesculap and Wagner containers. Compatible with steam sterilisation. Single use.

	Product code	Colour	Qty
J	01245	blue	1000 pcs
K	01246	green	1000 pcs
L	01247	white	1000 pcs
M	01248	red	1000 pcs
N	01249	yellow	1000 pcs

INDICATOR DOT LOCKING TAGS

Rectangular XL bar with universal seal. Compatible with steam sterilisation. Single use.

	Product code	Colour	Qty
O	02516	violet	1000 pcs



Tags – One-Loop Tags

ONE-LOOP TAGS


Designed to loop around tray handles to identify through colour-coding or providing a flat surface for barcode placement. Compatible with steam and EtO sterilisation. Single use.

	Product code	Colour	Size	Qty
A	PL-101	white	206 x 110mm	250 pcs
B	PL-201	white	236 x 50mm	250 pcs
C	PL-207	orange	236 x 50mm	250 pcs

ONE-LOOP TAGS - MESH BASKETS

Heat resistant up to 134°C. Suitable for washer-disinfectors. Single use.

	Product code	Colour	Size	Qty
D	09244	white	236 x 50mm	250 pcs
E	09264	blue	236 x 50mm	250 pcs
F	09265	green	236 x 50mm	250 pcs
G	09267	yellow	236 x 50mm	250 pcs
H	09637	red	160 x 35mm	250 pcs
I	09677	black	160 x 35mm	250 pcs
	09399	white	254 x 42mm	600 pcs
J	09570	green, urgent	240 x 50mm	250 pcs
K	09584	yellow, urgent	240 x 50mm	250 pcs
L	09752	red, sharpen	150 x 35mm	250 pcs
M	09753	orange, repair	150 x 35mm	250 pcs





Tags – Self-Tie Tags

SELF-TIE TAGS

Self-Tie Tags are designed for the identification of trays that are to be steam sterilised. Plain tag with dumbbell self-tie. Single use.

	Product code	Thickness	Colour	Size	Qty
A	SELFTIE01	195 mic	white	25 x 190.5mm	1000 per roll
B	SELFTIE01-GRN	195 mic	green	25 x 190.5mm	1000 per roll
C	SELFTIE01-YEL	195 mic	yellow	25 x 190.5mm	1000 per roll
D	SELFTIE01-RED	195 mic	red	25 x 190.5mm	1000 per roll
E	SELFTIE01-BLU	195 mic	blue	25 x 190.5mm	1000 per roll
F	SELFTIE01-ORA	195 mic	orange	25 x 190.5mm	1000 per roll





IDENTIFICATION CONTAINER CARRIER TAGS

Identification tags comply with different tracking and tracing systems and sterilisation methods. Plain tag with dumbbell self-tie.

IDENTIFICATION TAGS

	Product code	Description	Colour	Size	Qty
A	TAG-CONTAIN1	175 mic	white	75 x 228.6mm	250 per roll
B	TAG-CONTAIN1-GRE	175 mic	green	75 x 228.6mm	250 per roll
C	TAG-CONTAIN1-RED	175 mic	red	75 x 228.6mm	250 per roll
D	TAG-CONTAIN1-YEL	175 mic	yellow	75 x 228.6mm	250 per roll
	TAG-CONTAIN2	175 mic	white	96 x 280mm	250 per roll

IDENTIFICATION SECURITY TAGS

	Product code	Description	Colour	Size	Qty
	TAG-SECURE1	175 mic	white	152 x 101.6mm	500 per roll
E	TAG-SECURE1-BLU	175 mic	blue	152 x 101.6mm	500 per roll
F	TAG-SECURE1-GRE	175 mic	green	152 x 101.6mm	500 per roll
G	TAG-SECURE1-ORA	175 mic	orange	152 x 101.6mm	500 per roll
H	TAG-SECURE1-RED	175 mic	red	152 x 101.6mm	500 per roll
	TAG-SECURE1-YEL	175 mic	yellow	152 x 101.6mm	500 per roll

IDENTIFICATION AUTOCLAVE TAGS


	Product code	Description	Colour	Size	Qty
	TAG-TRAY1	175 mic, tray tags, cable tie	white	70 x 101.6mm	500 per roll
K	TAG-TRAY1-BLU	175 mic, tray tags, cable tie	blue	70 x 101.6mm	500 per roll

Tags – Washer Tags

WASHER TAGS

Washer tags can be fixed onto the handles of wash baskets to share information with the clean room.

Product code	Description	Colour	Size	Qty
TAG-WASH	self-tie, t-shape, polyplas 250 mic	white	50 x 228.6mm	500 per roll
A TAG-WASH-BLU	self-tie, t-shape, polyplas 250 mic	blue	50 x 228.6mm	500 per roll
B TAG-WASH-GRN	self-tie, t-shape, polyplas 250 mic	green	50 x 228.6mm	500 per roll
C TAG-WASH-PNK	self-tie, t-shape, polyplas 250 mic	pink	50 x 228.6mm	500 per roll
D TAG-WASH-RED	self-tie, t-shape, polyplas 250 mic	red	50 x 228.6mm	500 per roll
E TAG-WASH-ORA	self-tie, t-shape, polyplas 250 mic	orange	50 x 228.6mm	500 per roll
F TAG-WASH-YEL	self-tie, t-shape, polyplas 250 mic	yellow	50 x 228.6mm	500 per roll





WRITE-ON TAGS

Durable and economical. Identify and record specific instructions with these versatile tags. Tamper-resistant. Compatible with steam and gas plasma sterilisation. Customisable. Single use.

	Product code	Description	Colour	Qty
A	WT-0009	blank tag	black	100 pcs
B	WT-0203	loaner	green	100 pcs
C	WT-0304	priority	red	100 pcs
D	WT-0405	re-run	yellow	100 pcs
E	WT-0507	broken, do not use	orange	100 pcs
F	WT-1911	missing	pink	100 pcs
G	WT-2614	loaner	light green	100 pcs

MARKER PENS

Ink is waterproof, alcohol resistant, and compatible with steam and gas plasma sterilisation.

	Product code	Description	Colour	Qty
	PEN001	marker pen	black	10 pcs
	PEN002	marker pen	red	10 pcs
	PEN003	marker pen	blue	10 pcs
I	05514	marker Steriliner	black	10 pcs
J	RT-100P	permanent marker	black	10 pcs
	RT-209P	permanent marker	black	12 pcs
K	09170	permanent marker	red	10 pcs
L	09088	permanent marker	black	10 pcs
M	05505	permanent marker	black	1 pc



Tags – Metal Tray Tags



A



B




C



D



FRONT



BACK

E

PRODUCT SHOWN ACTUAL SIZE

METAL TRAY TAGS

Stainless steel tags easily clip onto the side of a tray or container. Flat surface allows for placement of barcode, label, or other identifying instructions. Available in five sizes. Compatible with steam sterilisation. Quantities available in 1 pc and 10 pcs.

Product code	Size	Qty
A MT-0303	30 x 30mm	1 pc
B MT-0553	55 x 30mm	1 pc
C MT-0703	70 x 30mm	1 pc
D MT-0903	90 x 30mm	1 pc
E MT-0905	90 x 50mm	1 pc

STAINLESS STEEL TRAY TAG

Product code	Size	Qty
TAGMET7530	75 x 30mm	1 pc



TRANSPORT AND DISPOSAL

Our selection of Transport Boxes and Lidded Crates ensure that contaminated or cleaned surgical instruments can be transported effectively, with a reduction in damage to wrapped sets.

Transportation Boxes (UN 3291) are suitable for the transportation of Medical Waste N.O.S (Class 6.2, Classification code I3, (UN 3291) in accordance with Packing Instruction P621 defined in ANNEX A of the European Agreement concerning the International Carriage of Dangerous Goods by Road. Also compliant with section 2.1.5 of the agreement and facilitate the transportation of all used medical instruments or equipment for the purposes of decontamination, cleaning, sterilisation, or repair.

Lidded Crates are for the transportation of medical waste N.O.S. These crates include an integrated lid that can be sealed using a pin or sealing strip to comply with transportation regulations. The range includes four sizes and boasts ergonomic handles for simple carrying. The two larger sizes can accommodate both wrapped and unwrapped DIN size trays.

Transportation Boxes

Lidded Crates

Transportation Bags

Transportation Trays

Cable Ties

Please note: Unless specified, all products are not made with natural rubber latex.





Transport and Disposal – Transportation Boxes

TRANSPORTATION BOXES

Stackable boxes that can be custom printed with text on the box and lid. They can withstand washer disinfectors for decontamination purposes. BTB005 must be used with a trolley BTB005TROL.

BTB001	Suitable for Transporting Packing Group I dangerous substances
BTB002	Suitable for Transporting Packing Group II dangerous substances
BTB003	Suitable for Transporting Packing Group I dangerous substances
BTB004A	Suitable for Transporting Packing Group II dangerous substances
BTB005	Suitable for Transporting Packing Group II dangerous substances

	Product code	Weight	Capacity	Colour	Size	Qty
A	BTB001	0.6kg	4l	blue	external: 300 x 200 x 120mm	1 pc
B	BTB002	0.9kg	7l	blue	external: 300 x 200 x 170mm	1 pc
C	BTB003	1.4kg	16l	blue	external: 400 x 300 x 220mm	1 pc
D	BTB004A	3.3kg	49l	blue	external: 600 x 400 x 320mm	1 pc
	BTB005	10.2kg	170l	blue	external: 800 x 600 x 600mm	1 pc
E	RTB002	0.9kg	7l	red	external: 300 x 200 x 170mm	1 pc
F	RTB003	1.4kg	16l	red	external: 400 x 300 x 220mm	1 pc
G	RTB004A	3.3kg	49l	red	external: 600 x 400 x 320mm	1 pc

INNER BASKETS

Powder coated aluminum inner baskets designed to fit Transportation Boxes.

	Product code	Description	Qty
H	BTB002BASK	inner basket to fit BTB002 and RTB002	1 pc
	BTB003BASK	inner basket to fit BTB003 and RTB003	1 pc
	BTB004BASK	inner basket to fit BTB004A and RTB004A	1 pc

Transport and Disposal – Transportation Boxes



TRANSPORTATION BOX OPENER

Stainless steel opener suitable for all BTB and RTB boxes.

Product code	Description	Qty
A BTBOPENER	box opener for BTB and RTB boxes, stainless steel	1 pc

TRANSPORTATION BOX TROLLEYS

Painted steel trolley with wheels for BTB and RTB boxes. Suitable for washer disinfectors.

Product code	Description	Qty
B BTB004TROL	trolley to fit BTB004A and RTB004A	1 pc
C BTB005TROL	trolley to fit BTB005	1 pc



TRANSPORTATION BOX LABELS

A selection of labels for 'BIOHAZARD' or 'USED MEDICAL DEVICE' to be inserted in plastic sleeves of transportation boxes. Our biohazard labels conform to the requirements for a 6.2 hazard (infectious substances).

Product code	Description	Size	Qty
D BTBLABEL001	biohazard, non-laminated label	125 x 40mm	20 pcs
BTBLABEL003	biohazard, laminated label	125 x 40mm	20 pcs
BTBLABEL004	biohazard, laminated label	150 x 70mm	20 pcs
E BTB232SML	used medical device, laminated label	125 x 40mm	100 pcs
E BTB232LRG	used medical device, laminated label	150 x 70mm	100 pcs



Transport and Disposal – Lidded Crates

A



LIDDED CRATES

Our range of boxes for the transportation of medical waste N.O.S. These crates include an integrated lid that can be sealed using a sealing strip to comply with transportation regulations. The range comprises of four sizes and boasts ergonomic handles for simple carrying. They accommodate wrapped and unwrapped DIN trays.

Product code	Weight	Capacity	Size	Qty
UNLC20	1.68kg	20l	external: 400 x 300 x 275mm	1 pc
A UNLC25	2.18kg	25l	external: 500 x 300 x 274mm	1 pc
B UNLC45	3.26kg	45l	external: 600 x 400 x 275mm	1 pc
C UNLC60	4.20kg	60l	external: 600 x 400 x 374mm	1 pc

B



Product code	Description	Qty
UNLC-X	sealing strip for packing groups i (goods of high risk), ii (goods of medium risk), and iii (goods of low risk)	100 pcs

C





Transport and Disposal – Transportation Bags and Trays

TRANSPORTATION BAGS

Our high-quality polyethylene transportation bags protect sterile packs from rigorous handling and provide an effective dust cover. With tamper proof seals, the bags also have ventilation holes to prevent condensate spreading in a humid environment.

	Product code	Size	Qty
A	TBAG001	380 x 760mm	250 bags
	TBAG002	585 x 765mm	200 bags
	TBAG003	380 x 710mm	250 bags
	TBAG004	385 x 490mm	250 bags
	TBAG005	664 x 762mm	200 bags
	TBAG006	585 x 765mm	150 bags

TRANSPORTATION TRAYS

Transportation Trays reduce the risk of damage to wrapped packs due to poor handling and racking during transport and storage. Can be used with a strap or interlocked to form a lidded box.


	Product code	Material	Size	Qty
B	TTRAY002	polypropylene	500 x 270 x 70mm	1 pc
C	TTRAY001	cardboard	500 x 270 x 70mm	1 pc
	TTRAY003	cardboard	500 x 540 x 70mm	1 pc

TRANSPORTATION TRAY STRAPS

These straps are extremely strong with double welded joints. They also offer a safe means of lifting heavy trays without risk of damage to wraps. The straps are held in place with hook and loop fastenings. Suitable for steam sterilisation.

	Product code	Colour	Size	Qty
D	STR004	blue	1400mm	10 pcs





Transport and Disposal – Cable Ties

SECURITY CABLE TIES

Plastic cable ties in a range of colours for use during transportation and identification.

	Product code	Colour	Size	Qty
A	BTBREDTIE	red	3.6 x 145mm	100 pcs
B	BTBYELTIE	yellow	3.6 x 145mm	100 pcs
C	TIE001	blue	3.5 x 145mm	100 pcs
	TIE002	clear	2.5 x 141mm	100 pcs
D	TIE004	clear	3.6 x 300mm	100 pcs
E	TIE006	clear	1.5 x 125mm	500 pcs
F	TIE007	yellow	4.8 x 300mm	100 pcs
G	TIE008	green	4.8 x 300mm	100 pcs
H	TIE010	red	4.8 x 300mm	100 pcs

CABLE TIES

Plastic cable ties in a range of colours for use during transportation. Coloured cable ties help identification, along with numbered cable ties.

	Product code	Description	Color	Size	Qty
I	TIE003	individually numbered	green	2 x 200mm	100 pcs
	TIE005	individually numbered	red	2 x 200mm	100 pcs
J	TIE009	resealable cable tie, not autoclavable	red	120mm	1000 pcs

A Adhesive Remover 138-139
Amcor Ultra Pouches and Reels 72-73

B Barrier Wrap 50-57
Clinipak Choice Wrap 50-53
Clinipak Choice Wrap Dry 54-55
Reliance Solo Wrap 56-57
Breakaway Tags 188-189
Brush Assortments 34-35

C Cable Ties 222-223
Security Cable Ties 222-223
CBM Baskets 110-111
Mini and Micro Baskets 110-111
Sterilisation Basket System 110-111
CBM Containers 130-131
Channel Cleaning Brushes 02-21
Bone Reamer Brushes 16-17
Polypropylene Brushes 18-21
Colour Coding 144-147
Colour Curls 144-145
Fixing Ring 146-147
Container Filters 130-131
Container Labels 132-133

Containers and Accessories 128-133

Crepe Wrap 58-59
Reliance Crepe Wrap 58-59
Reliance Non-Woven Wrap 58-59
Corner Protectors 118-119
Card and Crepe Corner Protectors 118-119
Surgical Tray Corner Protectors 118-119

E Endoscope Cleaning Brushes 36-47
Albarran Brush 44-45
Brush Assortments 46-47
Premium Brushes 26-47
Pull-Through 42-43
Single Ended 40-41
Valve Brushes 44-45

F Flat Base Baskets 100-101
Flat Base DIN Baskets 100-101
Flat Caps 150-153
Flat Cap Assortments 152-153
Precision Flat Caps 152-153

G General Instrument Cleaning Brushes 30-33

I Identification Tags 204-205
Identification Autoclave Tags 204-205
Identification Security Tags 204-205

Instrument Cleaning and Maintenance 00-47

Instrument Holders 166-171
Clear Point Protectors 170-171
Instrument Sleeves 166-167
Needle Protectors 170-171
Point Protectors 168-169
Tip Protectors 170-171
Webmesh Guard 168-169

Instrument Identification 134-147

Instrument Marking Tape 140-143

Instrument Protection 148-177

K KeyDot 136-137

L Labels 124-127
Aluminium Label 126-127
Insertion Label 126-127
Label Holder 126-127
Stainless Steel Sieve Tray Labels 124-125
Lidded Crates 218-219
Liners 60-65
Millennium Blue Liners 60-61
Tray Liners 62-65

M Mesh Baskets 102-105
Mesh Baskets and Lids 102-103
Perforated Plate Mesh Trays 102-103
Metal Tray Tags 210-211
Stainless Steel Tray Tag 210-211
Micro Mesh Trays 104-105
Minimally Invasive Surgical Instrument (MIS) Caps 162-163

N Nic-Nac Baskets 84-85

O One-Loop Tags 200-201
One-Loop Tags Mesh Baskets 200-201
Osteotome Caps 158-159
Precision Osteotome Caps 158-159


P Paper Seals 132-133
Pipe Cleaners 22-25
Flexistem 24-25
Microclean 24-25
Plastic Sterilisation Trays 112-113
Sterilisation Trays 112-113
Wash Trays 112-113
Plastic Tray Tags 180-187
Tray Tag Clip 184-185
Tray Tag Ring 184-185
Protection Boards 176-177
Sheets 176-177
Silicone Boards 176-177



Protective Mats 122-125
Absorbent Protective Mats 122-123
Grid Inserts and Mats 124-125
Silicone Mesh Cover 124-125
Silicone Nub Mats 124-125

R Retractor Caps 160-161
Gelpi Retractor Caps 160-161
Self-Retaining Retractor Caps 160-161
Roll Tape 138-139
Round Caps 154-157
Precision Round Caps 156-157

S Self-Sealed Pouches 76-77
Self-Tie Tags 202-203
Sieve Trays 86-99
Compartmental sieve for small parts 92-93
Sieve for small parts 86-91
Sieve Trays and Lids 96-97
Sieve Tray for small parts and lids 92-93
Sieve Tray Sets 98-99
Sieve Trays with Drop Handles 98-99





Silicone Inserts 116-117
 Skin Hook Cap 162-163
 Soaker Sheets 120-121
 Stainless Steel Baskets 106-109
 Premium Stainless Steel Baskets 108-109
 Strindberg Trays 108-109
 Trays and Cassettes 108-109
 Stainless Steel Scope Baskets 114-115

Sterilisation Packaging 48-79

Sterilisation Tape 66-69
 Adhesive Tapes 66-67
 Indicator Autoclave Tapes 66-67
 Latex-Free Tapes 68-69
 SMS Autoclave Tape 68-69
 Steriliser Cleaning Brushes 32-33
 Pistol-Grip Brush 32-33

Tags 178-211

Tapered Caps 164-165
 Toothbrush-Style Cleaning Brushes 26-31

Transport and Disposal 212-223

Transport Wrap 70-71
 Clinipak Transport Wrap 70-71
 Clinipak Transport Wrap Dry 70-71

Transportation Bags 220-221
 Transportation Boxes 214-217
 Inner Baskets 214-215
 Transportation Box Opener 216-217
 Transportation Box Labels 216-217
 Transportation Box Trolleys 216-217
 Transportation Trays 220-221
 Transportation Tray Straps 220-221

Trays and Accessories 80-127

Trinket Baskets 82-83
 Tube and Cable Binders 172-175
 Silicone Cables 174-175
 Silicone Tubes 172-173
 Silicone String 174-175
 Tube and Cable Binders 174-175
 Twin Caps 162-163
 Twin Cap Assortment 162-163

W Washer Tags 206-207
 Wipak Pouches and Reels 74-79
 Write-On Tags 208-209
 Marker Pens 208-209



