

**DIATHERMY CLEANING BRUSH****Product Code(s)**

- DIA001 Diathermy Cleaning Brush

Unique stainless-steel brushes mounted onto a strong plastic support, which may be assembled as a double brush with automatic depth guard or as two single brushes with guard, for cleaning the tips of diathermy forceps whilst protecting the insulation.

One pack contains two brushes in different widths (10mm and 15 mm) complete with guards.  
Produced to a patented design, this brush can be separated easily for safe cleaning and disinfection.

**Materials**

Guards are made from glass-reinforced Nylon 6.  
Brushes are made from Stainless steel.  
Glue made from Neoprene-based, Latex-Free adhesive.

**Sizes**

- Diathermy Brushes come in 2 widths – 10mm and 15mm.
- Each pack contains 1 x 10mm brush plus 1 x 15mm brush.

**Cleaning**

Diathermy Brushes can be wash-decontaminated in a standard washer-disinfector. These are not suitable for steam sterilisation.