

MICROCLENE DISPENSER**Product Code(s)**

- FLEX025

A simple-to-use dispenser for cutting lengths of Microclene small-bore cleaner. The dispenser ensures that lengths of cleaner are readily available for use.

**Contents**

A FLEX025 Microclene dispenser includes the following items:

1. Plastic Dispenser
2. 1 pair of cutters
3. FLEX021
4. FLEX022
5. FLEX023
6. FLEX024
7. Rawl Plugs and screws

Materials

Dispenser	– Perspex {Poly(methyl methacrylate)}
FLEX021/22/23/24	– Steel, Polyester and Nylon
Rawl Plugs	– Plastic + steel

Dimensions and weights

Dimensions : 200x 190 x 180mm

Weight : 2.1kg (Dispenser only)

Refill Pack

FLEX026 – contains 1 each of FLEX021, FLEX022, FLEX023, FLEX024