


**FLEXITEM DISPENSER****Product Code**

- FLEX008

A simple-to-use dispenser for cutting lengths of Flexistem small-bore cleaner. The dispenser ensures that lengths of cleaner are readily available for use.

**Contents**

A FLEX008 Flexistem dispenser includes the following items:

- Plastic Dispenser
- 1 pair of cutters
- FLEX001
- FLEX002
- FLEX003
- FLEX004
- FLEX005
- FLEX006
- FLEX007
- Rawl Plugs and screws

**Materials**

Dispenser	– Perspex (Poly(methyl methacrylate))
FLEX001/2/3/4/5/6/7	– Steel, Polyester and Nylon
Rawl Plugs	– Plastic + steel

**Dimensions and Weights**

Dimensions:	285 x 190 x 180mm
Weight:	2.4kg (Dispenser only)

**Refill Pack**

FLEX009 – contains 1 each of FLEX001, FLEX002, FLEX003, FLEX004, FLEX005, FLEX006, FLEX007