

FLEXISTEM AND MICROCLENE SMALL BORE CLEANERS

Product Code(s)

- FLEX001-033 (various types and sizes)



Flexistem and Microclene are unique single-use disposable cleaners for all small bore cannulae. Various diameters from 3mm to 15mms are available constructed from very soft to hard bristles inserted between two soft mild steel twisted wires which assist in following contours while resisting buckling.

PropertiesFlexistem

Code	Fibre	Diameter	Length	Type
FLEX001	Cotton	3mm	10m	Soft Cleaner
FLEX002	Cotton & Nylon	3mm	10m	Super Bristle
FLEX003	Cotton & Nylon	3mm	15cm	Double Ended
FLEX004	Polyester	9mm	30cm	Fleece
FLEX005	Polyester	12mm	30cm	Fleece
FLEX006	Polyester & Nylon	15mm	30cm	Bristle Stem
FLEX007	Polyester & Nylon	15mm	5m	Bristle Stem
FLEX007A	Polyester & Nylon	15mm	5m	Bristle Stem
FLEX007B	Polyester & Nylon	15mm	66.5cm	Bristle Stem
FLEX033	Polyester	3mm	30cm	Fleece

Microclene

Code	Fibre	Diameter	Length	Type
FLEX021	Polyester & Nylon	3mm	7.5m	Super Fine
FLEX022	Polyester & Nylon	6mm	7.5m	Medium Fine
FLEX023	Polyester & Nylon	3mm	7.5m	Super Brush
FLEX024	Polyester & Nylon	6mm	7.5m	Medium Brush

Flexistem FLEX007A and FLEX007B has a wavy bristle configuration to reduce buckling which can occur in long instruments.