

GOWN TIES

Product Code(s)

- GOW002 Crocodile Gown Tie

Gown Ties are an effective way to hold gowns together to facilitate aseptic donning procedures. The unique pattern allows for two ties to be put into both ends or through the centre for simple tie shortening.



Material

Crocodile Gown Ties are manufactured from a heavy weight, pure wet strengthened white card developed to withstand steam sterilisation and washing in the event of accidental re-processing. The card complies with the chloride and sulphate contents within EN 868 part 2 and contains no optical brightening agents. Crocodile Gown Ties are serrated to provide increased security, particularly when used on synthetic materials.

Standard Properties

| | | GOW002 |
|----------------------------|------------------|---------------|
| Grammage | g/m ² | 281 |
| Thickness | um | 350 |
| Size | mm | 50 x 160 |
| Weight | g | 2.2 |
| pH of Hot Extract | | 7.0 |
| Chlorides | % | 0.05 max |
| Sulphates | % | 0.25 max |
| Optical Brighteners | | nil |