

## INSTRUMENT/TRAY/CONTAINER LOCKS

## Product Code(s)

➤ LOCK994G	Instrument/Tray Lock ("Good")
➤ LOCK994R	Instrument/Tray Lock ("Repair")
➤ LOCK994S	Instrument/Tray Lock ("Sharpen")
➤ LOCK3030YELL1000	Container Lock, blank with no writing

A range of special locking tags designed for marking instruments, trays and containers.

The LOCK994 locks are printed with "Good", "Repair" or "Sharpen", while the container lock LOCK3030YELL1000 is supplied blank with no printing on.

The LOCK994 products are specifically designed for use with instruments, for circulation between the user, the CSSD and repair shops.

## Material

All Instrument/Tray/Container locks are manufactured from polypropylene. These locks are free of any natural rubber latex.


## Colours

LOCK994G – Green

LOCK994R – Red

LOCK994S – Yellow

LOCK3030YELL1000 - Yellow



## Sterilising

Instrument/Tray/Container locks may be sterilised with steam (up to 137°C/5 min), ethylene oxide and plasma sterilisation systems.