

**MARKER PENS****Product Code(s)**

- PEN001
- PEN002
- PEN003



Three different coloured permanent marker pens whose ink is suitable for marking items which are to be sterilised in steam or ethylene oxide. Containing instantly touch-dry ink, these pens are suitable for writing on the film surface of “see through” packaging and on Tray Tags.

**Material**

The ink is a mixture of dyes, solvents and resins.  
The solvent is 1-methoxy 2-propanol.

**Colours**

|        |       |
|--------|-------|
| PEN001 | Black |
| PEN002 | Red   |
| PEN003 | Blue  |

**Sterilising**

The inks will withstand repeated steam autoclave cycles (up to 137°C/4min), though a gradual fading of the red and blue inks will occur.

The black ink is practically unaffected, even up to 200 cycles.

Marker pens are also suitable for marking items sterilised by ethylene oxide.

Marker pens are not suitable for use in plasma sterilisation systems.