

INSTRUMENT MARKING TAPE

Product Code(s)

- | | |
|-------------------|-----------------------|
| ➤ MKT001 – MKT004 | Ring Binders |
| ➤ MKT005 – MKT015 | Plain colours |
| ➤ MKT016 – MKT025 | Zebra-striped colours |
| ➤ MKT031 & MKT033 | Tiger-striped colours |

Latex-free marking tapes which allow the effective marking of instruments so that they may be returned to the same department or surgeon.

Available in a range of colour-fast plain and zebra-striped colours.

**Material**

Tape material 66µ calendered, cadmium-free polyvinyl chloride (PVC) based film.
 Adhesive Permanent, acrylic adhesive.
 Backing 135 g/m² siliconised, white clay-coated kraft paper.

Properties

Tensile strength	>1.5 kN/m
Peel adhesion against stainless steel	>15 N/25mm
Solvent and chemical resistance	Resistance to most oils and greases, fuel Aliphatic solvents, weak acids and alkalis.

Colours

See attached table

Sterilisation

Instrument Marking Tape is compatible with steam sterilisation up to 137°C for up to 4 min.
 Instrument Marking Tape is not compatible with any other type of sterilisation modality.

MARKING TAPE COLOURS

CODE	PRODUCT
MKT001	Ring Binder with Index only
MKT002	Ring Binder with index and 1 sheet of each plain colour
MKT003	Ring Binder with index and 1 sheet of each zebra colour
MKT004	Ring Binder with index and 1 sheet of each plain and zebra colour
MKT005	Plain Purple
MKT006	Plain White
MKT007	Plain Beige/Tan
MKT008	Plain Yellow
MKT009	Plain Orange
MKT010	Plain Red
MKT011	Plain Blue
MKY012	Plain green
MKT013	Plain Silver
MKT014	Plain Brown
MKT015	Plain Black
MKT016	Zebra Beige/Tan
MKT017	Zebra Yellow
MKT018	Zebra Orange
MKT019	Zebra Red
MKT020	Zebra Blue
MKT021	Zebra Green
MKT022	Zebra Silver
MKT023	Zebra Brown
MKT024	Zebra Black
MKT025	Zebra Purple
MKT031	Tiger Orange
MKT033	Tiger Green