

DESCRIPTION: MEDICA SCRUB BRUSHES**Product Code:**

➤ MED001



Use: 1) To disinfect the hands of medical staff before a surgical procedure.
2) To clean the intact skin of a patient.

Material: Dry sterile (sterilised with ethylene oxide) and is latex-free.

Cleaning: MED001 brush is single use.

Standard Properties: 1. Store the product at 5 – 30 °C, protected from direct sunlight.
2. The expiry date on the brush is 5 years after manufacturing date.
3. Each carton contains 15 boxes of 25 brushes – a total of 375 brushes.

Sterilising: MED001 Dry Sterile Brush is sterilised with ethylene-oxide.