

OR SOAKER SHEETS

Product Code(s)

➤ ORSOAK1000

For use in operating theatres to absorb fluid spillage during the procedure.

The non-slip backing prevents movement when laid on the floor, while the upper surface absorbs up to 1 litre/m².

Particularly useful in urology and obstetrics cases.



Material

OR soaker sheets consist of a thermally-bonded 3-layer structure as follows:

Top Layer	Thermally-bonded printed PE/Kraft with non-slip coating
Middle layer	Proprietary blend of cellulose and bonding fibres
Bottom Layer	High strength porous non-woven material

Sheet size

ORSOAK1000 600 x 1000mm

Fluid Absorbency

ORSOAK1000 1 litres per m²
600 ml per sheet