

LATEX FREE SILICONE BANDS

Product Code(s)

- BAND001
- BAND002
- BAND003
- BAND004

Blue, latex-free, elastic bands can be used in numerous ways where natural latex is a problem. E.g. Hand surgery, Securing Instruments, Flap Retention, Bundling of Arterial forceps etc.



Materials

Latex free bands are manufactured from a silicone elastomer.

A very small amount of chalk (calcium carbonate) used as a mould-release agent may be present on the surface of the bands

Sizes

PRODUCT CODE	SIZE	PACK QUANTITY
BAND001	60 X 1.5 X 3 mm	1000/pack
BAND002	75 X 1.5 X 3 mm	800/pack
BAND003	60 X 1.5 X 3 mm	200/pack
BAND004	75 X 1.5 X 3 mm	200/pack

Sterilisation

Latex Free Bands are compatible with steam sterilisation (typically 134-137°C for 4 min) and ethylene oxide sterilisation.

Latex Free Bands are rated for single-use only and should not be sterilised more than once.