

Description: SILICONE INSERTS**Product Code(s)**

➤ See Table below

A range of blue silicone inserts designed to fit Plastic Sterilisation Trays. They increase the flexibility of the trays by allowing a range of instruments to be held in the tray during sterilisation

Material


All inserts are manufactured from medical grade silicone polymer.

All inserts are latex free.

Care should be taken during any cleaning process as the inserts can be stretched out of shape.

Sterilisation

Silicone inserts are suitable for use in steam, plasma and ethylene oxide sterilisation processes.

**SILICONE INSERTS**

| CODE | PRODUCT |
|----------|--|
| A237591K | Silicone bar, 1 slot open end round, 5mm, L 59mm x 15mm |
| A237592K | Silicone bar, 2 slot open end round, 5mm, L 59mm x 15mm |
| A237593 | Silicone bar, 3 slot straight cut, L 59mm x 15mm |
| A237593K | Silicone bar, 3 slot open end round, 5mm, L 59mm x 15mm |
| A237596 | Silicone bar, 6 slot straight cut, L 59mm x 15mm |
| A2616901 | Silicone bar, 2 slot open end round, 9mm and 12mm, L 67mm x 29mm |
| A2616902 | Silicone bar, 2 slot open end round, 9mm, L 67mm x 29mm |
| A327593 | Silicone bar, 3 slot open end round, 7mm, L 93mm x 31mm |
| A327593L | Silicone bar, 3 slot open end round, 15mm, 7mm and 12mm, L 93mm x 31mm |
| A32759P1 | Silicone bar, 1 slot open end round, 18mm, L 81mm x 31mm |
| A32759P2 | Silicone bar, 1 slot open end round, 7mm, L 81mm x 31mm |
| A3875912 | Silicone bar, 12 slot straight cut, L 96mm x 15mm |
| A387592 | Silicone bar, 2 slot open end round, 5mm, L 96mm x 15mm |
| A387593 | Silicone bar, 3 slot open end round, 5mm, L 96mm x 15mm |
| A387594 | Silicone bar, 4 slot open end round, 5mm, L 96mm x 15mm |
| A387595 | Silicone bar, 5 slot open end round, 5mm, L 96mm x 15mm |
| A387596 | Silicone bar, 6 slot straight cut, L 96mm x 15mm |