


**SILICONE PIN MAT****Product Code(s)**

➤ DINPIN001

A large silicone Pin Mat which may be cut up into a variety of sizes. Ideal for use in any tray where a bespoke mat is not available. Provided with holes in the bottom of the mat for rapid draining.

**Materials**

DINPIN001 mats are made of a silicone polymer.  
They are latex free.

**Size**

DINPIN001 -500 x 300mm  
PIN SIZE -13.3mm

**Washing and Sterilising**

Silicone Pin Mats are suitable for use in steam, ethylene oxide and plasma sterilisation systems.  
Silicone Pin Mats are suitable for use in all known wash decontamination procedures.