


**Description: CABLE TIES**

Single-use autoclavable cable ties.

**Product Code(s)**

- TIE001
- TIE002
- TIE003 individually numbered, Green
- TIE004
- TIE005 individually numbered, Red
- TIE006 Secure a tie
- TIE007 Yellow Security Tag Tie
- TIE008 Green Security Tag Tie
- TIE009 Re-sealable cable Tie
- TIE010 Red Security Tag Tie

**Material**

TIE001-TIE008 and TIE010 are manufactured from nylon.

TIE009 is manufactured from polyethylene

**Size**

TIE001 – 150mm x 3.5mm

TIE002 – 144mm x 2.5mm

TIE003 – 200mm x 2mm

TIE004 – 300mm x 3.6mm

TIE005 – 200mm x 2mm

TIE006 – 125mm x 1.5mm

TIE007 – 300mm x 4.8mm

TIE008 – 300mm x 4.8mm

TIE009 – 120mm x 3.5mm

TIE010 – 300mm x 4.8mm

The TIE003 and TIE005 versions are individually numbered.

TIE003 and TIE005 are fitted with a quick release mechanism.

**Colour**

TIE001 – Blue

TIE002 – Clear

TIE003 – Green

TIE004 - Clear

TIE005 – Red

TIE006 – Clear

TIE007 – Yellow

TIE008 – Green

TIE009 – Red

TIE010 – Red

**Sterilising Compatibility**

As TIE001-TIE008 are single-use products they are suitable for use in steam sterilisation processes up to 137°C and hydrogen peroxide (0.5%) plasma sterilisation systems.

TIE009 and TIE010 are not suitable for sterilisation.